

DENTALEZ

Connecting innovation to care



2025

Equipment & Utility

TABLE OF CONTENTS

EQUIPMENT

6400 Chair Mount Packages	8	Rear Mount Side Delivery Systems	67
Pivot Package Standard IC Control Head	8	Cabinetry Mounting Configurations	69
Sidebox Standard IC Control Head	9	7020PRO Duo Swing Delivery Systems	70
6400 Standard Package Ordering Guide	10	7021PRO Duo Swing Delivery Systems	72
6400 Designer Package Ordering Guide	14	7010 Duo Swing Delivery Systems	74
3900 Pivot Chair Mount Packages	22	3486 Swing Delivery Systems	76
Standard IC Control Head	22	Storage Organizer, Trays and Monitor Mount	78
Euro Control Head	23		
3900 Fixed Chair Mount with Sidebox Packages	24	Ancillary Ordering and Options	79
Standard IC Control Head	24	Standard IC Delivery Unit	79
Euro Control Head	25	Euro Delivery Unit	80
3900 Package Ordering Guide	26	Fusion Delivery Unit	81
Fusion Packages	38	Handpiece Bundles	82-83
Pivot Chair Mount	39	Control Head Accessories	84
Fixed Mount with Sidebox	40	Ancillary Products	84
Rear Mount Side Delivery	41	Simplified Cuspidors and Assistant's Instrumentation	87
Fusion Package Ordering Guide	42	Fixed Mount with Sidebox, External Umbilical	88
Surface Treatments	46	Fixed Chair Post Mount, External Umbilical	92
Designer Friendly Colors and Surfaces	47	Rear Beneath Shelf or Work Surface Mount	93
Chairs and Stools	49	Simplified Rear and Side Delivery Systems	
6400 Patient Chair	52	Side Cabinet and Wall Mount	94
3900 Hydraulic Patient Chair	53	Rear or Side Beneath Counter Mount	95
3900 Chair Upholstery Options	54	Rear Beneath Counter Mount	96
3800 Hydraulic Patient Chair	55	Panel Mounts	98
Chair Accessories	56	Rear or Side Mount Compact Controls	99
Economy and Deluxe Stools	57	Mobile Carts	100
Comfort Stools	58	Factory Options and Parts	101
Lights and Monitor Mounts	59	General Information	107
Post Mount Lights	61	Dimensions	107
Track Mount Lights	62	Equipment Maintenance	122
Ceiling Mount Lights	63	Equipment Terms and Warranty	123
Cabinet and Wall Mount Lights	64		
Monitor Mount Assemblies	65		

UTILITY

Compressors	128	Maintenance	137
Ramvac Compressors	128	Compressor and Vacuum Maintenance	137
Dry Vacuums	131	Utility Room Accessories	138
Dry Vacuums	131	Remote Panels	138
Badger LF	133	Master Water Control Valve	138
Tanks	134	Steel Frame Stands	138
15 Gallon and 30 Gallon Tanks	134	Amalgam Separator	139
InfiniTank	134	VACHECK / FLOWCHECK	139
		SlugBuster	139
Wet Vacuums	135	Quick Reference Guide	140
Barracuda	135	Utility Warranty	142
		DSX Sales Team	143



Connecting innovation to care

DENTALEZ is a proven leader, dedicated to providing innovative and trusted products, service, and support. We empower healthcare professionals to maximize their operations for increased patient care.

Forest equipment has reached new heights in innovation with the 6400 chair and its features. Real solutions are provided with everyday challenges in the dental world being "the new standard in smart". We maximize revenue opportunities with procedure efficiency data, minimize maintenance, and easily monitor the chair's performance with metrics from built-in Aeras Intel technology—all with no extra charge.

- [Forest 6400 Packages](#)
- [Forest 6400 Ordering Guides](#)
- [Forest 3900 Packages](#)
- [Forest 3900 Ordering Guide](#)
- [Forest Fusion Packages](#)
- [Forest Fusion Ordering Guide](#)
- [Surface Treatments](#)
- [Chairs and Stools](#)
- [Lights and Monitor Mounts](#)
- [Delivery Systems 7000 Series](#)
- [Delivery Systems](#)
- [Simplified](#)
- [Factory Options and Parts](#)
- [General Information](#)
- [Equipment Warranty](#)
- [Compressors](#)
- [Dry Vacuums](#)
- [Tanks](#)
- [Wet Vacuums](#)
- [Utility Room Accessories](#)
- [Utility Quick Reference Guide](#)
- [Utility Warranty](#)
- [Sales Reps](#)

Forest

PACKAGES



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

OPERATORY PACKAGES AVAILABLE WITH THREE CHAIR OPTIONS



*3800 Chair Package (Sidebox)
Forest Grey Base Color*

*3900 Chair Package (Pivot)
Shadow Base Color*




*6400 Chair Package (Sidebox)
Forest Grey Base Color*

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps


Pivot Chair Mount

Pivot Chair Mounts easily move left to right on large roller bearings to accommodate left and right-handed practitioners. Unlike a Fixed Chair Mount delivery system, the mounting arm comes from underneath the patient rather than over the top, making it less obtrusive and more patient-friendly. Pivot Assistant's Instrumentation, Operatory Lights and Monitor Mounts are popular options to support your style of practice.



Fixed Chair Mount With Sidebox

Fixed Chair Mounts With Sidebox allow for a Cuspidor, Assistant's Instrumentation, a Monitor and Light to mount on a sidebox, which is large enough to accommodate a host of ancillary support products. Finish the aluminum side covers in one of our nine powder coat paint colors and choose a complementary carbon fiber trim to finish the top cap. If you use a cuspidor, select our beautiful high-gloss ceramic black for a customized and beautiful patient-friendly delivery system.



6400 Pivot Chair Mount Packages

6400 Chair Package, Standard IC Delivery, Pivot

Base Package



Popular Options



Dr's Deluxe Stool,
Ultraleather Upholstery

Part number: 6159C



E-Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



LED Operatory Light, 9073

Part number: 9073



Assistants Instrumentation
rear Pivot with touchpad

Part number: 1096

Designer Friendly Package As Shown

Description	Part Number
6400C Patient Chair	6400C
Delivery System w/Standard IC, Pivot	4486PI
LED Operatory Light, 9073	9073
Light Post, Pivot	1113-654
Chrome Accents	9000-998-CS
Assistant Hygiene Instrumentation	1408-CS

Popular Options

Description	Part Number
Dr. Deluxe Stool, Ultraleather Upholstery	6159C
Assistants Deluxe Stool No Back, Ultraleather Upholstery	6160C
E-Torque Electric HP System	STAR-ETORQUE
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Quick Disconnect for Water	0016-454-CS
Quick Disconnect for Air	0016-456-CS
Routing Valve Kit	1115-315-CS
Additional HVE w/Dry Oral Cup	1400-113-CS
Seamed memory comfort w/Heat and Massage option	3625-962
CareGlide	3625-970
Assistants Instrumentation, Rear Pivot with Touchpad	1096

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

6400 Pivot Chair Mount Packages

6400 Chair Package, Standard IC Delivery, Sidebox

Base Package



Popular Options



Dr's Deluxe Stool,
Ultraleather Upholstery

Part number: 6159C



E-Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



LED Operatory Light, 9073

Part number: 9073

Designer Friendly Package As Shown

Description	Part Number
6400C Patient Chair	6400C
Delivery System w/Standard IC, Sidebox	4485SI
LED Operatory Light, 9073	9073
Light Post, Sidebox	1110-900
Chrome Accents	9000-998-CS
E-Torque Electric HP System	STAR-ETORQUE
Swivel Tray	1110-468
Assistants Instrumentation with Deluxe Cuspidor	1577

Popular Options

Description	Part Number
Dr. Deluxe Stool, Ultraleather Upholstery	6159C
Assistants Deluxe Stool No Back, Ultraleather Upholstery	6160C
E-Torque Electric HP System	STAR-ETORQUE
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Additional HVE w/Dry Oral Cup	1400-113-CS
Seamed memory comfort w/Heat and Massage option	3625-962
CareGlide	3625-970

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

6400 STANDARD PACKAGE ORDERING GUIDE

SURFACE FINISH OPTIONS

ULTRALEATHER® 12 CLASSIC COLORS



ULTRALEATHER 5 PEARLIZED COLORS



POWDER COAT PAINT 2 COLORS



Seamless Upholstery



Seamed Upholstery

*Attention to Detail
at an Affordable Price!*

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

6400 OPERATORY PACKAGE - DENTAL CHAIR SELECTION

STEP
1

STEP 1: SELECT YOUR CHAIR - Select 1 chair and 1 powder coat paint color

Description	Part #	Retail Price	Qty.	Total Price	Powder Coat Paint Color
6400 Patient Chair	6400				<input type="checkbox"/> Forest Grey <input type="checkbox"/> Shadow

Standard chair colors come with the baseplate painted Shadow, all Designer Friendly color options come standard with the baseplate painted the same color as the chair.

STEP
2

STEP 2: SELECT YOUR UPHOLSTERY - Select 1 upholstery option and 1 upholstery color

Description	Part #	Retail Price	Qty.	Total Price	Upholstery Color (used as inner color if dual color upholstery is selected)
Seamed, Memory Comfort	3625-962X				Primary Upholstery Color Standard Ultraleather® Colors
Seamless, Memory Comfort	3625-966X				<input type="checkbox"/> Autumn Leaf <input type="checkbox"/> Baltic
Seamed, Memory Comfort with Heat & Massage Option	3625-962				<input type="checkbox"/> Boysenberry <input type="checkbox"/> Charcoal
Seamless, Memory Comfort with Heat & Massage Option	3625-966				<input type="checkbox"/> Doe <input type="checkbox"/> Diplomat Blue
					<input type="checkbox"/> Dove Grey <input type="checkbox"/> Papyrus
					<input type="checkbox"/> Parrot <input type="checkbox"/> Raven Wing
					<input type="checkbox"/> St. John <input type="checkbox"/> Stone
					<input type="checkbox"/> Midnight* <input type="checkbox"/> Oz*
					<input type="checkbox"/> Pewter * <input type="checkbox"/> Pharaoh*
					<input type="checkbox"/> Tin Man* *Pearlized

STEP 3: CHAIR CONTROL LOCATION

Maximum of 1 control can be mounted on the delivery unit, maximum of 2 on the chair (left & right), 4 total maximum number of controls per operatory

STEP
3

First Chair Control (Standard) Integrated in chair	Included
Second Chair Control - Delivery Unit	6414-046
Footswitch	3625-858

STEP
4

STEP 4: CHAIR OPTIONS

CareGlide	3625-970
Duplex with USB	3658-495
Surgical Integrated Dual Rails must select this along with either the IV Armrest or Bar Style Armrest	3680-444
IV Armrest with clamp	3625-964
Bar Style Armrest with clamp	3625-965
Second IV Armrest with clamp	3951-306
IV Bag Stand	3625-812
22" Black Strap	3802-512
54" Black Strap	3802-513
80" Black Strap	3802-514

STEP
5

STEP 5: STOOLS

Description	Part #	Retail Price	Qty.	Total Price
Doctor's Deluxe Stool	6159C			
Assistant's Deluxe Stool	6160C			

TOTAL CHAIR PRICE:



6400 OPERATORY PACKAGE - STANDARD IC DELIVERY UNIT

STEP
6

STEP 6: SELECT YOUR DELIVERY SYSTEM - Select 1 delivery system

Description	Part #	Retail Price	Qty.	Total Price
Standard IC Delivery Head, Pivot Chair Mount	4486PI			
Standard IC Delivery Head, Fixed Chair Mount with Sidebox	4485SI			
Euro Delivery Head, Pivot Chair Mount	4988PI			
Euro Delivery Head, Fixed with Sidebox	4987SI			
Plastic Finishes	<input type="checkbox"/> Forest Grey <input type="checkbox"/> Greywood <input type="checkbox"/> Laser Leather <input type="checkbox"/> Lemans <input type="checkbox"/> White			

STEP
7

STEP 7: CONFIGURE YOUR DELIVERY SYSTEM

Note: Please refer to ancillary and options graph in the catalog to determine number of ancillaries and options possible

Handpiece Positions / Comments

Default position placement
 Added/optional positions

Description	Part #	Retail Price	Qty.	Total Price	0	1	2	3	4	5	6
Air/Water Syringe (included standard)											
4-hole handpiece tubing (included standard)											
5-hole system (1 tubing included, system accommodates a 2nd tubing)	1700-021-CS										
- Second 5-hole tubing	0017-615-CS										
6-pin system (1 tubing included, system accommodates a 2nd tubing)	1700-019-CS										
- Second 6-pin tubing	0017-603-CS										
Star® ETorque Electric System, Factory Installed*	STAR-ETORQUE										
Star Blis-sonic™ SW Scaler (with 4 line swivel)*	STAR-BLIS-SW										
Star Blis-sonic K Scaler (with fixed back-end)*	STAR-BLIS-K										
Satelec SP Scaler*	SATELEC-S										
Satelec SP LED Scaler*	SATELEC-S-LED										
Cavitron Scaler*	CAVITRON										
Additional Handpiece Position, Wet	1100-189-CS										
Additional Handpiece Position, Air Drive Only	1100-190-CS										
Assistant's hygiene instrumentation, control head-mounted Select 1 position for the SE and 1 position for the HVE	1408-CS				SE						
					HVE						
Ancillary Shelf	0010-310										
Star Handpiece Bundle (Write-in PN and price)											
Quick Disconnect w/Flow Control, 1/4" (for water)	0016-454-CS										Included with 4485SI
Quick Disconnect w/Flow Control, 3/8" (for air)	0016-456-CS										Included with 4485SI
Swivel Tray	1110-468										
Less sidebox telescoping assistant's instrumentation**	1400-000-CS										4485SI only
Chrome Accents	9000-998-CS										
Additional Umbilical, Per Foot (max of 9 feet)	1115-999-CS										
Routing Valve Kit w/Water shutoff & Regulator, Factory inst.	1115-315-CS										Included with 4485SI
Routing Valve Kit for Cuspidor, Factory Installed	1115-316-CS										

* Ancillaries utilize an existing handpiece position - some integrated options may require purchase of Ancillary Shelf PN 0010-310

** Sidebox mounted telescoping assistant's instrumentation included standard with 4485SI

STEP
8

STEP 8: ASSISTANT'S INSTRUMENTATION AND CUSPIDORS

Rear Pivot Instrumentation with Integrated Touchpad	1096
Rear Pivot Instrumentation with Deluxe Cuspidor**	1590
Deluxe Cuspidor only, Rear Pivot**	1580
Deluxe Cuspidor, Sidebox Mounted (4485SI only)	1577

**Not available with 4485SI

STEP
9

STEP 9: ASSISTANT'S INSTRUMENTATION AND CUSPIDOR OPTIONS

Quick Disconnect w/Flow Control, 1/4" (For 1095 instrumentation)	0016-455-CS	Specify One	<input type="checkbox"/> Air	<input type="checkbox"/> Water
Additional Autoclavable HVE with 5' Tubing & Dry Oral Cup	1400-113-CS			
Hands-Free HVE Holder with HVE and 5' Tubing (Select PN) \$Colors (holder/tubing): GRY (white/grey), STG (white/sterling), BLK (black/black), SHW (black/shadow)			<input type="checkbox"/> 1400-279-STG	<input type="checkbox"/> 1400-279-BLK
			<input type="checkbox"/> 1400-279-SHW	<input type="checkbox"/> 1400-279-GRY
Self-Contained Clean Water System w/Quick Disconnect, 1.5L bottle (For 1095 Instrumentation only)	1700-319			
Tray Assembly (For 1095 Instrumentation only)	1110-680			
Deluxe Cuspidor Unit-Mounted Activation (For 1577, 1580, 1590)	1100-480-CS			

TOTAL STANDARD IC DELIVERY UNIT PRICE:

STEP
10

STEP 10: OPERATORY LIGHT - Select Light Head Mounting Option - 1 light and mounting option required



Mount Type	9078 LED				9073 LED			
	Part #	Retail Price	Qty	Total Price	Part #	Retail Price	Qty	Total Price
Chair Mount	9078				9073			
Ceiling Mount	9078-CL				9073-CL			
Wall Mount	9079W				9074W			
Cabinet Mount	9079C				9074C			
Track Mount	9078T				9073T			

STEP 11: OPERATORY LIGHT CONFIGURATION AND OPTIONS

Configure based on the light-mounting option selected above

STEP
11

Chair Mounted Light Options (Pivot, Sidebox, and Fixed)

Description	Part #	Retail Price	Qty	Total Price
20" Rigid Arm in Lieu of Standard 32" Arm	1110-227-CS			
Light Posts (Select 1)				
Sidebox, 44" Post	1110-900			
Sidebox w/Monitor Cord slot, 44"	1110-170			
Pivot without chair adapter (adapter included with unit purchase)	1113-654			
Pivot with chair adapter (when chair-mounted unit not ordered)	1113-654-CA			
Monitor Mount Assemblies				
Ring Post Mount Assembly (Must pair with either PN 1110-170 or 1110-045 above) This option allows you to have both a light and monitor mounted to your chair	1300RM			
Top Post Mount Assembly (Must pair with one of the post options below) This option allows a monitor instead of a light on your chair	1300TPM			
Sidebox: Monitor Post, Dual Slot, Fixed Sidebox Mount, 30"	1331-998			
Pivot: Monitor in Lieu of Light Post, Pivot	1331-654			
Pivot: Monitor in Lieu of Light Post, Pivot with Chair Adapter	1331-654-CA			

Ceiling Mounted Light Configuration (Select Mounting Height and Post Length)

Mounting Height (Required) 8' 8'6" 9' 9'6" 10' 10'6" 11' Custom Height, Specify:Post Length (Required) 12" 18" 24" 30" 36" 42" 48" Custom Length (max 48"), Specify:

Cabinet and Wall Mounted Light Configuration and Options

Description	Part #	Retail Price	Qty	Total Price
Arm Height (Required): <input type="checkbox"/> Standard <input type="checkbox"/> Tall	N/A			
Specify Cabinet or Wall Mounting Height (write-in):	N/A			
Horizontal Top Cabinet Mount in Lieu of Internal Mount	1110-881-CS			
Vertical Exterior Cabinet Mount in Lieu of Internal Mount	1110-886-CS			

Track Mounted Light Configuration and Options

Monitor Mount Assembly (specify monitor post length in next section)	1305T
S Video Cable, Pre-Installed	0014-282
RCA Cable, Pre-Installed	0014-283

Select Track Mount Light Post Lengths

Mounting Height (Required) <input type="checkbox"/> 8' <input type="checkbox"/> 8'6" <input type="checkbox"/> 9' <input type="checkbox"/> 9'6" <input type="checkbox"/> 10' <input type="checkbox"/> 10'6" <input type="checkbox"/> 11' <input type="checkbox"/> Custom Height, Specify:	
Light Post Length (Required) <input type="checkbox"/> 7.25" <input type="checkbox"/> 13.75" <input type="checkbox"/> 19.75" <input type="checkbox"/> 25.75" <input type="checkbox"/> 31.75" <input type="checkbox"/> 37.75" <input type="checkbox"/> 43.75" <input type="checkbox"/> Custom Length (max 43.75"), Specify:	
Monitor Post Length (Required if 1305T selected) <input type="checkbox"/> 20.25" <input type="checkbox"/> 26.25" <input type="checkbox"/> 32.25" <input type="checkbox"/> 38.25" <input type="checkbox"/> 44.25" <input type="checkbox"/> 50.25" <input type="checkbox"/> 56.25" <input type="checkbox"/> Custom Length (max 56.25"), Specify:	

Forest

DENTALEZ

6400 DESIGNER PACKAGE ORDERING GUIDE




- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

SURFACE FINISH OPTIONS

ULTRALEATHER® 12 CLASSIC COLORS



ULTRALEATHER® 5 PEARLIZED COLORS



POWDER COAT PAINT 9 COLORS



* Designer Friendly Colors

CARBON FIBER 3 PATTERNS



Customize your treatment rooms with DENTALEZ!

Powder coat paint your drawer and door pulls, grab handles, trash receptacles, or other metal surfaces in your treatment room to match your delivery system paint color.

Source custom powder coat paint colors to match your provided sample.

Special order any UltraLeather® upholstery color shown on the Ultrafabrics® website (from original or pearlized palette).



Contact A Sales Rep



Seamless Upholstery



Seamed Upholstery

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

6400 OPERATORY PACKAGE - DENTAL CHAIR SELECTION

STEP 1: SELECT YOUR CHAIR - Select 1 chair and 1 powder coat paint color

STEP
1

Description	Part #	Retail Price	Qty.	Total Price	Powder Coat Paint Color														
6400 Patient Chair, Designer Friendly Includes any primary and secondary paint color on delivery unit and light	6400C				<p>Designer Friendly color options come standard with the baseplate painted the same color as the chair.</p> <p>Primary Color Secondary Color</p> <table><tbody><tr><td><input type="checkbox"/> Arctic Silver</td><td><input type="checkbox"/> Arctic Silver</td></tr><tr><td><input type="checkbox"/> Bronze</td><td><input type="checkbox"/> Bronze</td></tr><tr><td><input type="checkbox"/> Charcoal</td><td><input type="checkbox"/> Charcoal</td></tr><tr><td><input type="checkbox"/> Copper</td><td><input type="checkbox"/> Copper</td></tr><tr><td><input type="checkbox"/> Midnight</td><td><input type="checkbox"/> Midnight</td></tr><tr><td><input type="checkbox"/> Pewter</td><td><input type="checkbox"/> Pewter</td></tr><tr><td><input type="checkbox"/> White</td><td><input type="checkbox"/> White</td></tr></tbody></table> <p>Select the same color for primary and secondary if you want the whole chair to be the same color.</p>	<input type="checkbox"/> Arctic Silver	<input type="checkbox"/> Arctic Silver	<input type="checkbox"/> Bronze	<input type="checkbox"/> Bronze	<input type="checkbox"/> Charcoal	<input type="checkbox"/> Charcoal	<input type="checkbox"/> Copper	<input type="checkbox"/> Copper	<input type="checkbox"/> Midnight	<input type="checkbox"/> Midnight	<input type="checkbox"/> Pewter	<input type="checkbox"/> Pewter	<input type="checkbox"/> White	<input type="checkbox"/> White
<input type="checkbox"/> Arctic Silver	<input type="checkbox"/> Arctic Silver																		
<input type="checkbox"/> Bronze	<input type="checkbox"/> Bronze																		
<input type="checkbox"/> Charcoal	<input type="checkbox"/> Charcoal																		
<input type="checkbox"/> Copper	<input type="checkbox"/> Copper																		
<input type="checkbox"/> Midnight	<input type="checkbox"/> Midnight																		
<input type="checkbox"/> Pewter	<input type="checkbox"/> Pewter																		
<input type="checkbox"/> White	<input type="checkbox"/> White																		

STEP 2: SELECT YOUR UPHOLSTERY - Select 1 upholstery option and 1 upholstery color

STEP
2

Description	Part #	Retail Price	Qty.	Total Price	Upholstery Color (used as inner color if dual color upholstery is selected)
Seamed, Memory Comfort	3625-961				
Seamless, Memory Comfort	3625-823A				
Seamed, Memory Comfort with Heat & Massage Option	3625-962				
Seamless, Memory Comfort with Heat & Massage Option	3625-966				
Dual Upholstery Color	3907-999-DC				Primary Upholstery Color Standard Ultraleather® Colors
Secondary Color:					<input type="checkbox"/> Autumn Leaf <input type="checkbox"/> Baltic <input type="checkbox"/> Boysenberry <input type="checkbox"/> Charcoal <input type="checkbox"/> Doe <input type="checkbox"/> Diplomat Blue <input type="checkbox"/> Dove Grey <input type="checkbox"/> Papyrus <input type="checkbox"/> Parrot <input type="checkbox"/> Raven Wing <input type="checkbox"/> St. John <input type="checkbox"/> Stone <input type="checkbox"/> Midnight* <input type="checkbox"/> Oz* <input type="checkbox"/> Pewter * <input type="checkbox"/> Pharaoh* <input type="checkbox"/> Tin Man* *Pearlized
Contrast Stitch	3907-999-DS				

Thread Color:

Choose your thread color from available
upholstery color along with Gold, White or Red.

Upholstery:

Primary

Secondary



Designer Covers:

Primary

Secondary

6400 OPERATORY PACKAGE - DENTAL CHAIR SELECTION (CONTINUED)

**STEP
3**

STEP 3: CHAIR CONTROL LOCATION - Select at least 1 chair control

Maximum of 1 control can be mounted on the delivery unit, maximum of 2 on the chair (left & right), 4 total maximum number of controls per operatory

Description	Part #	Retail Price	Qty.	Total Price
First Chair Control (Required) Integrated in chair	Included			
Delivery Unit Control	6414-046			
Footswitch	3625-858			

**STEP
4**

STEP 4: CHAIR OPTIONS

CareGlide	3625-970
Duplex with USB	3658-495
Surgical Integrated Dual Rails must select this along with either the IV Armrest or Bar Style Armrest	3680-444
IV Armrest with clamp	3625-964
Bar Style Armrest with clamp	3625-965
Second IV Armrest with clamp	3951-306
IV Bag Stand	3625-812
22" Black Strap	3802-512
54" Black Strap	3802-513
80" Black Strap	3802-514

**STEP
5**

STEP 5: STOOLS

Description	Part #	Retail Price	Qty.	Total Price
Doctor's Deluxe Stool	6159C			
Assistant's Deluxe Stool	6160C			

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

6400 OPERATORY PACKAGE - STANDARD IC DELIVERY UNIT

STEP
6

STEP 6: SELECT YOUR DELIVERY SYSTEM - Select 1 delivery system

Description	Part #	Retail Price	Qty.	Total Price
Standard IC Delivery Head, Pivot Chair Mount	4486PI			
Standard IC Delivery Head, Fixed Chair Mount with Sidebox	4485SI			
Euro Delivery Head, Pivot Chair Mount	4988PI			
Euro Delivery Head, Fixed with Sidebox	4987SI			
Plastic Finishes	<input type="checkbox"/> Forest Grey <input type="checkbox"/> Greywood <input type="checkbox"/> Laser Leather <input type="checkbox"/> Lemans <input type="checkbox"/> White			

STEP 7: CONFIGURE YOUR DELIVERY SYSTEM

Note: Please refer to ancillary and options graph in the catalog to determine number of ancillaries and options possible

Description	Part #	Retail Price	Qty.	Total Price	0	1	2	3	4	5	6
					Default position placement	Added/optional positions					
Air/Water Syringe (included standard)											
4-hole handpiece tubing (included standard)											
5-hole system (1 tubing included, system accommodates a 2nd tubing)	1700-021-CS										
- Second 5-hole tubing	0017-615-CS										
6-pin system (1 tubing included, system accommodates a 2nd tubing)	1700-019-CS										
- Second 6-pin tubing	0017-603-CS										
Star® ETorque Electric System, Factory Installed*	STAR-ETORQUE										
Star Blis-sonic™ SW Scaler (with 4 line swivel)*	STAR-BLIS-SW										
Star Blis-sonic K Scaler (with fixed back-end)*	STAR-BLIS-K										
Satelec SP Scaler*	SATELEC-S										
Satelec SP LED Scaler*	SATELEC-S-LED										
Cavitron Scaler*	CAVITRON										
Additional Handpiece Position, Wet	1100-189-CS										
Additional Handpiece Position, Air Drive Only	1100-190-CS										
Assistant's hygiene instrumentation, control head-mounted Select 1 position for the SE and 1 position for the HVE	1408-CS				SE						
					HVE						
Ancillary Shelf	0010-310										
Star Handpiece Bundle (Write-in PN and price)											
Quick Disconnect w/Flow Control, 1/4" (for water)	0016-454-CS										Included with 4485SI
Quick Disconnect w/Flow Control, 3/8" (for air)	0016-456-CS										Included with 4485SI
Swivel Tray	1110-468										
Less sidebox telescoping assistant's instrumentation**	1400-000-CS										4485SI only
Chrome Accents	9000-998-CS										
Additional Umbilical, Per Foot (max of 9 feet)	1115-999-CS										
Routing Valve Kit w/Water shutoff & Regulator, Factory inst.	1115-315-CS										Included with 4485SI
Routing Valve Kit for Cuspidor, Factory Installed	1115-316-CS										

* Ancillaries utilize an existing handpiece position - some integrated options may require purchase of Ancillary Shelf PN 0010-310

** Sidebox mounted telescoping assistant's instrumentation included standard with 4485SI

STEP 8: ASSISTANT'S INSTRUMENTATION AND CUSPIDORS - Must select 1 option for order to be valid

No Assistant's Instrumentation or Cuspidor*	N/A
Rear Pivot Instrumentation with Integrated Touchpad	1096
Rear Pivot Instrumentation with Deluxe Cuspidor**	1590
Deluxe Cuspidor only, Rear Pivot**	1580
Basic Cuspidor, Sidebox Mounted (4485SI only)	1679
Deluxe Cuspidor, Sidebox Mounted (4485SI only)	1577

**Not available with 4485SI

STEP 9: ASSISTANT'S INSTRUMENTATION AND CUSPIDOR OPTIONS

Quick Disconnect w/Flow Control, 1/4" (For 1095 instrumentation)	0016-455-CS	Specify One	<input type="checkbox"/> Air	<input type="checkbox"/> Water
Additional Autoclavable HVE with 5' Tubing & Dry Oral Cup	1400-113-CS			
Hands-Free HVE Holder with HVE and 5' Tubing (Select PN) Colors (holder/tubing): GRY (white/grey), STG (white/sterling), BLK (black/black), SHW (black/shadow)		<input type="checkbox"/> 1400-279-STG <input type="checkbox"/> 1400-279-BLK <input type="checkbox"/> 1400-279-SHW <input type="checkbox"/> 1400-279-GRY		
Self-Contained Clean Water System w/Quick Disconnect, 1.5L bottle (For 1095 Instrumentation only)	1700-319			
Tray Assembly (For 1095 Instrumentation only)	1110-680			
Deluxe Cuspidor Unit-Mounted Activation (For 1577, 1580, 1590)	1100-480-CS			
Chair Control Touchpad - Mounted to 1095 Instrumentation	1700-088-CS			

STEP
10**STEP 10: OPERATORY LIGHT - Select Light Head Mounting Option - 1 light and mounting option required**

Mount Type	9078 LED				9073 LED			
	Part #	Retail Price	Qty	Total Price	Part #	Retail Price	Qty	Total Price
Chair Mount	9078				9073			
Ceiling Mount	9078-CL				9073-CL			
Wall Mount	9079W				9074W			
Cabinet Mount	9079C				9074C			
Track Mount	9078T				9073T			

STEP 11: OPERATORY LIGHT CONFIGURATION AND OPTIONS

Configure based on the light-mounting option selected above

STEP
11**Chair Mounted Light Options (Pivot, Sidebox, and Fixed)**

Description	Part #	Retail Price	Qty	Total Price
20" Rigid Arm in Lieu of Standard 32" Arm	1110-227-CS			
Light Posts (Select 1)				
Sidebox, 44" Post	1110-900			
Sidebox w/Monitor Cord slot, 44"	1110-170			
Pivot without chair adapter (adapter included with unit purchase)	1113-654			
Pivot with chair adapter (when chair-mounted unit not ordered)	1113-654-CA			
Monitor Mount Assemblies				
Ring Post Mount Assembly (Must pair with either PN 1110-170 or 1110-045 above) This option allows you to have both a light and monitor mounted to your chair	1300RM			
Top Post Mount Assembly (Must pair with one of the post options below) This option allows a monitor instead of a light on your chair	1300TPM			
Sidebox: Monitor Post, Dual Slot, Fixed Sidebox Mount, 30"	1331-998			
Pivot: Monitor in Lieu of Light Post, Pivot	1331-654			
Pivot: Monitor in Lieu of Light Post, Pivot with Chair Adapter	1331-654-CA			

Ceiling Mounted Light Configuration (Select Mounting Height and Post Length)

Mounting Height (Required)	<input type="checkbox"/> 8'	<input type="checkbox"/> 8'6"	<input type="checkbox"/> 9'	<input type="checkbox"/> 9'6"	<input type="checkbox"/> 10'	<input type="checkbox"/> 10'6"	<input type="checkbox"/> 11'	<input type="checkbox"/> Custom Height, Specify:
Post Length (Required)	<input type="checkbox"/> 12"	<input type="checkbox"/> 18"	<input type="checkbox"/> 24"	<input type="checkbox"/> 30"	<input type="checkbox"/> 36"	<input type="checkbox"/> 42"	<input type="checkbox"/> 48"	<input type="checkbox"/> Custom Length (max 48"), Specify:

Cabinet and Wall Mounted Light Configuration and Options

Description	Part #	Retail Price	Qty	Total Price
Arm Height (Required): <input type="checkbox"/> Standard <input type="checkbox"/> Tall	N/A			
Specify Cabinet or Wall Mounting Height (write-in):	N/A			
Horizontal Top Cabinet Mount in Lieu of Internal Mount	1110-881-CS			
Vertical Exterior Cabinet Mount in Lieu of Internal Mount	1110-886-CS			

Track Mounted Light Configuration and Options

Monitor Mount Assembly (specify monitor post length in next section)	1305T
S Video Cable, Pre-Installed	0014-282
RCA Cable, Pre-Installed	0014-283

Select Track Mount Light Post Lengths

Mounting Height (Required)	<input type="checkbox"/> 8'	<input type="checkbox"/> 8'6"	<input type="checkbox"/> 9'	<input type="checkbox"/> 9'6"	<input type="checkbox"/> 10'	<input type="checkbox"/> 10'6"	<input type="checkbox"/> 11'	<input type="checkbox"/> Custom Height, Specify:
Light Post Length (Required)	<input type="checkbox"/> 7.25"	<input type="checkbox"/> 13.75"	<input type="checkbox"/> 19.75"	<input type="checkbox"/> 25.75"	<input type="checkbox"/> 31.75"	<input type="checkbox"/> 37.75"	<input type="checkbox"/> 43.75"	<input type="checkbox"/> Custom Length (max 43.75"), Specify:
Monitor Post Length (Required if 1305T selected)	<input type="checkbox"/> 20.25"	<input type="checkbox"/> 26.25"	<input type="checkbox"/> 32.25"	<input type="checkbox"/> 38.25"	<input type="checkbox"/> 44.25"	<input type="checkbox"/> 50.25"	<input type="checkbox"/> 56.25"	<input type="checkbox"/> Custom Length (max 56.25"), Specify:

Forest 6400
PackagesForest 6400
Ordering GuidesForest 3900
PackagesForest 3900
Ordering GuideForest Fusion
PackagesForest Fusion
Ordering GuideSurface
TreatmentsChairs
and StoolsLights and
Monitor MountsDelivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and PartsGeneral
InformationEquipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
AccessoriesUtility Quick
Reference GuideUtility
Warranty

Sales Reps

Notes

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

3900 CHAIR MOUNT PACKAGE

Forest

DENTALEZ



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

3900 Pivot Chair Mount Packages

3900 Chair Package, Standard IC Delivery, Pivot

Base Package



Note: Package shown in shadow paint with raven wing ultraleather.

Base Package As Shown

Description	Part Number
3900 Patient Chair	3900
Delivery System w/Standard IC, Pivot	4486PI
Assistants Instrumentation, Rear Pivot	1095
LED Operatory Light, 9078	9078
Light Post, Pivot	1113-654
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Equipment dimensions, see page 109

Shipping weight estimate:

Delivery unit & chair	450 lb / 204 kg
Base package	580 lb / 264 kg
As shown package	630 lb / 286 kg



Popular Options



Dr's Deluxe Stool,
Ultraleather Upholstery

Part number: 6159C



E-Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



LED Operatory Light, 9073

Part number: 9073



Assistants Instrumentation
rear Pivot with touchpad

Part number: 1096

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Popular Options

Description	Part Number
-------------	-------------

Dr. Deluxe Stool,
Ultraleather Upholstery

6159C

Assistants Deluxe Stool No Back,
Ultraleather Upholstery

6160C

E-Torque Electric HP System

STAR-ETORQUE

6-PIN Handpiece Illumination w/tubing

1700-019-CS

Quick Disconnect for Water

0016-454-CS

Quick Disconnect for Air

0016-456-CS

Routing Valve Kit

1115-315-CS

Additional HVE w/Dry Oral Cup

1400-113-CS

Soft Touch Back Massage

3900-100

LED Operatory Light, 9073

9073

Assistants Instrumentation,
Rear Pivot with Touchpad

1096

3900 Pivot Chair Mount Packages

3900 Chair Package, Euro Pivot Delivery

Base Package



Popular Options



Dr's Deluxe Stool,
Ultraleather Upholstery

Part number: 6159C



E- Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



LED Operatory Light, 9073

Part number: 9073



Assistants Instrumentation
rear Pivot with touchpad

Part number: 1096

Base Package As Shown

Description	Part Number
3900 Patient Chair	3900
Delivery System w/Euro, Pivot	4988PI
Assistants Instrumentation, Rear Pivot	1095
LED Operatory Light, 9078	9078
Light Post, Pivot	1113-654
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Equipment dimensions, see page 109

Shipping weight estimate:

Delivery unit & chair 445 lb / 203 kg

Base package 533 lb / 264 kg

As shown package 630 lb / 286 kg

Popular Options

Description	Part Number
Dr. Deluxe Stool, Ultraleather Upholstery	6159C
Assistants Deluxe Stool No Back, Ultraleather Upholstery	6160C
E-Torque Electric HP System	STAR-ETORQUE
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Quick Disconnect for Water	0016-454-CS
Quick Disconnect for Air	0016-456-CS
Routing Valve Kit	1115-315-CS
Additional HVE w/Dry Oral Cup	1400-113-CS
Soft Touch Back Massage	3900-100
LED Operatory Light, 9073	9073
Assistants Instrumentation, Rear Pivot with Touchpad	1096



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

3900 Fixed Chair Mount Packages

3900 Chair Package, Standard IC Delivery, Sidebox

Base Package



Popular Options



Dr's Deluxe Stool,
Ultraleather Upholstery

Part number: 6159C



E- Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



Swivel Tray

Part number: 1110-468



LED Operatory Light, 9073

Part number: 9073

Base Package As Shown

Description	Part Number
3900 Patient Chair	3900
Delivery System, Sidebox w/Standard IC	4485SI
Routing Valve	Included
Quick Disconnect for Water	Included
Quick Disconnect for Air	Included
Deluxe Cuspidor	1577
LED Operatory Light, 9078	9078
Light Post, Sidebox	1110-900
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Equipment dimensions, see page 110

Shipping weight estimate:

Delivery unit & chair 440 lb / 202 kg

Base package 556 lb / 253 kg

As shown package 561 lb / 255 kg



Popular Options

Description	Part Number
Dr. Deluxe Stool, Ultraleather Upholstery	6159C
Assistants Deluxe Stool No Back, Ultraleather Upholstery	6160C
E-Torque Electric HP System	STAR-ETORQUE
Satelec SP	SATELEC-S
Satelec SP, LED lit	SATELEC-S-LED
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Swivel Tray	1110-468
Additional Handpiece Position, Wet	1100-189-CS
Deluxe Cuspidor Remote Activation	1100-480-CS
LED Operatory Light, 9073	9073
Soft Touch Back Massage	3900-100

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

3900 Fixed Chair Mount Packages

3900 Chair Package, Euro Sidebox Delivery

Base Package



Popular Options



Dr's Deluxe Stool,
Ultraleather Upholstery

Part number: 6159C



Assistants Deluxe Stool,
No Back,
Ultraleather Upholstery

Part number: 6160C



E-Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



LED Operatory Light, 9073

Part number: 9073

Base Package As Shown

Description	Part Number
3900 Patient Chair	3900
Delivery System, Sidebox w/Euro	4987SI
Routing Valve	Included
Quick Disconnect for Water	Included
Quick Disconnect for Air	Included
Deluxe Cuspidor	1577
LED Operatory Light, 9078	9078
Light Post, Sidebox	1110-900
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Equipment dimensions, see page 110

Shipping weight estimate:

Delivery unit & chair 445 lb / 203 kg

Base package 551 lb / 251 kg

As shown package 551 lb / 251 kg



Popular Options

Description	Part Number
Dr. Deluxe Stool, Ultraleather Upholstery	6159C
Assistants Deluxe Stool No Back, Ultraleather Upholstery	6160C
E-Torque Electric HP System	STAR-ETORQUE
Satelec SP	SATELEC-S
Satelec SP, LED lit	SATELEC-S-LED
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Additional Handpiece Position, Wet	1100-189-CS
Deluxe Cuspidor Remote Activation	1100-480-CS
LED Operatory Light, 9073	9073
Soft Touch Back Massage	3900-100

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



SURFACE FINISH OPTIONS

ULTRALEATHER® 12 CLASSIC COLORS



Autumn Leaf




Baltic



Boysenberry



Charcoal



Diplomat Blue



Doe



Midnight



Oz



Pewter




Dove Grey



Papyrus



Parrot



Raven Wing



St. John



Stone

ULTRALEATHER® 5 PEARLIZED COLORS




Pharaoh




Tin Man

POWDER COAT PAINT 9 COLORS



Forest Grey




Shadow



Arctic Silver*



Bronze*




Charcoal*

CARBON FIBER 3 PATTERNS



Greywood



Laser Leather



Lemans

* Designer
Friendly Colors

DESIGN YOUR 3900 PATIENT CHAIR

**STEP
1**

STEP 1: SELECT YOUR CHAIR - Select 1 chair and 1 powder coat paint color

Description	Part #	Powder Coat Paint Color
<input type="checkbox"/> 3900 Patient Chair	3900	<input type="checkbox"/> Forest Grey* <input type="checkbox"/> Shadow*
<input type="checkbox"/> 3900 Patient Chair, Designer Friendly	3900C	<input type="checkbox"/> Arctic Silver <input type="checkbox"/> Bronze <input type="checkbox"/> *White <input type="checkbox"/> Charcoal <input type="checkbox"/> Copper <input type="checkbox"/> Midnight <input type="checkbox"/> Pewter
*White, Shadow, and Forest Grey chairs come standard with the baseplate painted Shadow, all other colors come standard with the baseplate painted the same color as the chair.		

**STEP
2**

STEP 2: SELECT YOUR UPHOLSTERY - Select 1 option and 1 upholstery color

Description	Part #	Upholstery Color (used as inner color if dual color upholstery is selected)
<input type="checkbox"/> Wide Back, Plush	3907-PL-WB	Standard Ultraleather® Colors
<input type="checkbox"/> Wide Back, Seamless	3907-AS-WB	<input type="checkbox"/> Autumn Leaf <input type="checkbox"/> Baltic <input type="checkbox"/> Boysenberry <input type="checkbox"/> Charcoal
<input type="checkbox"/> Wide Back, Memory Comfort	3907-MC-WB	<input type="checkbox"/> Doe <input type="checkbox"/> Diplomat Blue <input type="checkbox"/> Dove Grey <input type="checkbox"/> Papyrus
<input type="checkbox"/> Narrow Back, Plush	3907-PL-NB	<input type="checkbox"/> Parrot <input type="checkbox"/> Raven Wing
<input type="checkbox"/> Narrow Back, Seamless	3907-AS-NB	<input type="checkbox"/> St. John <input type="checkbox"/> Stone
<input type="checkbox"/> Narrow Back with Slings, Plush	3907-PL-NS	<input type="checkbox"/> Midnight* <input type="checkbox"/> Pewter* <input type="checkbox"/> Oz* <input type="checkbox"/> Tin Man* <input type="checkbox"/> Pharaoh* *Pearlized
<input type="checkbox"/> Narrow Back with Slings, Seamless	3907-AS-NS	<input type="checkbox"/> Special color, select PN 3907-999-SP from Step 3 and write color choice here: _____

**STEP
3**

STEP 3: ADDITIONAL UPHOLSTERY OPTIONS - Select at least 1 of the below options

<input type="checkbox"/> No Additional Options	N/A	
<input type="checkbox"/> Special Color Ultraleather Upholstery	3907-999-SP	Color:
<input type="checkbox"/> Dual Color, Plush - specify outer color (Not available with seamless or memory comfort)	3907-999-DC	Outer Color: (inner color selected in Step 2 above)
<input type="checkbox"/> Dual Color for Memory Comfort specify outer color	3907-MC-DC	Outer Color: (inner color selected in Step 2 above)
<input type="checkbox"/> Contrast/Deviated Stitching, Plush (Not available with seamless or memory comfort)	3907-999-DS	Thread Color:
<input type="checkbox"/> Neck Pillow	3907-318	Color:
<input type="checkbox"/> Knee Pillow	3907-317	Color (inner if dual color):
<input type="checkbox"/> Dual Color Upholstery, Knee Pillow	3907-317-DC	Outer Color:
<input type="checkbox"/> Contrast Stitch Upholstery, Knee Pillow	3907-317-DS	Thread Color:
<input type="checkbox"/> Pedo Booster Set	3907-180	
<input type="checkbox"/> 29" Toe Guard in Lieu of 21"	3907-028-CS	

**STEP
4**

STEP 4: MASSAGE OPTION - Select at least 1 of the below options

<input type="checkbox"/> No Massage Options	N/A	
<input type="checkbox"/> Soft Touch Massage*	3900-100	Not available with Memory Comfort Upholstery

*If a chair-mounted delivery unit is being ordered then a utility center housing is included with purchase. Otherwise we recommend adding a utility center for the massage air compressor (Part numbers below).

**STEP
5**

STEP 5: UTILITY CENTER SELECTION - Select at least 1 of the below options

<input type="checkbox"/> No Additional Options	N/A	Select this option if you are buying a chair-mount delivery unit and no massager
Contoured Utility Center Housing and Cover		
<input type="checkbox"/> Forest Grey	3900-101-G	
<input type="checkbox"/> Shadow	3900-101-SH	
<input type="checkbox"/> White	3900-101	
<input type="checkbox"/> Designer Friendly®	3900-101-C	

**STEP
6**

STEP 6: CHAIR CONTROL LOCATION - Select at least 1 chair control

Maximum of 1 control can be mounted on the delivery unit, maximum of 2 on the chair (left & right), 4 total maximum number of controls per operatory

<input type="checkbox"/> First Chair Control (Required) First control included with chair purchase	3914-998-CS	<input type="checkbox"/> Delivery Unit <input type="checkbox"/> Footswitch <input type="checkbox"/> Chair, right side <input type="checkbox"/> Chair, left side
<input type="checkbox"/> Second Chair Control (Optional)	3914-999-CS	<input type="checkbox"/> Delivery Unit <input type="checkbox"/> Footswitch <input type="checkbox"/> Chair, right side <input type="checkbox"/> Chair, left side

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

PATIENT CHAIR

STEP

1


Chair paint color



STEP

2

Upholstery color (inner color if dual color selected)



STEP

3

Outer upholstery color (if dual color selected, PN 3907-99-DC)
Contrast stitch color (if color stitch selected, PN 3907-999-DS) -----



CHAIR IMAGES



- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

DESIGN YOUR STANDARD IC DELIVERY SYSTEM




SURFACE FINISH OPTIONS

POWDER COAT PAINT 9 COLORS




Arctic Silver



Forest Grey



Charcoal



Copper



Bronze




Midnight



Pewter



Shadow




White


CARBON FIBER 3 PATTERNS



Greywood



Laser Leather



Lemans

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

STANDARD IC DELIVERY UNIT

STEP

8


Powder coat
paint color



STEP

8


Powder coat paint with upgrade
option to chrome/polished finish
(PN 9000-998-CS, Step 9)




STEP

8

Plastic finish



UNIT IMAGES



Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty


Sales Reps

DESIGN YOUR EURO DELIVERY SYSTEM




SURFACE FINISH OPTIONS

POWDER COAT PAINT 9 COLORS



Forest Grey



Shadow




Arctic Silver*



Bronze*



Charcoal*



Copper*



Midnight*



Pewter*




White*

* Designer
Friendly Colors

CARBON FIBER 3 PATTERNS



Greywood



Laser
Leather



Lemans

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

EURO DELIVERY UNIT

STEP

8


Powder coat
paint color



STEP

8




Powder coat paint with upgrade
option to chrome/polished finish
(PN 9000-998-CS, Step 9)



STEP

8

Plastic finish



UNIT IMAGES



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

DESIGN YOUR LED OPERATORY LIGHT




SURFACE FINISH OPTIONS


POWDER COAT PAINT 9 COLORS



Forest Grey




Shadow




Arctic Silver*



Bronze*




Charcoal*



Copper*



Midnight*



Pewter*



White*

* Designer
Friendly Colors

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty


Sales Reps

OPERATORY LIGHT

STEP

12

Powder coat paint



LIGHT IMAGES



Forest 6400
Ordering Guide

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

THE FUSION PACKAGE

Forest

DENTALEZ



LUXURY MEETS AFFORDABILITY

The Fusion Package integrates thoughtful design with premium materials for years of beauty at an exceptional value.

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Fusion Packages

Fusion, Pivot

Base Package



Popular Options



Dr's Economy Stool,
Naugahyde Upholstery

Part number: 6156V



Assistants Economy Stool,
Naugahyde Upholstery

Part number: 6157V



E-Torque Electric
Handpiece System

Part number:
STAR-ETORQUE



Assistants Instrumentation
Rear Pivot with Touchpad

Part number: 1296

Base Package As Shown

Description	Part Number
3800 Patient Chair, Pretubed	3800PT
Delivery System, w/Fusion, Pivot	4296PI
Assistants Instrumentation, Rear Pivot	1295
LED Operatory Light, 9078	9078
Light Post, Pivot	1113-654
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Equipment dimensions, see page 109

Shipping weight estimate:

Delivery Unit & Chair 445 lb / 203 kg
Base package 445 lb / 203 kg
As shown package 585 lb / 266 kg



Popular Options

Description	Part Number
Dr. Economy Stool, Naugahyde Upholstery	6156V
Assistants Economy Stool, Naugahyde Upholstery	6157V
E-Torque Electric HP System	STAR-ETORQUE
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Quick Disconnect for Water	0016-454-CS
Quick Disconnect for Air	0016-456-CS
Routing Valve Kit	1115-315-CS
Additional HVE w/Dry Oral Cup	1400-113-CS
Ultraleather Upgrade	3807-999-UL
Assistants Instrumentation, Rear Pivot with Touchpad	1296

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Fusion Packages

Fusion, Sidebox

Base Package



Popular Options



Dr's Economy Stool,
Naugahyde Upholstery

Part number: 6156V



Assistants Economy Stool,
Naugahyde Upholstery

Part number: 6157V



E-Torque Electric
Handpiece System

Part number:
STAR-ETORQUE

Base Package As Shown

Description	Part Number
3800 Patient Chair, Pretubed	3800PT
Delivery System, Sidebox w/Fusion	4295SI
Routing Valve	Included
Quick Disconnect for Water	Included
Quick Disconnect for Air	Included
Deluxe Cuspidor	1577
LED Operatory Light, 9078	9078
Light Post, Sidebox	1110-900
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Equipment dimensions, see page 110

Shipping weight estimate:

Delivery unit & chair 440 lb / 202 kg
Base package 440 lb / 202 kg
As shown package 561 lb / 255 kg



Popular Options

Description	Part Number
Dr. Economy Stool, Naugahyde Upholstery	6156V
Assistants Economy Stool, Naugahyde Upholstery	6157V
E-Torque Electric HP System	STAR-ETORQUE
Satelec SP	SATELEC-S
Satelec SP, LED lit	SATELEC-S-LED
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Swivel Tray	1110-468
Additional Handpiece Position, Wet	1100-189-CS
Deluxe Cuspidor Remote Activation	1100-480-CS

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Fusion Packages

Fusion, 2186BC Rear Delivery

Base Package



Popular Options



Dr's Economy Stool,
Naugahyde Upholstery

Part number: 6156V



Assistants Economy Stool,
Naugahyde Upholstery

Part number: 6157V

Base Package As Shown

Description	Part Number
3800 Patient Chair	3800
Delivery System, Beneath Counter Duo	2186BC
LED Operatory Light, 9078	9078
Light Post, Pivot w/Chair Adapter	1113-654-CA
Chair Control (select from unit-mounted, chair-mounted, or foot control)	3914-998-CS

Popular Options

Description	Part Number
Dr. Economy Stool, Naugahyde Upholstery	6156V
Assistants Economy Stool, Naugahyde Upholstery	6157V
6-PIN Handpiece Illumination w/tubing	1700-019-CS
Quick Disconnect for Water	0016-454-CS
Quick Disconnect for Air	0016-456-CS
Routing Valve Kit	1115-315-CS
Additional HVE w/Dry Oral Cup	1400-113-CS
Ultraleather Upgrade	3807-999-UL
LED Operatory Light, Ceiling Mounted	9078-CL
LED Operatory Light, Wall Mounted	9079W
Deluxe Utility Center Housing and Cover	1115-007-G

Shipping weight estimate:

Delivery Unit & Chair 445 lb / 203 kg
Base package 445 lb / 203 kg
As shown package 585 lb / 266 kg




- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages**
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series

- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors


- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps



FUSION PACKAGE ORDERING GUIDE



NAUGAHYDE® UPHOLSTERY OPTIONS



Black



Cocoa



Deep
Sapphire



Saddle



Sand



Taupe

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

DESIGN YOUR FUSION OPERATORY CHAIR

**STEP
1**

STEP 1: SELECT YOUR CHAIR - Select 1 chair

Description	Part #	Qty.	Comments, color
<input type="checkbox"/> 3800 Patient Chair, non-tubed	3800		Forest Grey w/Shadow baseplate only
<input type="checkbox"/> 3800 Patient Chair, pre-tubed	3800PT		Forest Grey w/Shadow baseplate only

**STEP
2**

STEP 2: SELECT YOUR UPHOLSTERY - Select 1 option and 1 upholstery color

<input type="checkbox"/> Standard Naugahyde® Upholstery, Seamless (Included)	N/A	<input type="checkbox"/> Black <input type="checkbox"/> Cocoa <input type="checkbox"/> Deep Sapphire <input type="checkbox"/> Saddle <input type="checkbox"/> Sand <input type="checkbox"/> Taupe
<input type="checkbox"/> Upgrade to standard Ultraleather® color, Seamless	3807-999-UL	Standard Ultraleather Colors <input type="checkbox"/> Autumn Leaf <input type="checkbox"/> Baltic <input type="checkbox"/> Boysenberry <input type="checkbox"/> Charcoal <input type="checkbox"/> Diplomat Blue <input type="checkbox"/> Doe <input type="checkbox"/> Dove Grey <input type="checkbox"/> Papyrus <input type="checkbox"/> Parrot <input type="checkbox"/> Raven Wing <input type="checkbox"/> St. John <input type="checkbox"/> Stone
<input type="checkbox"/> Upgrade, special color Ultraleather, Seamless	3807-999-UL-SP	Standard Pearlized Ultraleather Colors <input type="checkbox"/> Midnight <input type="checkbox"/> Pewter <input type="checkbox"/> Oz <input type="checkbox"/> Pharaoh <input type="checkbox"/> Tin Man

**STEP
3**

STEP 3: ADDITIONAL UPHOLSTERY OPTIONS - Select at least 1 of the below options

Description	Part #	Qty.	Comments, color
<input type="checkbox"/> No Additional Options	N/A		
<input type="checkbox"/> Neck Pillow (only available in Ultraleather)	3907-318		Color:
<input type="checkbox"/> Knee Pillow (only available in Ultraleather)	3907-317		Color (inner if dual color):
<input type="checkbox"/> Dual Color Upholstery, Knee Pillow	3907-317-DC		Outer Color:
<input type="checkbox"/> Contrast Stitch Upholstery, Knee Pillow	3907-317-DS		Thread Color:
<input type="checkbox"/> Pedo Booster Seat (only available in Ultraleather)	3907-180		Color:

STEP 4: UTILITY CENTER OPTIONS

Description	Part #	Qty.	Comments, color
<input type="checkbox"/> No Utility Center*	N/A		Select this option if you are buying a chair-mounted delivery unit
<input type="checkbox"/> Deluxe Utility Center Housing & Cover, Rear Umbilical Access, Forest Grey	1115-007-G		
<input type="checkbox"/> Deluxe Utility Center Housing & Cover, Right Umbilical Access, Forest Grey	1115-009-G		

*If you are ordering a chair-mounted delivery unit then select "No Utility Center" because one will be included with your delivery unit

**STEP
5**

STEP 5: CHAIR CONTROL LOCATION - Select at least 1 chair control

Maximum of 1 control can be mounted to the delivery unit, maximum of 2 on the chair (left and right),
4 total maximum number of controls per operatory

<input type="checkbox"/> First Chair Control (Required) First control included with chair purchase	3914-998-CS	<input type="checkbox"/> Delivery Unit <input type="checkbox"/> Footswitch <input type="checkbox"/> Chair, right side <input type="checkbox"/> Chair, left side
<input type="checkbox"/> Second Chair Control (Optional)	3914-999-CS	<input type="checkbox"/> Delivery Unit <input type="checkbox"/> Footswitch <input type="checkbox"/> Chair, right side <input type="checkbox"/> Chair, left side

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

DESIGN YOUR OPERATORY LIGHT SELECTION

**STEP
10**

STEP 10: OPERATORY LIGHT

Select Light Head Mounting Option - 1 light and mounting option required

	9083 Halogen	9078 LED		
	Part #	Qty	Part #	Qty
Lights with Fusion packages come in Forest Grey Only	<input type="checkbox"/> 9083		<input type="checkbox"/> 9078	
Chair Mount	<input type="checkbox"/> 9083-CL		<input type="checkbox"/> 9078-CL	
Ceiling Mount	<input type="checkbox"/> 9093W		<input type="checkbox"/> 9079W	
Wall Mount	<input type="checkbox"/> 9093C		<input type="checkbox"/> 9079C	
Cabinet Mount	<input type="checkbox"/> 9083T		<input type="checkbox"/> 9078T	
Track Mount				

**STEP
11**

STEP 11: OPERATORY LIGHT CONFIGURATION AND OPTIONS

Chair Mounted Light Options (Pivot, Sidebox, and Fixed)

Description	Part #	Qty
<input type="checkbox"/> 20" Rigid Arm in Lieu of Standard 32" Arm	1110-227-CS	
Light Posts (Select 1)		
<input type="checkbox"/> Sidebox, 44" Post	1110-900	
<input type="checkbox"/> Sidebox w/Monitor Cord slot, 44"	1110-170	
<input type="checkbox"/> Pivot without chair adapter (adapter included with unit purchase)	1113-654	
<input type="checkbox"/> Pivot with 3800/3900 chair adapter (when chair-mounted unit not ordered)	1113-654-CA	
Monitor Mount Assemblies		
<input type="checkbox"/> Ring Post Mount Assembly (Must pair with either PN 1110-170 or 1110-045 above) This option allows you to have both a light and monitor mounted to your chair	1300RM	
<input type="checkbox"/> Top Post Mount Assembly (Must pair with one of the post options below) This option allows a monitor instead of a light on your chair	1300TPM	
<input type="checkbox"/> Sidebox: Monitor Post, Dual Slot, Fixed Sidebox Mount, 30"	1331-998	
<input type="checkbox"/> Pivot: Monitor in Lieu of Light Post, Pivot	1331-654	
<input type="checkbox"/> Pivot: Monitor in Lieu of Light Post, Pivot with Chair Adapter	1331-654-CA	

Ceiling Mounted Light Configuration (Select Mounting Height and Post Length)

Mounting Height (Required) 8' 8'6" 9' 9'6" 10' 10'6" 11' Custom Height, Specify:

Post Length (Required) 12" 18" 24" 30" 36" 42" 48" Custom Length (max 48"), Specify:

Cabinet and Wall Mounted Light Configuration and Options

Description	Part #	Qty
Arm Height (Required): <input type="checkbox"/> Standard <input type="checkbox"/> Tall	N/A	
Specify Cabinet or Wall Mounting Height (write-in):	N/A	
<input type="checkbox"/> Horizontal Top Cabinet Mount in Lieu of Internal Mount	1110-881-CS	
<input type="checkbox"/> Vertical Exterior Cabinet Mount in Lieu of Internal Mount	1110-886-CS	

Track Mounted Light Configuration and Options

<input type="checkbox"/> Monitor Mount Assembly (specify monitor post length in next section)	1305T
<input type="checkbox"/> S Video Cable, Pre-Installed	0014-282
<input type="checkbox"/> RCA Cable, Pre-Installed	0014-283

Select Track Mount Light Post Lengths

Mounting Height (Required) 8' 8'6" 9' 9'6" 10' 10'6" 11' Custom Height, Specify:

Light Post Length (Required) 8.25" 14.25" 20.25" 26.25" 32.25" 38.25" 44.25" Custom Length (max 44.25"), Specify:

Monitor Post Length (Required if 1305T selected) 20.25" 26.25" 32.25" 38.25" 44.25" 50.25" 56.25" Custom Length (max 56.25"), Specify:

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface Treatments

Chairs and Stools

Lights and Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

Surface Treatments

Forest is the only dental equipment with a WOW factor. It really is a kick watching a dentist react to Forest equipment when they first see it. A common reaction is “wow, that is beautiful.”

Carbon Fiber

Control Heads

- Standard IC - Chassis
- Euro - Lid

Sidebox - Cover

Basic Cuspidor - Cover

Assistant's Instrumentation

- 1095 - Arm Cover

Powder Coat Paint

Dentists and interior designers love being able to put something other than white into their beautiful patient friendly decors.

Control Heads

- Standard IC - Lid, Brake Handles and Holder Bar
- Euro - Chassis, Whip Arms, Brake Handles and Swivel Tray
- 7000 Series Rear Mount, Side Delivery and Mobile Cart
- 7020 PRO Series Rear Mount, Side Delivery and Mobile Cart

Sidebox - Housing

Basic Cuspidor - Housing and Holder Bar

Assistant's Instrumentation

- 1095 - Arm
- 1096 - Arm and Integrated Holder

Ultraleather® (standard)

Once you see a Forest chair, delivery system and light package with color coordinated paint and upholstery in a well designed treatment room, you wonder why it's taken so long for a manufacturer to offer dental equipment with a full range of surface treatments.

3900 Chair

- Plush Upholstery
- Seamless Upholstery
- Knee Pillow, Neck Pillow and Pedo Booster Seat
- Deluxe and Economy Stools

Solid Surface (standard)

Control Heads

- 7010 Series Rear Mount Side Delivery and Mobile Cart
- 7020 PRO Series Rear Mount Side Delivery and Mobile Cart
- 3486 Rear Mount Side Delivery and Mobile Cart
- 2184BC/2186BC Rear and Side Delivery

Laminate (optional)

Control Heads

- 3486 Rear Mount Side Delivery
- 2184BC/2186BC Rear and Side Delivery
- 7010 Series Rear Mount Side Delivery and Mobile Cart

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments

Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series

Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty




Compressors
Dry Vacuums

Tanks
Wet Vacuums
Utility Room Accessories

Utility Quick Reference Guide
Utility Warranty

Sales Reps

POWDER COAT PAINT (9)



CARBON FIBER (3) FORMICA® (4)



SOLID SURFACES (10)



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories


Utility Quick
Reference Guide

Utility
Warranty


Sales Reps

NAUGAHYDE®* (6)

*For The Fusion Packages only. Pages 22-29




BLACK
NaugaHyde®




COCOA
NaugaHyde®



DEEP SAPPHIRE
NaugaHyde®



SADDLE
NaugaHyde®



SAND
NaugaHyde®



TAUPE
NaugaHyde®

ULTRALEATHER®* (12 CLASSIC) (5 PEARLIZED)

*Pearlized




AUTUMN LEAF
Ultraleather®



BALTIC
Ultraleather®



BOYSENBERRY
Ultraleather®



CHARCOAL
Ultraleather®




DIPLOMAT BLUE
Ultraleather®



DOE
Ultraleather®




DOVE GREY
Ultraleather®




PAPYRUS
Ultraleather®



PARROT
Ultraleather®



RAVEN WING
Ultraleather®



ST. JOHN
Ultraleather®



STONE
Ultraleather®




MIDNIGHT*
Ultraleather®



OZ*
Ultraleather®



PEWTER*
Ultraleather®



PHARAOH*
Ultraleather®



TIN MAN*
Ultraleather®

CUSTOMIZE YOUR TREATMENT ROOMS WITH DENTALEZ

Powder coat paint your drawer and door pulls, grab handles, trash receptacles, or other metal surfaces in your treatment room to match your delivery system paint color

Source custom powder coat paint colors to match your provided sample

Special order any Ultraleather® upholstery color shown on the Ultrafabrics® website (from original or pearlized palette)

Forest

CHAIRS AND STOOLS



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

**Chairs
and Stools**

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks


Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



* Not a standard color

* Not a standard color

3800 CHAIR

- Material: hydraulic cantilever mechanism
Benefit: smooth and steady heavy lifting
- Material: cast aluminum chair back, base, armrest and headrest
Benefit: strong, impervious to fading or cracking
- Material: flip down armrests
Benefit: easy patient entry and exit
- Material: optional soft-touch back massage with integrated compressor
Benefit: saves central compressor capacity for delivery systems

3900 CHAIR

- Material: cast aluminum chair back, chair base and utility center base
Benefit: strong, impervious to fading or cracking
- Material: pressure-formed plastic chair and utility center covers
Benefit: warp, impact and stain resistant
- Material: aluminum hydraulic manifold and speed adjustment screws
Benefit: independent speed adjustment of base and back both up and down
- Material: optional soft-touch back massage with integrated compressor
Benefit: saves central compressor capacity for delivery systems
- Material: Ultraleather® upholstery standard
Benefit: withstands heavy use and cleaning, with glove-soft leather feel
- Material: Forest Grey or Shadow powder coat standard on 3900 chair
Benefit: seven optional Designer Friendly® colors
- Material: complements treatment room interior design color palette
Benefit:



A Beautiful Blend of Ergonomics & Intelligence

6400 Chair

Forest
DENTALEZ

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools**
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

Chairs and Stools - Chairs

6400 Hydraulic Patient Chair

Model 6400

Part Number 6400



Plush

Seamless

Standard Features

Hydraulic Chair 6400

Ultraleather® upholstery with wide back for elbow support
Plastic toe guard 21"
Double articulating pinch-style headrest
Right/left patient access
Swing-down locking armrests
Hydraulic cantilever mechanism
Emergency shut-off
30 degree base rotation w/lock
Multifunction touchpad
w/4 programmable preset positions
Seat height: minimum 17" maximum 36"
Smart technology monitors chair performance and operatory efficiency

Factory Options

Part No

CHAIR (Forest Grey or Shadow) 6400

CHAIR Designer Friendly® 6400C

Delivery Unit Control 6414-046

Footswitch 3625-858

! UPHOLSTERY SELECTIONS:

Seamed, Memory Comfort 3625-961

Seamless, Memory Comfort 3625-823A

Seamed, Memory Comfort with Heat & Massage Option 3625-962

Seamless, Memory Comfort with Heat & Massage Option 3625-966

Dual Upholstery Color 3907-999-DC

Secondary Color:

Contrast Stitch 3907-999-DS

Thread Color: Choose your thread color from available upholstery color along with Gold, White or Red

! CHAIR OPTIONS:

CareGlide 3625-970

Duplex with USB 3658-495

Surgical Integrated Dual Rails: must select this along with either the IV Armrest or Bar Style Armrest 3680-444

IV Armrest with clamp 3625-964

Bar Style Armrest with clamp 3625-965

Second IV Armrest with clamp 3951-306

IV Bag Stand 3625-812

22" Black Strap 3802-512

54" Black Strap 3802-513

80" Black Strap 3802-514

! LIGHT POSTS (for unit mount lights/post mount monitor assemblies see pages 62)

Fixed (chair adaptor sold separately)

44" light post without slot, chair mount 1110-794

44" light post with slot for monitor cord 1110-050

Pivot (includes pivot cast chair adaptor, bushing, bearing kit - mounts on Forest chairs only) 1113-654-CA

Pivot monitor post, (includes pivot cast chair adapter, bushing, bearing kit - mounts on Forest chairs only) (monitor mount assembly not included) 1331-654-CA

Monitor mount assembly, top post mount (see page 65) 1300TPM

Monitor mount assembly, ring mount (see page 65) 1300RM

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options

and Parts

General

Information

Equipment

Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room

Accessories

Utility Quick

Reference Guide

Utility

Warranty

Sales Reps

! Designer Friendly® Surface Treatments

STEP 1 SELECT STANDARD POWDER COAT PAINT (CHAIR):

Forest Grey Shadow

UPGRADE: Chair Base/Back Powder Coat Paint:

Arctic Silver Bronze Charcoal

Copper Midnight Pewter

White

STEP 2 SELECT ULTRALEATHER® UPHOLSTERY COLOR CHOICES

FOR CHAIR AND DELUXE /ECONOMY STOOLS:

Autumn Leaf Baltic Charcoal

Diplomat Blue Doe Dove Grey

Boysenberry Papyrus Parrot

Raven Wing Stone St. John

Pearlized Colors: Midnight Pewter Oz

Pharaoh Tin Man

! OR SELECT Special Color Ultraleather (Call for leadtime)

Chair

3907-999-SP

! LIGHT POSTS (for unit mount lights/post mount monitor assemblies see pages 62)

Fixed (chair adaptor sold separately)

44" light post without slot, chair mount 1110-794

44" light post with slot for monitor cord 1110-050

Pivot (includes pivot cast chair adaptor, bushing, bearing kit - mounts on Forest chairs only)

Pivot monitor post, (includes pivot cast chair adapter, bushing, bearing kit - mounts on Forest chairs only) (monitor mount assembly not included)

Monitor mount assembly, top post mount (see page 65)

Monitor mount assembly, ring mount (see page 65)

Chairs and Stools

3900 Hydraulic Patient Chair



Plush


Model 3900

Part Number

3900



Seamless



Standard Features

Hydraulic Chair 3900

Ultraleather® upholstery with wide/thin back for elbow support
Plastic toe guard 21"
Double articulating pinch-style headrest
Thin back w/full lumbar support
Right/left patient access
Swing-down locking armrests
Hydraulic cantilever mechanism
Emergency shut-off
30 degree base rotation w/lock
Multifunction touchpad
w/4 programmable preset positions
Seat height: minimum 17" maximum 29"

Factory Options

Part No

CHAIR (Forest Grey or Shadow)	3900
CHAIR Designer Friendly®	3900C
CHAIR PRO Designer Friendly Includes: Plush dual color and deviated/contrast stitching upholstery, massage, contoured utility center housing and cover, knee pillow and neck pillow	3900C-PRO
CHAIR PRO Designer Friendly less utility center housing and cover	3900C-PRO-LUC
One chair control included (specify type below)	3914-998-CS*
<input type="radio"/> Footswitch control <input type="radio"/> Touchpad, chair-mounted	
! Additional chair control (limit 2 on chair, limit 1 on unit)	3914-999-CS
<input type="radio"/> Footswitch control <input type="radio"/> Touchpad, chair-mounted	
Chair adapter, fixed (for external umbilical)	1115-195
Contoured utility center housing & cover	3900-101
Contoured utility center housing & cover Designer Friendly	3900-101-C
! UPHOLSTERY SELECTIONS:	
Pedo booster seat (see page 56)	3907-180
Knee pillow (see page 56)	3907-317
Neck pillow (see page 56)	3907-318
Soft-touch back massage option (not available with Memory Comfort)	3900-100
Plastic toe guard 29" in lieu of 21"	3907-028-CS
! LIGHT POSTS (for unit mount lights/post mount monitor assemblies see pages 62)	
Fixed (chair adaptor sold separately) 44" light post without slot, chair mount	1110-794
44" light post with slot for monitor cord	1110-050
Pivot (includes pivot cast chair adaptor, bushing, bearing kit - mounts on Forest chairs only)	1113-654-CA
Pivot monitor post, (includes pivot cast chair adapter, bushing, bearing kit - mounts on Forest chairs only) (monitor mount assembly not included)	1331-654-CA
Monitor mount assembly, top post mount (see page 65)	1300TPM
Monitor mount assembly, ring mount (see page 65)	1300RM

- ! Often overlooked
- Δ Lead time varies

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

3900 Chair Upholstery Options



Memory Comfort



Narrow Back or Narrow Back with Slings

Standard Features	Factory Options	Part No
Memory Comfort Upholstery Designed for all day/ultra-long procedures or greater general patient comfort 3 layers of foam topped with 1.5 inches of memory foam for extra comfort Side-bolsters on seat cushion support legs Memory foam headrest	UPHOLSTERY SELECTIONS: Plush selections: Standard back, wide/thin Narrow back Narrow back with slings Plush additional options: Contrast/deviated stitching Dual color Dual color with contrast/deviated stitching Call customer service for information (1.800.423.3555) Massage option Seamless selections: Standard back, wide/thin Narrow Back Narrow Back with Slings Seamless additional options: Massage option Memory comfort selections: Wide/thin Dual color	3907-PL-WB 3907-PL-NB 3907-PL-NS 3907-999-DS 3907-999-DC 3900-100 3907-AS-WB 3907-AS-NB 3907-AS-NS 3900-100 3907-MC-WB 3907-MC-DC
Narrow Back Perfect for Pediatric or Orthodontic dentistry Unrestricted patient access		
Narrow Back with Slings Same as above, with additional features: Slings provide alternative for additional patient elbow support Slings create a barrier between patient and practitioner		



3900 Chair Equipment dimensions, see page 109

Shipping weight estimate:
335 lb / 152 kg

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools**
- Lights and Monitor Mounts
- Delivery Systems 7000 Series

- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

3800 Hydraulic Patient Chair



Model 3800

Part Number 3800

Standard Features

Hydraulic Chair 3800

Naugahyde® upholstery
w/wide back for elbow support
Plastic toe guard 21"
Double articulating pinch-style headrest
Right/left patient access
Swing-down locking armrests
Hydraulic cantilever mechanism
Aluminum casting chair back, base,
armrest and headrest
Emergency shut-off
30-degree swivel rotation w/lock

Factory Options

Part No

CHAIR (pre-tubed)	3800PT
CHAIR (non-tubed)	3800
One chair control included (specify type below)	3914-998-CS
<input type="radio"/> Footswitch control	
<input type="radio"/> Touchpad, chair-mounted	
! Additional chair control (limit 2 on chair, limit 1 on unit)	3914-999-CS
<input type="radio"/> Footswitch control	
<input type="radio"/> Touchpad, chair-mounted	
Chair adapter, fixed (for external umbilical)	1115-195
Deluxe utility center housing and cover, rear umbilical access	1115-007-G
Deluxe utility center housing and cover, right umbilical access	1115-009-G
! UPHOLSTERY SELECTIONS (only available in ultraleather):	
Pedo booster seat (see page 56)	3907-180
Knee pillow (see page 56)	3907-317
Neck pillow (see page 56)	3907-318
! LIGHT POSTS (for unit mount lights/post mount monitor assemblies see pages 62)	
Fixed (chair adaptor sold separately)	1110-794
44" light post without slot, chair mount	
44" light post with slot for monitor cord	1110-050
Pivot (includes pivot cast chair adaptor, bushing, bearing kit - mounts on Forest chairs only)	1113-654-CA
Pivot monitor post, (includes pivot cast chair adapter, bushing, bearing kit - mounts on Forest chairs only) (monitor mount assembly not included)	1331-654-CA
Monitor mount assembly, top post mount (see page 65)	1300TPM
Monitor mount assembly, ring mount (see page 65)	1300RM

- !** Often overlooked
- Δ** Lead time varies

Designer Friendly® Surface Treatments

STEP 1 SELECT STANDARD POWDER COAT PAINT (CHAIR):

Forest Grey

STEP 2 SELECT NAUGAHYDE® UPHOLSTERY COLOR FOR YOUR CHAIR:

Black Cocoa Deep Sapphire
 Saddle Sand Taupe

UPGRADE: Standard Ultraleather® upholstery (select from

Forest standard options) 3807-999-UL

Δ OR SELECT Special Color upholstery (Call for leadtime)

Ultraleather 3807-999-UL-SP

Naugahyde 3807-999-SP



3800 Chair Equipment dimensions, see page 108
305 lb / 152 kg

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Chairs and Stools Accessories

Chair Accessories

Neck Pillow



3907-318

Provides ergonomic support for patients with neck discomfort and/or during long procedures.

Neck support offers patients a comfortable relaxed experience.

Specify Ultraleather® upholstery color

Knee Pillow



3907-317

Knee support ideal for patients with joint discomfort and/or long procedures.

Zero gravity posture, eliminates pressure on knees.

Specify Ultraleather upholstery color

3900 dual color upholstery option (upcharge) 3907-317-DC

3900 contrast/deviated stitching option (upcharge) 3907-317-DS

Pedo Booster Seat (to accommodate smaller patients)



3907-180

A cushion that takes the lumbar out of the chair, essentially turning it into a bed, ideal for children.

Specify Ultraleather upholstery color

3900/3800 Chair Control, Footswitch

Additional Chair Control, Footswitch



Factory installed 3914-998-CS

Factory installed 3914-999-CS

After Market 3914-010-SH

4 programmable preset positions

3900/3800 Remote Mount Touchpad with 25' Cable



FIXED HORIZONTAL MOUNT

White 1700-079

Forest Grey 1700-079-G

Designer Friendly® (specify color) 1700-079-C

SWING HORIZONTAL MOUNT

White 1700-089-HM

Forest Grey 1700-089-HM-G

Designer Friendly (specify color) 1700-089-HM-C

SWING VERTICAL MOUNT

White 1700-089-VM

Forest Grey 1700-089-VM-G

Designer Friendly (specify color) 1700-089-VM-C

• 4 programmable preset positions

3900 Contoured Utility Center Housing & Cover



Factory option: 3900-101

Forest Grey 3900-101-G

Shadow 3900-101-SH

Designer Friendly (specify color) 3900-101-C

Included w/4486PI, 4485SI, 4987SI and 4988PI chair mounted packages

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options
and PartsGeneral
InformationEquipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
AccessoriesUtility Quick
Reference Guide

Utility

Warranty

Sales Reps

Chairs and Stools - Stools

Economy and Deluxe Stools



Model 6156



Model 6157



Model 6159



Model 6160



Model 6161

Part Number 6156C

Part Number 6157C

Part Number 6159C

Part Number 6160C

Part Number 6161C

Part Number 6156V

Part Number 6157V

Part Number 6159V

Part Number 6160V

Part Number 6161V

Standard Features

Economy and Deluxe Stools

Ultraleather® or Naugahyde® Upholstery
Injection molded foam
5-leg fully polished base
Twin wheel hard floor casters
Accommodates maximum of 300 lbs.

Doctor's Economy 6156C/6156V

Backrest, non-adjusting
Round seat with single lever height adjustment
Pneumatic cylinder height adjustment, 18" - 23"

Assistant's Economy 6157C/6157V

Round seat with single lever height adjustment
Torso support, non-ratcheting, with height adjustment
Pneumatic cylinder height adjustment, 21 1/2" - 28"
Foot ring affixed to base

Doctor's Deluxe 6159C/6159V

Backrest with height adjustment
Contoured seat with triple lever adjustment
Pneumatic cylinder height adjustment, 18 1/2" - 23"

Assistant's Deluxe 6160C/6160V

Contoured seat with single lever adjustment
Torso support with ratcheting and height adjustment
Pneumatic cylinder height adjustment, 21 1/2" - 28"
Foot ring with adjustment

Assistant's Deluxe with Back 6161C/6161V

Backrest with height adjustment
Contoured seat with triple lever adjustment
Torso support with ratcheting and height adjustment
Pneumatic cylinder height adjustment, 21 1/2" - 28"
Foot ring with adjustment

Equipment dimensions, see page 109

Shipping weight estimate:

6156/6159 35 lb / 15 kg

6161 45 lb / 21 kg

6157/6160 40 lb / 18 kg

Nomenclature Key

C = Ultraleather

6157C

V = Naugahyde

6156V

Example

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Designer Friendly® Surface Treatments

STEP 1 SELECT UTRALEATHER® UPHOLSTERY COLOR CHOICES FOR CHAIR AND DELUXE/ECONOMY STOOLS:

- Autumn Leaf Baltic Charcoal Diplomat Blue Doe Dove Grey
- Boysenberry Papyrus Parrot Raven Wing Stone St. John

Pearlized Colors: Midnight Pewter Pharaoh Oz Tin Man

OR SELECT NAUGAHYDE® UPHOLSTERY COLOR:

- Black Cocoa Deep Sapphire Saddle Sand Taupe

OR SELECT Special Color Ultraleather® (Call for leadtime)

- Deluxe/Economy Stool 6007-999-SP

Comfort Stools



Model 6176

Part Number 6176



Model 6187

Part Number 6187

Standard Features

Comfort Stools

Air circulating weave construction
Active tilt® seat (actively moves with user)
Tilt seat lever function
Contoured seat
5-leg base
Twin wheel hard floor casters
Accommodates maximum of 300 lbs.

Doctor's Comfort 6176

Triple lever adjustment - independent back and seat
Backrest with angle & height adjustment
Pneumatic cylinder height adjustment,
19 1/4"– 24 1/2"

Assistant's Comfort 6187

Triple lever adjustment - independent back and seat
Backrest with angle & height adjustment
Fully adjustable ratcheting assistant arm
Foot ring with height adjustment
Pneumatic cylinder height adjustment, 23"– 30 1/2"

Factory Options

Part No

STOOL	MODEL CHOICE	Part No
! Cylinder choices:		
Doctor's/Office short cylinder, 16 " – 21" (6176)	6000-628-CS	
Doctor's/Office tall cylinder, 21 1/2" – 28" (6176) foot ring with support highly recommended	6000-604-CS	
Assistant's short cylinder, 18" – 23" (6187)	6000-605-CS	
Foot ring with support (6176)	6000-602	

! Often overlooked

Equipment dimensions, see page 109

Shipping weight estimate:

6176 48 lb / 23 kg

6187 60 lb / 28 kg

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools

Lights and Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

Forest

LIGHTS AND MONITOR MOUNTS



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

**Lights and
Monitor Mounts**

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks


Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

**Lights and
Monitor Mounts**

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

HALOGEN

Material: reflector material embedded into glass | Lexan light shield

Benefit: reflector will not peel or crack | shield is scratch resistant

Material: oil impregnated bronze bushing | 3-axis light head joint

Benefit: no light drift | 3-axis light keeps light pattern parallel w/the long axis of mouth while freeing up space above patient's head

LED

Material: energy efficient LED lamps

Benefit: creates a light pattern with feathered edges reducing eye fatigue | consumes approximately 75% less power than halogen lights | no heat | 22 -year projected LED life span

Material: cast aluminum light head

Benefit: allows for nine Designer Friendly® powder coat paint colors

Lights and Monitor Mounts

Post Mount Lights



Model 9083

Halogen Light

Part Number 9083



Model 9078

LED Light

Part Number 9078



Model 9073

LED Light

Part Number 9073

Standard Features

Halogen Light 9083

Automatic on/off activation w/override
3-Axis light head movement
Additional bulb
Color shade matching temperature setting
4200 to 5500 Kelvin
Designer Friendly® powder coat paint light head

LED Light 9078

Automatic on/off activation w/override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

LED Light 9073

Automatic on/off activation w/override
3-Axis light head movement
Six-mode, polycarbonate button keyboard
Easy release autoclavable handles
Light source 8 1-watt LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

Factory Options

LIGHT

Replacement bulb (9083/9093) 0014-300
Short rigid arm 20" in lieu of standard 32" 1110-227-CS

LIGHT POSTS:

44" with slot for sidebox access 1110-900
44" without slot, chair mount 1110-794

35" with post clamp 1110-820

Pivot (mounts on Forest 3900/3800 chair only)
(does not include cast chair adapter) 1113-654

Pivot (includes 3900/3800 cast chair adapter) 1113-654-CA

LIGHT POSTS, WITH MONITOR CORD SLOT:

44" with dual slot for sidebox access 1110-170
44" with single slot, chair mount 1110-050

35" with single slot & post clamp 1110-045

Part No

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Designer Friendly® Surface Treatments

STEP 1 SELECT POWDER COAT PAINT FOR UNIT & LIGHT

Arctic Silver

Bronze

Charcoal

Copper

Forest Grey

Midnight

Pewter

Shadow

White



Equipment dimensions, see page 110
Shipping weight estimate:
35 lb / 16 kg

Lights and Monitor Mounts

Track Mount Lights



Model 9083T

Halogen Light/Track
Shown with Optional 1305T Monitor Mount
Part Number 9083T



Model 1305T

Monitor Mount
Monitor Screen and Track Not Included
Part Number 1305T



Model 9078T

LED Light/Track
Part Number 9078T



Model 9073T

LED Light/Track
Part Number 9073T

Standard Features

Light and Track Assembly

Pre-wired through flex arm
Track/Trolley assembly
Vertical light tube w/6 inch horizontal side offset for clearance of optional monitor
Flex arm & light head assembly
Trolley mounting plate and shroud
! Drop post (specify mounting height)

Monitor Mount Assembly

Pre-wired from track through flex arm
Low voltage wiring (2 conductor 18 awg wire)
HDMI shielded cable
Lower positioning handle

Monitor Requirements:

12 volt w/power supply
weight: up to 18 lbs max, width of 19"



Equipment dimensions, see page 110

Shipping weight estimate:

Light/track 115 lb / 53 kg
1305T 140 lb / 64 kg (includes mount/track light weight)

Factory Options

LIGHT / TRACK

Monitor mount assembly
Replacement bulb (9083/9093)
S Video, cable pre-installed
RCA cables (3 conductor), pre-installed

MODEL CHOICE

1305T
0014-300
0014-282
0014-283

! Often overlooked

Standard Mounting Options

STEP 1 CHOOSE YOUR STANDARD MOUNTING OPTION:

MOUNTING HEIGHT:

- 8'
- 8' 6"
- 9'
- 9' 6"
- 10'
- 10' 6"
- 11'

○ _____ Custom mounting height ○ _____ * Custom Post Length ○ _____ * Custom Post Length

* Cannot exceed 44.25" and 56.25"

LIGHT POST LENGTH:

- 8.25"
- 14.25"
- 20.25"
- 26.25"
- 32.25"
- 38.25"
- 44.25"
- 50.25"
- 56.25"

MONITOR POST LENGTH:

- 20.25"
- 26.25"
- 32.25"
- 38.25"
- 44.25"
- 50.25"
- 56.25"

Designer Friendly® Surface Treatments

STEP 2 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:

- | | | | | |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|-----------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze | <input type="radio"/> Charcoal | <input type="radio"/> Copper | |
| <input type="radio"/> Forest Grey | <input type="radio"/> Midnight | <input type="radio"/> Pewter | <input type="radio"/> Shadow | |
| | | | | <input type="radio"/> White |

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Lights and Monitor Mounts

Ceiling Mount Lights



Model 9083-CL

Halogen Light

Part Number 9083-CL



Model 9078-CL

LED Light

Part Number 9078-CL



Model 9073-CL

LED Light

Part Number 9073-CL

Standard Features

Lights

Ceiling mounting plate and shroud w/drop post
! Specify mounting height

Halogen Light 9083-C

Automatic on/off activation w/override
3-Axis light head movement
Additional bulb
Color temperature, 5,500 Kelvin
Designer Friendly® powder coat paint light head

LED Light 9078-C

Automatic on/off activation w/override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

LED Light 9073-C

Automatic on/off activation w/override
3-Axis light head movement Six-mode
Polycarbonate button keyboard
Easy release autoclavable handles
Light source 8 1-watt LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

Factory Options

LIGHT

Replacement bulb (9083/9093)

! Often overlooked

MODEL CHOICE

0014-300

Part No

Standard Mounting Options

STEP 1 CHOOSE YOUR STANDARD MOUNTING OPTION: MOUNTING HEIGHT:

- 8'
- 8' 6"
- 9'
- 9' 6"
- 10'
- 10' 6"
- 11'
- _____ Custom mounting height
* Cannot exceed 48"

POST LENGTH:

- 12"
- 18"
- 24"
- 30"
- 36"
- 42"
- 48"
- _____ * Custom post length

Designer Friendly® Surface Treatments

STEP 2 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:

- Arctic Silver Bronze Charcoal Copper
- Forest Grey Midnight Pewter Shadow
- White



Equipment dimensions, see page 110

Shipping weight estimate:
50 lb / 23 kg

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Lights and Monitor Mounts

Cabinet and Wall Mount Lights

**Model 9079C**(Standard) 9079C-T (Tall)
LED Cabinet

Part Number 9079C

**Model 9079W**(Standard) 9079W-T (Tall)
LED Wall

Part Number 9079W

**Model 9074C**(Standard) 9074C-T (Tall)
LED Cabinet

Part Number 9074C

**Model 9093C**(Standard) 9093C-T (Tall)
Halogen Cabinet

Part Number 9093C

**Model 9093W**(Standard) 9093W-T (Tall)
Halogen Wall

Part Number 9093W

**Model 9074W**(Standard) 9074W-T (Tall)
LED Wall

Part Number 9074W

Standard Features

Lights

Cabinet mount, internal vertical mount w/post
Wall mount, vertical mounting plate
Standard rigid arm
Tall rigid arm (additional 1' drop)

Halogen Light 9093

Automatic on/off activation w/override
3-Axis light head movement
Additional bulb
Color shade matching temperature
setting 4200 to 5500 Kelvin
Designer Friendly® powder coat paint light head

LED Light 9079

Automatic on/off activation w/override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

LED Light 9074

Automatic on/off activation w/override
3-Axis light head movement Six-mode
Polycarbonate button keyboard
Easy release autoclavable handles
Light source 8 1-watt LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

Factory Options

LIGHT

Replacement bulb (9083/9093)
Cabinet mounting options: (for dimensions see page 95)
Horizontal top cabinet mount in lieu of internal mount
Vertical exterior cabinet mount in lieu of internal mount

MODEL CHOICE

0014-300
1110-881-CS
1110-886-CS

! **SPECIFY CABINET/WALL MOUNTING HEIGHT** _____

! Often overlooked

Part No

Forest 6400
PackagesForest 6400
Ordering GuidesForest 3900
PackagesForest 3900
Ordering GuideForest Fusion
PackagesForest Fusion
Ordering GuideSurface
TreatmentsChairs
and StoolsLights and
Monitor MountsDelivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and PartsGeneral
InformationEquipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
AccessoriesUtility Quick
Reference GuideUtility
Warranty

Sales Reps

Designer Friendly® Surface Treatments

STEP 1 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:

- | | | | |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze | <input type="radio"/> Charcoal | <input type="radio"/> Copper |
| <input type="radio"/> Forest Grey | <input type="radio"/> Midnight | <input type="radio"/> Pewter | <input type="radio"/> Shadow |
| <input type="radio"/> White | | | |



Equipment dimensions, see page 111

Shipping weight estimate:

Cabinet 60 lb / 28 kg

Wall 65 lb / 30 kg

Monitor Mount Assemblies



Model 1300RM Ring Post Mount

Part Number 1300RM



Model 1300TPM Top Post Mount

Part Number 1300TPM

Monitor Screens And Chair-Mounted Posts Not Included

Standard Features

Mount Assembly, Ring Post Mount 1300RM

ICW Quicklink monitor mount
Rex positioning handles
22" ring mount rigid arm

Monitor Mount Assembly, Top Post Mount 1300TPM

ICW Quicklink monitor mount
Rex positioning handles
22" top post mount rigid arm

Monitor Requirements:
12 volt with power supply
Weight: 8 1/2 to 13 lbs max, width of 19"

Factory Options

MONITOR MOUNT ASSEMBLY

- ! Light Post (1300RM):
 - 44" with dual slot, fixed sidebox mount, with Forest light bushing
 - 35" with slot & post clamp, fixed mount, with Forest light bushing
 - Pivot light posts & requirements
- ! MONITOR POST (1300TPM):
 - 30" with dual slot, fixed sidebox mount, with Oilite bushing
 - Pivot monitor posts & requirements
 - 30" with slot & post clamp, fixed mount, with Oilite bushing

Part No

MODEL CHOICE	Part No
	1110-170
	1110-045
(call to order)	
	1331-998
(call to order)	
	1331-997

! Often overlooked

Designer Friendly® Surface Treatments

STEP 1 SELECT POWDER COAT PAINT FOR UNIT & LIGHT:

- | | | | |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze | <input type="radio"/> Charcoal | <input type="radio"/> Copper |
| <input type="radio"/> Forest Grey | <input type="radio"/> Midnight | <input type="radio"/> Pewter | <input type="radio"/> Shadow |
| <input type="radio"/> White | | | |

Equipment dimensions, see page 111

Shipping weight estimate:

1300RM, 1300TPM 15 lb / 7 kg

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

Notes

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest

DELIVERY SYSTEMS



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

**Delivery Systems
7000 Series**

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories


Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



MIDNIGHT
Powdercoat Paint



GREYWOOD
Carbon Fiber



BRONZE
Powdercoat Paint



LASER LEATHER
Carbon Fiber

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

STANDARD IC

Material:	aluminum hinged lid 1/8" ABS plastic chassis with durable DR acrylic protective coating
Benefit:	hinge allows quick access to controls and ancillary devices / lid and chassis are resistant to UV, impact and compatible surface cleaners / aluminum lid is strong and lightweight
Material:	electro-less nickel plated brass air and water Ts, routing valves, pneumatic control and distribution blocks
Benefit:	20 year limited warranty on all parts "under the hood" (on all control heads); refer to page 107 for warranty details
Material:	handpiece holder bar with glass-filled nylon injection molded handpiece holders
Benefit:	allows for up to four handpiece positions, a syringe and ancillary touchpad at your desired angle and spacing
Material:	hinged aluminum control head lid
Benefit:	available in nine two-coat powder coat paint colors
Material:	plastic chassis
Benefit:	available in three carbon fiber look and two solid color finishes


EURO

Material:	whip arms rotate on nylon wheels with aluminum and brass components
Benefit:	feather-touch whip movement / ultralight handpiece feel / long, trouble free life
Material:	control head swivels at its center point
Benefit:	240-degrees of control head rotation allows for beside-or over-the-patient placement
Material:	forward mounted aluminum whip arm system
Benefit:	allows for greater reach of handpiece and syringe expanding control head placement options.
Material:	plastic control head lid
Benefit:	available in three carbon fiber looks and two solid color finishes
Feature:	autoclavable silicone handpiece pad available in black and white blends with Designer Friendly® choices

Cabinetry Mounting Configurations

Made for Each Other

From the ground up, Forest delivery systems are designed for easy cabinet integration.



3486 Beneath Counter Mount (BC)

Start with a strong foundation

- Your choice of 3/4" hardened steel 11" x 14-1/2" or optional 10" x 24" mounting plate that bolts to the bottom of your counter top.
- Followed by a 6" x 6" leveling plate with friction adjustment.

3486 Under Cabinet Floor Mount (UC)

- Low profile 1/2" steel floor mount plate allows room for your stool casters.
- Needle roller bearings with friction adjustment eliminates drift.
- Durable 2-1/2" steel tube post and arm assembly with needle roller bearings installed between pivoting sections.

7000 Series Beneath Counter Mount (BC)


Start with a strong foundation

- Your choice of 3/4" hardened steel 11" x 14-1/2" or optional 10" x 24" mounting plate that bolts to the bottom of your counter top
- Followed by a 6" x 6" leveling plate and a set of needle roller bearings with friction adjustment.
- Durable 2-1/2" steel tube dual arm assembly with needle roller bearings installed between pivoting arm sections.

7000 Series Under Cabinet Floor Mount (UC)

Start with a strong foundation

- Low profile 1/2" steel floor mount plate allows room for your stool casters.
- Needle roller bearings with friction adjustment eliminates drift.
- Durable 2-1/2" steel tube post and arm assembly with needle roller bearings installed between pivoting sections.



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories


Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Rear Mount Side Delivery Systems - 7000 Series

7020PRO Duo Swing Delivery System



Model 7020PRO-BC

Beneath Counter Mount

Part Number

7020PRO-BC

Universal heavy-duty mounting plate 11" x 14 1/2"
included (mounts to under counter)

Specify Cabinet Manufacturer

Standard Features

Duo Delivery System 7020PRO Base Model

Duo control

- automatic control for three handpieces
- handpiece lockouts
- 3-way air water syringe w/tubing, 5'
- 4-hole handpiece tubing, 5'
- pneumatic master control block w/diaphragm
- master on/off
- individual handpiece pressure adjustment
- individual water coolant adjustment
- non-retraction water coolant
- air coolant adjustment
- handpiece flush
- handpiece exhaust canister/oil collector
- wet/dry foot control
- Self-contained clean water system w/quick disconnect,
single bottle, 1 1/2 liter
- "Burr Buddy" tray holder and covered Zirc tray
- Chrome accents
- Central vacuum accessory package
 - central vacuum canister with disposable solids collector
 - 3-way air/water syringe w/tubing, 5'
 - autoclavable standard HVE lever valve w/tubing, 5'
 - autoclavable saliva ejector lever valve w/tubing, 5'
 - dual assistant's quad holder assembly
- Solid work surface, 17" x 30"
- Pneumatic work surface easy lift/close mechanism
- Height adjustable work surface, 5"
- Large and small organizer bins w/frame, dividers
and sharps container
- Procedure shelf
- Keyboard tray
- Stainless trash/CPU tray (MC)
- Universal heavy duty mounting plate 11" x 14 1/2",
beneath counter (BC)
- Wall mount panel 16" x 18" (WM)
- Standard utility center housing and cover (MC)
- Air pressure regulator w/filter
- Master air shut-off valve
- Umbilical, 6'



Beneath Counter Mount



Mobile Cart



Under Cabinet Floor Mount



Wall Mount

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Mounting plate dimensions, see pages 112

Shipping weight estimate:

7020PRO BC, MC, UC, WM 100 lb / 46 kg


Rear Mount Side Delivery Systems

7021PRO Duo Swing Delivery System

Handpieces Not Included



7021PRO



Model 7021PRO-BC

Beneath Counter Mount

Part Number

7021PRO-BC

Universal heavy-duty mounting plate 11" x 14 1/2"
included (mounts to under counter)

Specify Cabinet Manufacturer

Standard Features

Duo Delivery System 7021PRO

Duo control

automatic control for three handpieces
handpiece lockouts
3-way air water syringe w/tubing, 5'
4-hole handpiece tubing, 5'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
wet/dry foot control
Self-contained clean water system w/quick disconnect,
single bottle, 1 1/2 liter
"Burr Buddy" tray holder and covered Zirc tray
Chrome accents
Central vacuum accessory package
central vacuum canister with disposable solids collector
3-way air/water syringe w/tubing, 5'
autoclavable standard HVE lever valve w/tubing, 5'
autoclavable saliva ejector lever valve w/tubing, 5'
dual assistant's quad holder assembly
Solid work surface, 17" x 30"
Pneumatic work surface easy lift/close mechanism
Height adjustable work surface, 5"
Large and small organizer bins w/frame, dividers
and sharps container
Universal heavy duty mounting plate 11" x 14 1/2",
beneath counter (BC)
Wall mount panel 16" x 18" (WM)
Standard utility center housing and cover (MC)
Air pressure regulator w/filter
Master air shut-off valve
Umbilical, 6'



Beneath Counter Mount



Mobile Cart



Under Cabinet Floor Mount



Wall Mount

Mounting plate dimensions, see pages 112-121

Shipping weight estimate:

7021PRO BC, MC, UC, WM 100 lb / 46 kg

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories


Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Rear Mount Side Delivery Systems

7021PRO Duo Swing Delivery System



Model 7021PRO-MC

Mobile Cart Mount

Part Number 7021PRO-MC

Model 7021PRO-UC

Under Cabinet Floor Mount

Part Number 7021PRO-UC

Model 7021PRO-WM

Wall Mount

Part Number 7021PRO-WM

! Designer Friendly® Surface Treatments

STEP 1 SELECT POWDER COAT PAINT:

- Arctic Silver Bronze Charcoal
- Copper Forest Grey Midnight
- Pewter Shadow White

STEP 2 SOLID WORK SURFACE COLORS:

Corian® Antarctica

Wilsonart®: Arctic Melange Baja Melange

Meganite®: Cinder Granite Goldrush Gemstone
 Moonstone Platinum Granite Rocky Road Granite
 Silver Mist Starry Starry Night

Cabinet mounting configurations, see page 69

Factory Options

DELIVERY SYSTEM

Factory Installed STAR ETorque Electric System

Part No

STAR-ETORQUE

Handpiece illumination system

- One handpiece (5-hole), accommodates 1 add'l tubing 1700-021-CS
- Add'l handpiece illumination tubing (5-hole) 0017-615-CS
- One handpiece (6-pin), accommodates 1 add'l tubing 1700-019-CS
- Add'l handpiece illumination tubing (6-pin) 0017-603-CS

! Ancillaries use existing handpiece position (tubing lengths may vary) (see page 85)

Write in selection here _____

Write in selection here _____

! Handpiece Bundle option (see pages 82-83)

Write in selection here _____

Additional handpiece position, wet 1100-189-CS

Additional handpiece position, air drive only 1100-190-CS

Additional standard autoclavable HVE with tubing,
5' & autoclavable oral dry cup 1400-113-CS

Hands-Free HVE, see page 103, Write in selection here _____

Quick disconnect with flow control 1/4" (for water) 0016-454-CS

Quick disconnect with flow control 3/8" (for air) 0016-456-CS

! Procedure shelf

0010-864-CS

! Keyboard tray

0010-862-CS

! Stainless trash/CPU tray, H-frame (MC)

0010-863

! Monitor mount assembly (monitor screen not included)

base model requires purchase of 0010-864-CS 1300-025-CS

! Swivel tray

1110-909-CS

Large bin divider kit (white or black)

0010-860 / 0010-860-BLK

Small bin divider kit (white or black)

0010-861 / 0010-861-BLK

Sharps container kit (white or black)

0010-870 / 0010-870-BLK

! Touchpad for 3900 chair (25 ft. Cable)

1700-081-1-CS

Standard utility center housing & cover (BC/UC/WM)

1115-020

In floor recessed utility center housing & cover

1115-305

Additional umbilical (per foot)

1115-999-CS

Routing valve kit with water shutoff &

regulator/factory installation

Universal mounting plate, 10" x 24"
in lieu of 11" x 14 1/2" (BC)

0010-084-CS

△ Special color solid work surface

(call for price and leadtime) 0023-856-SSV

! Often overlooked

△ Lead time varies

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps


Rear Mount Side Delivery Systems

7010 Duo Swing Delivery System



1110-185-CS and 1110-908-CS

Shown with optional storage organizer
and swivel tray



Shown with optional
handpiece lockouts 0010-748
and additional assistant's quad
holder assembly 0024-156

Model 7010BC

Beneath Counter Mount

Part Number

7010BC

Universal heavy-duty mounting plate 11" x 14 1/2"
included (mounts to under counter)

Specify Cabinet Manufacturer

Standard Features

Duo Swing Delivery System 7010

Duo control

automatic control for three handpieces
3-way air water syringe w/tubing, 4'
4-hole handpiece tubing, 4'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
self-contained clean water system w/quick disconnect,
single bottle, 1 1/2 liter
wet/dry foot control
Height adjustable work surface, 5" (BC/UC/WM)
"Burr Buddy" tray holder and covered Zirc tray
Pneumatic work surface easy lift/close mechanism
Central vacuum accessory package
central vacuum canister w/disposable
solids collector
3-way air/water syringe w/tubing, 5'
autoclavable standard HVE lever valve w/tubing, 5'
autoclavable saliva ejector lever valve w/tubing, 5'
Solid work surface, 19" x 24"
Universal heavy duty mounting plate
11" x 14 1/2", beneath counter (BC)
Wall mount panel 16" x 18" (WM)
Adjustable height frame (MC)
Twin wheel hard floor casters (MC)
Standard utility center housing and cover (MC)
Air pressure regulator w/filter
Master air shut-off valve
Umbilical, 7'
External umbilical, 6' (MC)




Beneath Counter Mount



Mobile Cart



Under Cabinet Floor Mount



Wall Mount

Mounting plate dimensions, see pages 116-119

Shipping weight estimate:

7010BC, MC, UC, WM 130 lb / 59 kg

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories




Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Rear Mount Side Delivery Systems

7010 Duo Swing Delivery System



Model 7010MC

Mobile Cart Mount

Part Number 7010MC

Model 7010UC

Under Cabinet Floor Mount

Part Number 7010UC

Model 7010WM

Wall Mount

Part Number 7010WM

! Designer Friendly® Surface Treatments

STEP 1 SELECT POWDER COAT PAINT:

- | | | |
|-------------------------------------|-----------------------------------|--------------------------------|
| <input type="radio"/> Arctic Silver | <input type="radio"/> Bronze | <input type="radio"/> Charcoal |
| <input type="radio"/> Copper | <input type="radio"/> Forest Grey | <input type="radio"/> Midnight |
| <input type="radio"/> Pewter | <input type="radio"/> Shadow | <input type="radio"/> White |

STEP 2 SOLID WORK SURFACE COLORS:

Corian® Antarctica

Wilsonart®: Arctic Melange
 Baja Melange

Meganite®: Cinder Granite
 Goldrush Gemstone
 Moonstone
 Platinum Granite
 Rocky Road Granite
 Silver Mist
 Starry Night

Cabinet mounting configurations, see page 69

Factory Options

DELIVERY SYSTEM

Factory Installed STAR ETorque Electric System

Handpiece illumination system

One handpiece (5-hole), accommodates 1 add'l tubing

Part No

STAR-ETORQUE

1700-021-CS

0017-615-CS

1700-019-CS

0017-603-CS

! Ancillaries use existing handpiece position (tubing lengths may vary) (see page 85)

Write in selection here _____

Write in selection here _____

! Handpiece Bundle option (see pages 82-83)

Write in selection here _____

Additional handpiece position, wet

1100-189-CS

Additional handpiece position, air drive only

1100-190-CS

Lockout feature for 3 handpieces

0010-748

Additional standard autoclavable HVE with tubing,
5' & autoclavable oral dry cup

1400-113-CS

Hands-Free HVE, see page 103, Write in selection here _____

0024-156

Additional assistant's quad holder assembly

0016-454-CS

Quick disconnect with flow control 1/4" (for water)

0016-456-CS

Quick disconnect with flow control 3/8" (for air)

1110-185-CS

! Storage organizer with bins & frame

1110-908-CS

! Swivel tray (requires purchase of 1110-185-CS)

1110-079-CS

! Swivel tray (mounts without storage organizer)

1300-025-CS

Monitor mount assembly (monitor screen not included)
base model requires purchase of 0010-864-CS

1110-909-CS

Swivel tray

1700-081-1-CS

Touchpad for 3900 chair (25 ft. Cable)

9000-998-CS

Chrome accents

Standard utility center housing & cover (BC/UC/WM)

1115-020

In floor recessed utility center housing & cover

1115-305

Additional umbilical (per foot)

1115-999-CS

Routing valve kit with water shutoff &

1115-315-CS

regulator/factory installation

Universal mounting plate, 10" x 24"
in lieu of 11" x 14 1/2" (BC)

0010-084-CS

Standard color laminate work surface

0023-776-LAM

in lieu of solid surface

△ Special color solid or laminate work surface
(call for price and leadtime)

0023-776-SSV/

0023-776-X

! Often overlooked

△ Lead time varies

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories


Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Rear Mount Side Delivery Systems

3486 Swing Delivery System



1110-179-CS and 1110-908-CS

Shown with optional storage organizer and swivel tray

Model 3486BC

Beneath Counter Mount

Part Number

3486BC

Universal heavy-duty mounting plate 11" x 14 1/2" included (mounts to under counter)

Specify Cabinet Manufacturer

Standard Features

Assistant's Swing Delivery System 3486 Series

Central vacuum accessory package

arm-mounted central vacuum canister
w/disposable solids collector

3-way air/water syringe w/tubing, 5'
autoclavable standard HVE lever valve w/tubing, 5'
autoclavable saliva ejector lever valve w/tubing, 5'
self-contained clean water system w/quick disconnect,
single bottle, 1 1/2 liter

Assistant's quad holder assembly

Solid work surface, 19" x 24"

Beneath counter rigid arm assembly

w/heavy duty mounting plate, 11" x 14 1/2" (BC)

Adjustable height frame (MC)

Twin wheel hard floor casters (MC)

Under cabinet floor mount delivery post (UC)

Wall mount panel, 16" x 18" (WM)

Standard utility center housing and cover (MC)

Air pressure regulator w/filter

Master shut-off valve

Umbilical, 7'

electrical outlet/UL (MC)

Note: Arms can be reversed on the beneath counter and wall mount to adjust work surface height.



Beneath Counter Mount



Mobile Cart



Under Cabinet Floor Mount



Wall Mount

Mounting plate dimensions, see pages 118-119

Shipping weight estimate:

3486 115 lb / 53 kg

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories




Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Rear Mount Side Delivery Systems

3486 Swing Delivery System



Model 3486MC

Mobile Cart Mount

Part Number 3486MC

Model 3486UC

Under Cabinet Floor Mount

Part Number 3486UC

Model 3486WM

Wall Mount

Part Number 3486WM

! Designer Friendly® Surface Treatments

STEP 1 SELECT POWDER COAT PAINT:

- Arctic Silver Bronze Charcoal
- Copper Forest Grey Midnight
- Pewter Shadow White

STEP 2 SELECT UNIT PLASTIC FINISHES:

- Forest Grey Greywood Laser Leather
- Lemans White

STEP 3 SOLID WORK SURFACE COLORS:

Corian® Antarctica

Wilsonart®: Arctic Melange Baja Melange

Meganite®: Cinder Granite Goldrush Gemstone
 Moonstone Platinum Granite
 Rocky Road Granite Silver Mist
 Starry Night

Factory Options

Part No

DELIVERY SYSTEM

Additional standard autoclavable HVE with tubing,
5' & autoclavable oral dry cup

1400-113-CS

Hands-Free HVE, see page 103, Write in selection here

Additional assistant's quad holder assembly

0024-149

Quick disconnect with flow control 1/4" (for water)

0016-454-CS

Quick disconnect with flow control 3/8" (for air)

0016-456-CS

! Storage organizer with bins & frame

1110-179-CS

! Swivel tray (requires purchase of 1110-179-CS)

1110-908-CS

! Swivel tray (mounts without storage organizer)

1110-482-CS

Monitor mount assembly (monitor screen not included)
requires purchase of 1110-179-CS

1300-025-CS

Touchpad for 3900 chair (25 ft. Cable)

1700-081-1-CS

Standard utility center housing & cover (BC/UC/WM)

1115-020

In floor recessed utility center housing & cover

1115-305

Additional umbilical (per foot)

1115-999-CS

Routing valve kit with water shutoff &
regulator/factory installation

1115-315-CS

Standard color laminate work surface
in lieu of solid surface

0023-003-LAM

△ Special color solid or laminate work surface
(call for price and leadtime)

0023-003-X/

0023-003-SSV

! Often overlooked

△ Lead time varies

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty


Sales Reps

Cabinet mounting configurations, see page 69

Rear Mount Side Delivery Systems

Storage Organizer, Trays and Monitor Mount

Storage Organizer with Bins and Frame



Factory option model: 7010 1110-185-CS
Factory option models: 2184BC/2186BC/3486 1110-179-CS

- Metal powder coated organizer frame
- Three storage bins with removable covers
- Storage dividers

Swivel Tray



Factory option models: 1110-908-CS

2184BC/2186BC/3486/7010 requires purchase of 1110-179-CS or 1110-185-CS

- Inserts inside storage organizer tower
- Metal powder coated tray holder
- Plastic instrument tray

Swivel Tray with Left or Right Mounting Bracket



Factory option model: 7010 1110-079-CS
Factory option models: 2184BC/2186BC/3486 1110-482-CS

- Metal powder coated organizer frame
- Three storage bins with removable covers
- Storage dividers

Swivel Tray with Left and Right Mounting Brackets



Factory option models: 1110-909-CS

- Metal powder coated tray holder towers and brackets
- Plastic instrument tray

PRO Series and 7010 third tray addition mounts with organizer

Monitor Mount



Factory option models: 1300-025-CS

- Inserts inside storage organizer or PRO Series procedure shelf tower
- Forest's unique locking tilt mechanism with knob 1/2" pin

PRO Series 7021 requires purchase of 0010-864-CS
7010 requires purchase of 1110-185-CS
2184BC/2186BC/3486 requires purchase of 1110-179-CS

Monitor requirements (not included):
12 volt with power supply
Weight 8-1/4 to 13 lbs. max, width of 19"

Forest 6400 Packages	Forest 6400 Ordering Guides
Forest 3900 Packages	Forest 3900 Ordering Guide
Forest Fusion Packages	Forest Fusion Ordering Guide
Surface Treatments	Lights and Monitor Mounts
Chairs and Stools	Delivery Systems 7000 Series
Simplified	Delivery Systems
Factory Options and Parts	General Information
Equipment Warranty	Compressors
Compressors	Dry Vacuums
Tanks	Wet Vacuums
Utility Room Accessories	Utility Quick Reference Guide
Utility Warranty	Sales Reps

Ancillary Ordering and Options

Standard IC Delivery Unit

Standard IC


Please use graph below to determine number of ancillaries and options possible for the Standard IC.



Standard Features

Automatic control for three handpieces
3-way air/water syringe with tubing, 5'
4-hole handpiece tubing, 5';
pneumatic master control block with diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air/water coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
handle brake release function

% Standard IC Control Head Space Absorbed per Option



Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

* Uses existing handpiece position

★ Standard on SE and SI deliveries

+ 3900 Chair touchpad cable standard on PI and SI deliveries

◆ Two or more may require purchase of 0010-310 ancillary shelf (see page 86)

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options

and Parts

General

Information

Equipment

Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room

Accessories

Utility Quick

Reference Guide

Utility

Warranty

Sales Reps

Euro Delivery Unit


Euro

Please use graph below to determine number of ancillaries and options possible for the Euro.



Standard Features

Automatic control for four handpieces
3-way air/water syringe with tubing, 3'
4-hole handpiece tubing, 3'
pneumatic master control block with diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air/water coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
dual handle brake release function
swivel instrument tray



Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

* Uses existing handpiece position

★ Standard on SE and SI deliveries

+ 3900 Chair touchpad cable standard on PI and SI deliveries

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps

Fusion Delivery Unit

Fusion


Please use graph below to determine number of ancillaries and options possible for the delivery system.



Standard Features

Automatic control for three handpieces
3-way air/water syringe w/tubing, 5'
4-hole handpiece tubing, 5'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
handle brake release function

% Control Head Space Absorbed per Option



Please note:

If your ancillary and option choices exceed 100%, please call customer service for additional information.

⌘ Uses existing handpiece position

★ Standard on SI deliveries

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps

Handpiece Bundles

Save time and money with handpiece bundles

Our Star® handpiece bundles make it easy to be set up for success. Each bundle package is designed to complete a new operatory setup with easy ordering and tremendous savings.

Bundle Packages	Bundle Details	Part No	Quantity	
Star 500 fiber optic Bundle				
Part Number	268925	Star 500 Elite LubeFree, fiber optic Star 500 Swivel, fiber optic	266040 266064	3 1
430 Torque fiber optic Bundle				
Part Number	268927	430SWL Torque LubeFree, fiber optic LED Plug 'N Play Swivel	265890 265646	3 1
430 Series fiber optic Bundle				
Part Number	268928	430SWL Satin LubeFree, fiber optic LED Plug 'N Play Swivel	268840 265646	3 1
430K Series fiber optic Bundle	430K LubeFree, fiber optic	262052	3	
Part Number	268929			

Star fiber optic bundles require 6-pin tubing configuration and fiber optic light source.
All other fiber optic bundles require fiber optic light source.

NOTES:

Cannot be combined with other Star promotions or discounts.

*Savings calculated off the total retail price for the individual items in each package.

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps

Bundle Packages	Bundle Details	Part No	Quantity
Star Titan® Low-Speed Bundle			
Part Number	268934		
	Star Titan 20K Motor	266271	1
	Titan 3 Motor To Angle Adaptor	264073	3
	Titan LubeFree Auto Latch Angle	264146	3
	Titan 3 Straight Nose Cone	257506	1
	Titan 4-Line Swivel	260144	1
Titan 3 Low-Speed Bundle			
Part Number	268930		
	Titan 3 20K Motor	262262	1
	Titan 3 Motor to Angle Adaptor	264073	3
	Titan LubeFree Auto Latch Angle	264146	3
	Titan 3 Straight Nose Cone	257506	1
	Titan 4-Line Swivel	260144	1
Star 500 non-fiber optic Bundle			
Part Number	268931		
	Star 500 Elite LubeFree, non-fiber optic	266041	3
	Star 500 Swivel, non-fiber optic	266065	1
430 Torque non-fiber optic Bundle			
Part Number	268932		
	430SW Torque LubeFree	265891	3
	4-line non-fiber optic Swivel	261547	1
430 Series non-fiber optic Bundle			
Part Number	268933		
	430SW Satin LubeFree	268841	3
	4-line non-fiber optic Swivel	261547	1

Star fiber optic bundles require 6-pin tubing configuration and fiber optic light source.
All other fiber optic bundles require fiber optic light source.

NOTES:

Cannot be combined with other Star promotions or discounts.

*Savings calculated off the total retail price for the individual items in each package.

Control Head Accessories

Ancillary Products

Star ETorque™ Electric System



Factory Option

STAR-ETORQUE

Quiet and powerful, with a brushless motor reaching up to 200,000 rpm, this electric system can perform all procedures. It has a user-friendly display module that can be customized as needed and has a fully autoclavable motor. Handpiece attachments not included in price.

To order StarETorque and accessories, see the Star handpiece catalog.

Star® Titan® Blis-sonic™ Scaler



Blis-sonic SW

Blis-sonic K

To order Titan Blis-sonic scalers and accessories, see the Star handpiece catalog.

Blis-sonic SW Scaler
(Includes 4 line swivel, Star PN 261547)

STAR-BLIS-SW

Blis-sonic K Scaler (Fixed Backed)

STAR-BLIS-K

Universal scaler tip

- Use for all Surfaces, medium and heavy calculus

Sickle scaler tip

- Recommended for sensitive patients
- Less energetic than the Universal Tip
- Less suitable for deep root scaling

Perio scaler tip

- Recommended for lingual supragingival scaling of mandibular incisors and deep local pockets
- Provides same energy as Universal Tip

Satelec SP



Factory option

SATELEC-S

Factory option

SATELEC-S-LED

- Uses existing handpiece position
- Satelec SP Newtron LED lit Piezo Ultrasonic Scaler w/power supply
- Color coding system ensures proper tip and power setting
- Three tips: #1, #10Z, H3 (Intro Scaling Kit)
- Tips available for periodontics, endodontics, scaling, restorative, and cavity preparation
- Cruise Control® automatically adjusts frequency amplitude
- On/off controlled by rheostat
- First ever autoclavable LED lit handpiece – 100,000 LUX
- Newtron technology means power, control and reliability
- Satelec warranty 1 year / tips 30 days from date of purchase
- Please call 800-289-6367 for warranty and technical support of the Scaler

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Ancillary Products

Cavitron Scaler



Factory option

CAVITRON

- Uses existing handpiece position Dentsply Cavitron Magnetostrictive Ultrasonic scaler
- Blue Zone® for patient comfort
- SPS® (Sustained Performance System) constant balance between scaling efficiency and patient comfort
- Handpiece w/330 degree swivel and lavage control (reduces cord drag)
- Autoclavable Steri-mate® (infection control)
- PowerLINE 1000 30K insert
- Frequency 30 kHz
- Cavitron warranty 2 years from date of purchase
- Please call 800-289-6367 for warranty and technical support of the Scaler

HP Illumination System, 5 Hole with NEW integrated LED bulb

HP Illumination System, 6 Pin



5 Hole

6 Pin

Factory option

1700-021-CS

Factory option

1700-019-CS

- Uses existing handpiece position
- Dual channel outputs individually adjustable from 2.5VDC to 4.0VDC
- Accommodates one additional handpiece tubing
- Illumination delay (adjustable between 0-15 seconds) can be set individually for 2 HP
- Automatic activation of HP illumination w/the dental unit foot control
- Simple installation into any new or existing dental unit
- Indicator LEDs confirms correct operation of the control module
- Powers all brands of optic handpieces and air-motors
- ISO standard handpiece terminal with stainless steel nut
- Handpiece tubing, (4' to 5' installed)
- Power transformer, output 24VAC, 1.05 Amps
- Installation hardware and instructions

HP Illumination Tubing, 5 Hole with NEW integrated LED bulb



Factory option

0017-615-CS

- Used to replace power optic handpiece tubing
- ISO standard power optic handpiece terminal with stainless steel nut
- Handpiece tubing, (4' to 5' factory installed)

HP Illumination Tubing, 6 Pin



Factory option

0017-603-CS

- Used to replace power optic handpiece tubing
- ISO standard power optic handpiece terminal with stainless steel adapter and nut
- Handpiece tubing, (4' to 5' factory installed)

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Control Head Accessories

Ancillary Products

Chair Control, Touchpad, Unit-Mounted



Factory option

3914-998-CS

- 4 programmable preset positions

Additional Chair Control, Touchpad, Unit-Mounted (3900/3800)

Factory option

3914-999-CS

Assistant's Hygiene Instrumentation



Factory option

1408-CS

- Canister with 1/2" drain
- Mounting bracket
- Drain tubing
- Autoclavable standard HVE w/tubing, 5'
- Autoclavable saliva ejector w/tubing, 5'
- Standard vacuum holder
- Universal holder
- Installed on Standard IC and Value control heads (not available on Euro)

Ancillary Shelf



Factory option

0010-310

- Accommodates additional ancillary products
- Designed for easy access to installed devices
- Powder coated aluminum
- Installed on Standard IC control head only

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Ordering Forest Equipment Made Easy **Simplified**

Too many options can sometimes be overwhelming. That's why this select product offering has just two color choices. You can make upgrading a treatment room or matching existing equipment easy.

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps

Simplified Delivery Systems, Chair Mounted

Fixed Mount with Sidebox, External Umbilical

Choose Your Control Head



Fusion Control Head

Pnuematic 3 HP Syringe

Part Number

4295SE

Standard IC Control Head

Pnuematic 3 HP Syringe

Part Number

4485SE

Euro Control Head

Pnuematic 4 HP Syringe

Part Number

4987SE

Delivery Systems can be mounted on any 2" fixed post mount brand of chair

Optional Equipment Choose One:

1577

Cuspidor, Deluxe

(Sidebox comes with delivery)



1679

Cuspidor, Basic

(Sidebox comes with delivery)



9083

Light, Halogen

1110-900

Light Post



9078

Light, LED

1110-900

Light Post



9073

Light, LED

1110-900

Light Post



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Standard Features	Factory Options	Part No	
Delivery System	DELIVERY SYSTEM	MODEL CHOICE	Forest 6400 Packages
Choice of Control Head	Factory Installed STAR ETorque Electric System	STAR-ETORQUE	Forest 6400 Ordering Guides
Automatic control for handpieces 3-way air/water syringe with tubing, 5' (Euro tubing 3') 4-hole handpiece tubing, 5'; (Euro tubing 3') pneumatic master control block w/diaphragm master on/off	Handpiece illumination system One handpiece (5-hole), accommodates 1 add'l tubing Add'l handpiece illumination tubing (5-hole) One handpiece (6-pin), accommodates 1 add'l tubing Add'l handpiece illumination tubing (6-pin)	1700-021-CS 0017-615-CS 1700-019-CS 0017-603-CS	Forest 3900 Packages
individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust canister/oil collector single handle brake release function dual handle brake release function (Euro) swivel instrument tray (Euro) wet/dry foot control - Shadow powder coat standard	! Ancillaries use existing handpiece position (see page 85) Write in selection here _____ Write in selection here _____		Forest 3900 Ordering Guide
Flexible control arm with pneumatic brake vertical travel 28" / horizontal travel 62"	! Handpiece Bundle option (see pages 82-83) Write in selection here _____		Forest Fusion Packages
Sidebox with self-contained clean water system w/ quick disconnect, single bottle, 1 1/2 liter electrical outlet/UL	Additional handpiece position, wet (4485) Additional handpiece position, air drive only (4485) Swivel tray, arm-mounted (over control head) Less sidebox telescoping arm and assistant's instrumentation Chrome accents Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup	1100-189-CS 1100-190-CS 1110-468 1400-000-CS 9000-998-CS 1400-113-CS	Forest Fusion Ordering Guide
Water quick disconnect 1/4" w/flow control Air quick disconnect 3/8" w/flow control Toggle valve for city/bottle water routing Telescoping arm with swivel quad holder bar Central vacuum accessory package central vacuum canister with disposable solids collector 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5'	Hands-Free HVE, see page 103, Write in selection here _____ Additional umbilical (maximum 9 feet) (per foot) 3800/3900 chair fixed 2" post mount chair adapter	1115-999-CS 1115-195	Surface Treatments
External umbilical, 5' Deluxe utility center air pressure regulator w/filter master air shut-off valve ! Chair adapter sold separately	△ non 3800/3900 2" post mount chair adapter, fixed (Only sold with unit) (specify chair make and model; non-cancellable/non-returnable)	1115-099-CA	Chairs and Stools
Deluxe Cuspidor 1577 One-piece ceramic bowl with touch lever water activator Automatic cup fill Timed bowl rinse	CUSPIDOR , (sold only with delivery) DELUXE <input type="radio"/> 1577 OR BASIC <input type="radio"/> 1679 Deluxe cuspidor unit mount activation (1577) LIGHT, HALOGEN (mounting options, see page 61) <input type="radio"/> 9083 OR LED <input type="radio"/> 9078 OR LED <input type="radio"/> 9073	1100-480-CS	Lights and Monitor Mounts
Basic Cuspidor 1679 Removable ceramic bowl Manual cup fill Timed bowl rinse	•! Light post (specify light make/model) Light post, 44" (Forest) Light post, 44" (Pelton & Crane) MONITOR MOUNT ASSEMBLY (monitor screen not included) (see page 65) Ring mount (used with light post) Top post mount (used with monitor post) • Monitor post (for top post mount)	1110-900 1110-890 1300RM 1300TPM 1331-998	Delivery Systems 7000 Series
	• Choose only one post when sold with unit ! Often overlooked △ Lead time varies		Delivery Systems
			Simplified
			Factory Options and Parts
			General Information
			Equipment Warranty
			Compressors
			Dry Vacuums
			Tanks
			Wet Vacuums
			Utility Room Accessories
			Utility Quick Reference Guide
			Utility Warranty
			Sales Reps

!

Simplified Surface Treatments**4295SE AVAILABLE IN FOREST GREY ONLY****UNIT PLASTIC FINISHES:** Forest Grey White**POWDER COAT PAINT FINISHES:** Forest Grey White**Equipment dimensions, see page 108**

Shipping weight estimate:

4295SE & 4485SE Delivery Units 105 lb / 48 kg

4987SE Delivery Unit 110 lb / 50 kg

Fixed Mount without Sidebox, External Umbilical

Choose Your Control Head



Fusion Control Head

Pnuematic 3 HP Syringe

Part Number

4295FE

Standard IC Control Head

Pnuematic 3 HP Syringe

Part Number

4485FE

Euro Control Head

Pnuematic 4 HP Syringe

Part Number

4987FE

Delivery Systems can be mounted on any 2" fixed post mount brand of chair

Optional Equipment Choose One:

1075FE

Assistant's Telescoping Arm with Instrumentation



1636FE

Cuspidor, Basic



1637FE

Cuspidor, Basic with Assistant's Instrumentation



9083

Light, Halogen

1110-820

Light Post



9078

Light, LED

1110-820

Light Post



9073

Light, LED

1110-820

Light Post



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Standard Features	Factory Options	Part No	Forest 6400 Packages
Delivery System Choice of Control Head Automatic control for handpieces asepsis 3-way air/water syringe w/tubing, 5' (Euro tubing 3') handpiece tubing, 4-hole Midwest, asepsis style, 5' (Euro tubing 3') master control block master on/off individual handpiece pressure adjustment individual water coolant adjustment non-retraction water coolant air coolant adjustment handpiece flush handpiece exhaust canister/oil collector single handle brake release function self-contained clean water system w/ quick disconnect, single bottle, 1 1/2 liter wet/dry foot control Flexible control arm w/pneumatic brake vertical travel 28" / horizontal travel 62" External umbilical, 5' * Deluxe utility center air pressure regulator w/filter master air shut-off valve ! Chair adapter sold separately	DELIVERY SYSTEM Factory Installed STAR ETorque Electric System Handpiece illumination system One handpiece (5-hole), accommodates 1 add'l tubing Add'l handpiece illumination tubing (5-hole) One handpiece (6-pin), accommodates 1 add'l tubing Add'l handpiece illumination tubing (6-pin) ! Ancillaries use existing handpiece position (see page 85) Write in selection here _____ Write in selection here _____ ! Handpiece Bundle option (see pages 82-83) Write in selection here _____ Additional handpiece position, wet (4485) Additional handpiece position, air drive only (4485) Assistant's Hygiene Instrumentation (control head mount) Swivel tray, arm-mounted (over control head) Chrome accents Quick disconnect with flow control, 1/4" (specify air or water) Additional umbilical (maximum 9 feet) (per foot) Routing valve kit with water shutoff & regulator/factory installation Routing valve kit for use with cuspidor/factory installation 3800/3900 chair fixed 2" post mount chair adapter △ non 3800/3900 2" post mount chair adapter, fixed (Only sold with unit) (specify chair make and model; non-cancellable/non-returnable)	STAR-ETORQUE 1700-021-CS 0017-615-CS 1700-019-CS 0017-603-CS 1100-189-CS 1100-190-CS 1408-CS 1110-468 9000-998-CS 0016-455-CS 1115-999-CS 1115-315-CS 1115-316-CS 1115-195 1115-099-CA 1075FE 1400-113-CS ○ 1636FE ○ 1637FE 1400-113-CS ○ 9083 ○ 9078 ○ 9073 1110-820 1110-822 1300RM 1300TPM 1331-997	Forest 6400 Ordering Guides Forest 3900 Packages Forest 3900 Ordering Guide Forest Fusion Packages Forest Fusion Ordering Guide Surface Treatments Chairs and Stools Lights and Monitor Mounts Delivery Systems 7000 Series Delivery Systems
Assistant's Instrumentation 1075FE Central vacuum accessory package central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5' External umbilical, 5' Post mount, 2"	ASSISTANT'S INSTRUMENTATION Additional autoclavable standard HVE with tubing, 5' and autoclavable oral dry cup Hands-Free HVE, see page 103, Write in selection here _____	1075FE 1400-113-CS	Simplified
Cuspidor, Basic 1636FE Gravity drain cuspidor w/swing away arm timed bowl rinse manual cup fill removable ceramic bowl Post mount adapter, 2" External 5' umbilical Water pressure regulator w/filter Master water shut-off valve	* CUSPIDOR, BASIC with assistant's instrumentation Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup (1637FE) LIGHT, HALOGEN (mounting options, see page 61) OR LED OR LED •! LIGHT POST with 2" post clamp (specify light make/model) Light post, 35" (Forest) Light post, 35" (Pelton & Crane) MONITOR MOUNT ASSEMBLY (monitor screen not included) (see page 65) Ring mount (used with light post) Top post mount (used with monitor post) • Monitor post (for top post mount)	1400-113-CS ○ 1636FE ○ 1637FE 1400-113-CS ○ 9083 ○ 9078 ○ 9073 1110-820 1110-822 1300RM 1300TPM 1331-997	Factory Options and Parts General Information Equipment Warranty Compressors Dry Vacuums
Cuspidor, Basic 1637FE Same as 1636FE above, w/additional features: Individual assistant's instrumentation holders Central vacuum accessory package central vacuum canister w/disposable solids collector asepsis 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5'	• Choose only one post when sold with unit * When optional cuspidor added, standard utility center housing and cover is used to accommodate two umbilicals. ! Often overlooked △ Lead time varies		Tanks Wet Vacuums

! Simplified Surface Treatments

4295FE AVAILABLE IN FOREST GREY ONLY

UNIT PLASTIC FINISHES: ○ Forest Grey ○ White

POWDER COAT PAINT FINISHES: ○ Forest Grey ○ White

Equipment dimensions, see page 108

Shipping weight estimate:

4295FE & 4485FE Delivery Units 95 lb / 44 kg

4987FE Delivery Unit 100 lb / 46 kg

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility

Warranty

Sales Reps

Simplified Cuspidors and Assistant's Instrumentation

Fixed Chair Post Mount, External Umbilical



Model 1636

Part Number 1636



Model 1637

Part Number 1637



Model 1075

Part Number 1075

Standard Features

Cuspidor, Basic 1636

Gravity drain cuspidor w/swing away arm

timed bowl rinse

manual cup fill

removable ceramic bowl

Post mount adapter, 2"

External umbilical, 5'

! Does not include post or chair adapter

Cuspidor, Basic 1637

Same as 1636 above, w/additional features:

Individual assistant's instrumentation holders

Central vacuum accessory package

central vacuum canister w/disposable solids collector

3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

! Does not include post or chair adapter

Post Mount Vacuum Telescoping Arm 1075

Post mount adapter, 2"

Central vacuum accessory package

post mounted central vacuum canister w/disposable solids collector

3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5'

autoclavable saliva ejector lever valve w/tubing, 5'

External umbilical, 8'

! Does not include post or chair adapter

Factory Options

CUSPIDOR, BASIC

MODEL CHOICE _____
Additional autoclavable standard HVE with tubing,
5' & autoclavable oral dry cup (1637) 1400-113-CS

Hands-Free HVE, see page 103, Write in selection here _____

Quick disconnect with flow control, 1/4"
(specify air or water) 0016-455-CS

! Self-contained clean water system with quick disconnect
Standard utility center housing & cover 1700-310-CS
1115-020

► Standard utility center
Additional umbilical (per foot) 1115-001
1115-999-CS

! MOUNTING POSTS 2" diameter:
Control post with cap, 18" 1110-487

! LIGHT POSTS: 44" 2" diameter
Forest model 1110-794
Pelton & Crane 1110-785

ASSISTANT'S INSTRUMENTATION 1075
Additional standard autoclavable HVE with tubing,
5' & autoclavable oral dry cup 1400-113-CS

! Self-contained clean water system with quick disconnect
Assistant's air shut-off & regulator utility center
w/o housing & cover 1700-323-CS
1115-008

Standard utility center housing & cover
Additional umbilical (per foot) 1115-020
1115-998-CS

► Standard utility center includes: housing, cover, air & water pressure regulator with filter, master air and water shut-off valves.

! Often overlooked

WARNING: Isolated water sources (as opposed to city water hook-up) are highly recommended to allow full control of the source and quality of water circulating through the syringe and cuspidor water lines. If you do not have a water quality protocol, the optional self-contained water system allows you to isolate this product from the municipal water supply; however, it requires monitoring, maintenance and disinfection of the water lines, reservoir water bottle and source water.

Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

White

POWDER COAT PAINT FINISHES:

Forest Grey

White

Equipment dimensions, see page 108

Shipping weight estimate:

1636 17 lb / 8 kg

1637 20 lb / 9 kg

1075 15 lb / 7 kg

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility

Warranty

Sales Reps

Simplified Cuspidors and Assistant's Instrumentation

Rear Beneath Shelf or Work Surface Mount



Model 1284 - Mounts under floating work surface

Shown with optional Clean Water System 1700-312 not included in price Assistant's Utilities 1115-008 sold separately

Part Number 1284



Model 1286 - Mounts beneath shelf

Optional Clean Water System 1700-312 and Assistant's Utilities 1115-008 sold separately

Part Number 1286

Standard Features

Horizontal Assistant's Delivery System 1284/1286

Central vacuum accessory package, swing mount arm-mounted central vacuum canister w/disposable solids collector 3-way air/water syringe w/tubing, 5' autoclavable standard HVE lever valve w/ tubing, 5' autoclavable saliva ejector lever valve w/ tubing, 5' Horizontal surface mount Link assembly, arm to plate (1286) External umbilical, 7'

Factory Options

DELIVERY SYSTEM	MODEL CHOICE	Part No
Additional autoclavable standard HVE with tubing, 5' & autoclavable oral dry cup	1400-113-CS	
Hands-Free HVE, see page 103, Write in selection here		
Slide mount, 7 5/8" extension (1284)	1700-142	
! Standard utility center housing & cover	1115-020	
! Self-contained clean water system with quick disconnect, single bottle	1700-312	
! Assistant's air shut-off & regulator utility center w/o housing & cover	1115-008	

! Often overlooked

Part No

Forest 6400 Packages

Forest 6400 Ordering Guides

Forest 3900 Packages

Forest 3900 Ordering Guide

Forest Fusion Packages

Forest Fusion Ordering Guide

Surface Treatments

Chairs and Stools

Lights and Monitor Mounts

Delivery Systems 7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

⚠ WARNING: Isolated water sources (as opposed to city water hook-up) are highly recommended to allow full control of the source and quality of water circulating through the syringe and cuspidor water lines. If you do not have a water quality protocol, the optional self-contained water system allows you to isolate this product from the municipal water supply, however, it requires monitoring, maintenance and disinfection of the water lines, reservoir water bottle and source water.

!

Simplified Surface Treatments

Forest Grey

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Mounting plate dimensions, see page 111

Shipping weight estimate:

1284 12 lb / 5 kg

1286 14 lb / 6 kg

Simplified Rear and Side Delivery Systems

Side Cabinet & Wall Mount



Model 3235

12" Rigid Arm Cabinet Mount

Part Number 3235



Model 3230

22" Rigid Arm Wall Mount

Part Number 3230

Standard Features

Doctor Standard IC Delivery System 3235/3230

Standard IC control head
automatic control for three handpieces
3-way air/water syringe w/tubing, 5'
4-hole handpiece tubing, 5'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
single handle brake release function
self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
wet/dry foot control - Shadow powder coat standard
Flexible control arm w/pneumatic brake (3235)
cabinet vertical mount 12" rigid arm & 18" flex arm
Flexible control arm w/pneumatic brake (3230)
wall vertical mount 22" rigid arm & 18" flex arm
External umbilical, 5'
Standard utility center housing and cover (3230)
Air pressure regulator w/filter
Master air shut-off valve

Factory Options

DELIVERY SYSTEM

MODEL CHOICE	Part No
Factory Installed STAR ETorque Electric System	STAR-ETORQUE
Handpiece Illumination System	
One handpiece (5-hole), accommodates	
1 add'l tubing	1700-021-CS
Add'l handpiece illumination tubing (5-hole)	0017-615-CS
One handpiece (6-pin), accommodates	
1 add'l tubing	1700-019-CS
Add'l handpiece illumination tubing (6-pin)	0017-603-CS
Assistant's hygiene instrumentation (control head mount)	1408-CS
Quick disconnect with flow control, 1/4" (specify air or water)	0016-455-CS
Swivel tray, arm-mounted (over control head)	1110-468
Standard utility center housing & cover (3235)	1115-020
Additional umbilical (per foot)	1115-999-CS
Routing valve kit with water shutoff & regulator/factory installation	1115-315-CS
Wall mount panel, 16" x 18"	1110-586-CS

! Often overlooked

Part No

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

!

Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

Mounting plate dimensions, see page 112

Shipping weight estimate:
3235, 3230 50 lb / 23 kg

Simplified Rear and Side Delivery Systems

Rear or Side Beneath Counter Mount



Model 2400

Beneath Counter Mount with swivel included
Countertop/work surface not included

Part Number 2400

Optional Beneath Counter Slide 1700-140

Optional Folding Arm Assembly 1110-115-CS
(Not available with STAR-ETORQUE option)

Standard Features

Doctor Standard IC Delivery System 2400

Beneath counter swivel mount
(work surface not included)
Standard IC control head
automatic control for three handpieces
3-way air/water syringe w/tubing, 5'
4-hole handpiece tubing, 5'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
wet/dry foot control - Shadow powder coat standard
External braided umbilical, 7'
Air pressure regulator w/ filter
Master air shut-off valve

Factory Options

DELIVERY SYSTEM	MODEL CHOICE
Factory Installed STAR ETorque Electric System	STAR-ETORQUE
Handpiece Illumination System	
One handpiece (5-hole), accommodates	
1 add'l tubing	1700-021-CS
Add'l handpiece illumination tubing (5-hole)	0017-615-CS
One handpiece (6-pin), accommodates	
1 add'l tubing	1700-019-CS
Add'l handpiece illumination tubing (6-pin)	0017-603-CS
Assistant's hygiene instrumentation (control head mount)	1408-CS
Quick disconnect with flow control, 1/4" (specify air or water)	0016-455-CS
Slide mount with swivel, 7 5/8" extension	1700-140
Folding arm assembly, 8" x 12" with wall mount	1110-115-CS
Standard utility center housing & cover	1115-020
Routing valve kit with water shutoff & regulator/factory installation	1115-315-CS

Part No

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

! Often overlooked

!

Simplified Surface Treatments

Mounting plate dimensions, see page 112

Shipping weight estimate:
2400 28 lb / 13kg

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

Simplified Rear and Side Delivery Systems

Rear Beneath Counter Mount



Panel Mount is recessed in cabinet separate from duo articulating instrument delivery 2184.



Model 2184

Mounts under floating work surface.

Part Number 2184



Model 2186

Mounts under floating work surface.

Part Number 2186

Standard Features

Duo Delivery System 2184/2186

Automatic panel control for three handpieces
3-way air/water syringe w/tubing, 5'
4-hole handpiece tubing, 5'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
wet/dry foot control - Shadow powder coat standard
Central vacuum accessory package, swing mount
arm-mounted central vacuum canister
w/disposable solids collector
3-way air/water syringe w/tubing, 5'
autoclavable standard HVE lever
valve w/tubing, 5'
autoclavable saliva ejector lever valve
w/tubing, 5'
Articulating left/right instrument positioning
Self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
External umbilical, 7'
Air pressure regulator w/filter
Master air shut-off valve

Factory Options

DELIVERY SYSTEM

Factory Installed STAR ETorque Electric System

STAR-ETORQUE

Handpiece Illumination System

One handpiece (5-hole), accommodates

1 add'l tubing

1700-021-CS

Add'l handpiece illumination tubing (5-hole)

0017-615-CS

One handpiece (6-pin), accommodates

1 add'l tubing

1700-019-CS

Add'l handpiece illumination tubing (6-pin)

0017-603-CS

Additional standard autoclavable HVE position
with tubing, 5' & autoclavable oral dry cup

1400-113-CS

Hands-Free HVE, see page 103, Write in selection here

Slide mount, 7 5/8" extension (2184)

1700-142

Touchpad for 3900 chair (25 ft. Cable)

1700-081-1-CS

Standard utility center housing & cover

1115-020

! Quick disconnect for air or water - dealer installed (see page 106)

! Routing valve kit for city water - dealer installed (see page 106)

! Often overlooked

Part No

Mounting plate dimensions, see pages 120

Panel cut-out, see page 120

Shipping weight estimate:

2184 30 lb / 14 kg

2186 30 lb / 14 kg

!

Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Simplified Rear and Side Delivery Systems

Rear Beneath Counter Mount



Panel Mount is recessed in cabinet separate from duo articulating instrument delivery 2184BC.



Model 2184BC

Universal Heavy-Duty Mounting Plate 11" x 14 1/2" included (mounts to under counter)
Specify cabinet manufacturer

Part Number

2184BC



Model 2186BC

Universal Heavy-Duty Mounting Plate 11" x 14 1/2" included (mounts to under counter)
Specify cabinet manufacturer

Part Number

2186BC

Standard Features

Duo Delivery System 2184BC/2186BC

Automatic panel control for three handpieces
3-way air/water syringe w/tubing, 5'
4-hole handpiece tubing, 5'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
wet/dry foot control - Shadow powder coat standard
Central vacuum accessory package, swing mount
arm-mounted central vacuum canister
w/disposable solids collector
3-way air/water syringe w/tubing, 5'
autoclavable standard HVE lever valve w/tubing, 5'
autoclavable saliva ejector lever valve w/tubing, 5'
Articulating left/right instrument positioning
Self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
Solid work surface, 19" x 24"
► Beneath counter rigid arm assembly w/heavy duty
mounting plate, 11" x 14 1/2"
External umbilical, 7'
Air pressure regulator w/filter
Master air shut-off valve



► Beneath Counter Mount

Mounting plate dimensions, see pages 120

Panel cut-out, see page 120

Shipping weight estimate:

2184BC 85 lb / 36 kg

2186BC 85 lb / 36 kg

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

FACTORY OPTIONS	MODEL CHOICE	Part No
DELIVERY SYSTEM	MODEL CHOICE	Part No
Factory Installed STAR ETorque Electric System	STAR-ETORQUE	
Handpiece Illumination System		
One handpiece (5-hole), accommodates 1 add'l tubing	1700-021-CS	
Add'l handpiece illumination tubing (5-hole)	0017-615-CS	
One handpiece (6-pin), accommodates 1 add'l tubing	1700-019-CS	
Add'l handpiece illumination tubing (6-pin)	0017-603-CS	
Additional standard autoclavable HVE position with tubing, 5' & autoclavable oral dry cup	1400-113-CS	
Hands-Free HVE, see page 103, Write in selection here		
Quick disconnect with flow control 1/4" (for water)	0016-454-CS	
Quick disconnect with flow control 3/8" (for air)	0016-456-CS	
Storage organizer with bins & frame	1110-179-CS	
Swivel tray (requires purchase of 1110-179-CS)	1110-908-CS	
Swivel tray (mounts without storage organizer)	1110-482-CS	
Monitor mount assembly (monitor screen not included) (requires purchase of 1110-179-CS)	1300-025-CS	
Touchpad for 3900 chair (25 ft. Cable)	1700-081-1-CS	
Standard utility center housing & cover	1115-020	
Routing valve kit with water shutoff & regulator/factory installation	1115-315-CS	
Standard laminate work surface in lieu of solid surface	0023-003-SS-LAM	
Special color laminate or solid work surface upgrade △ (call for price and leadtime)	0023-003-X / 0023-003-SSV	

! Often overlooked

△ Lead time varies

!

Simplified Surface Treatments

UNIT PLASTIC FINISHES: Forest Grey

POWDER COAT PAINT FINISHES: Forest Grey

SOLID WORK SURFACE:

CORIAN®: Antarctica

WILSONART®: Arctic Melange Baja Melange

MEGANITE®: Cinder Granite Goldrush Gemstone Moonstone Platinum Granite

Rocky Road Granite Silver Mist Starry Starry Night

Simplified Rear and Side Delivery Systems

Panel Mounts



Model 2500

Part Number 2500



Optional Assistant's Hygiene
Instrumentation 1408-CS

Model 2520

Part Number 2520

Standard Features

Doctor Control 2500

Automatic control for three handpieces
3-way air/water syringe w/tubing, 7'
4-hole handpiece tubing, 7'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
wet/dry foot control - Shadow powder coat standard
Remote mount holder assembly w/control tubing, 6'
External umbilical, 6'
Air pressure regulator w/filter
Master air shut-off valve

Doctor Control 2520

Automatic control for two handpieces
3-way air/water syringe w/tubing, 7'
4-hole handpiece tubing, 7'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
non-retraction water coolant
air coolant adjustment
self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
wet/dry foot control - Shadow powder coat standard
Remote mount holder assembly w/control tubing, 6'
External umbilical, 6'
Air pressure regulator w/filter
Master air shut-off valve

Factory Options

DELIVERY SYSTEM

Handpiece Illumination System

- | | |
|---|-------------|
| One handpiece (5-hole), accommodates 1 add'l tubing | 1700-021-CS |
| Add'l handpiece illumination tubing (5-hole) | 0017-615-CS |
| One handpiece (6-pin), accommodates 1 add'l tubing | 1700-019-CS |
| Add'l handpiece illumination tubing (6-pin) | 0017-603-CS |
| Standard utility center housing & cover | 1115-020 |

Assistant's hygiene instrumentation with extended 10 1/2"
holder bar (Drain Tubing 12') standard autoclavable
HVE, saliva ejector with 2 holders and canister

- ! Quick disconnect for air or water - dealer installed (see page 106)
! Routing valve kit for city water - dealer installed (see page 106)

! Often overlooked

Part No

MODEL CHOICE

- | | |
|---|-------------|
| Handpiece Illumination System | _____ |
| One handpiece (5-hole), accommodates 1 add'l tubing | 1700-021-CS |
| Add'l handpiece illumination tubing (5-hole) | 0017-615-CS |
| One handpiece (6-pin), accommodates 1 add'l tubing | 1700-019-CS |
| Add'l handpiece illumination tubing (6-pin) | 0017-603-CS |
| Standard utility center housing & cover | 1115-020 |

Assistant's hygiene instrumentation with extended 10 1/2"
holder bar (Drain Tubing 12') standard autoclavable
HVE, saliva ejector with 2 holders and canister

1408-CS

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

!

Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

Wall/cabinet cut-out dimensions, see page 119

Shipping weight estimate:

2500 18 lb / 8 kg

2520 15 lb / 7 kg

Simplified Rear and Side Delivery Systems

Rear or Side Mount Compact Controls



Model MT600

Work surface not included

Part Number MT600



Model MT610

Work surface not included

Part Number MT610

Standard Features

Doctor Control with 3 HP

Top Shelf Mount
Beneath Shelf Mount w/Swivel

Doctor Control with 2 HP

Top Shelf Mount
Beneath Shelf Mount w/Swivel

All Units Include the Following Features:

Automatic control for three handpieces (MT600), or
Automatic control for two handpieces (MT610)
3-way air/water syringe w/tubing, 4'
4-hole handpiece tubing, 4'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
non-retraction water coolant
air coolant adjustment
self-contained clean water system w/quick
disconnect, single bottle, 1 1/2 liter
wet/dry foot control - Shadow powder coat standard
External umbilical, 12'
Air pressure regulator w/filter
Master air shut-off valve

Factory Options

DELIVERY SYSTEM

Handpiece Illumination System

One handpiece (5-hole), accommodates 1 add'l tubing

1700-021-CS

Add'l handpiece illumination tubing (5-hole)

0017-615-CS

One handpiece (6-pin), accommodates 1 add'l tubing

1700-019-CS

Add'l handpiece illumination tubing (6-pin)

0017-603-CS

Slide mount with swivel, 7 5/8" extension (order BS)

1700-140

Folding arm assembly, 8" x 12" with wall mount (order TS)

1110-115-CS

Standard utility center housing & cover

1115-020

Assistant's hygiene instrumentation with extended 10 1/2"
holder bar (Drain Tubing 12') standard autoclavable HVE,
saliva ejector with 2 holders & canister

1408-CS

! Quick disconnect for air or water - dealer installed (see page 106)

! Routing valve kit for city water - dealer installed (see page 106)

! Often overlooked

Part No

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

!

Simplified Surface Treatments

Mounting plate dimensions, see page 121

Shipping weight estimate:

MT600, MT610 20 lb / 10 kg

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

Simplified Rear and Side Delivery Systems - Carts

Mobile Carts




Model 5350

Part Number 5350



Shown with optional swivel tray,
arm mounted 1110-731



Model 5355

Part Number 5355

Standard Features

Doctor Standard IC Cart 5350

Standard IC Control Head
automatic control for three handpieces
3-way air/water syringe w/tubing, 4'
4-hole handpiece tubing, 4'
pneumatic master control block w/diaphragm
master on/off
individual handpiece pressure adjustment
individual water coolant adjustment
non-retraction water coolant
air coolant adjustment
handpiece flush
handpiece exhaust canister/oil collector
self-contained clean water system w/quick
disconnect, single bottle 1 1/2 liter
wet/dry foot control - Shadow powder coat standard
Adjustable height frame
Twin wheel hard floor casters
External 7' umbilical
Standard utility center housing and cover
Air pressure regulator w/filter
Master air shut-off valve

Duo Standard IC Cart 5355

Same as above, w/additional features:
Central vacuum canister w/disposable solids collector
autoclavable standard HVE lever valve w/tubing, 5'
autoclavable saliva ejector lever valve w/tubing, 5'

Factory Options

DELIVERY SYSTEM

Factory Installed STAR ETorque Electric System

STAR-ETORQUE

Handpiece Illumination System

One handpiece (5-hole), accommodates 1 add'l tubing

1700-021-CS

Add'l handpiece illumination tubing (5-hole)

0017-615-CS

One handpiece (6-pin), accommodates 1 add'l tubing

1700-019-CS

Add'l handpiece illumination tubing (6-pin)

0017-603-CS

Additional standard autoclavable HVE

with tubing, 5' & autoclavable oral dry cup (5355)

1400-113-CS

Hands-Free HVE, see page 103, Write in selection here

Quick disconnect with flow control, 1/4"
(specify air or water)

0016-455-CS

Swivel tray, arm-mounted

1110-731

Foot control support shelf, clamps to 2" base

0023-665

Additional umbilical (per foot)

1115-999-CS

Routing valve kit with water shutoff &
regulator/factory installation

1115-315-CS

Part No

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

!

Simplified Surface Treatments

UNIT PLASTIC FINISHES:

Forest Grey

POWDER COAT PAINT FINISHES:

Forest Grey

Equipment dimensions, see page 122

Shipping weight estimate:
5350 50 lb / 23 kg
5355 55 lb / 25 kg

Standard Chair Adapters - for 2" Post

3900/3800 fixed accommodates external umbilical 1115-195

Monitor Posts - 2" Diameter with Oilitel Bushing (used with 1300TPM)

30" with dual slot, fixed sidebox mount 1331-998
30" with slot & post clamp, fixed mount 1331-997

Light Posts - 2" with Light Bushing

35" fixed (includes post clamp)
Forest 1110-820
Pelton & Crane 1110-822

44" fixed
Forest 1110-794
Pelton & Crane 1110-785

44" fixed with slot (for monitor cord)
Forest 1110-050

44" fixed with slot (sidebox access)
Forest 1110-900
Pelton & Crane 1110-890

44" fixed with dual slot (for monitor cord & sidebox access)
Forest 1110-170

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps



250 Watt (Universal Power Supply)

Used when two or more ancillary devices are ordered with a Forest Operatory Light (excluding the 3900 soft-touch massager power source)

- Provides power for MULTIPLE electrical accessories such as: fiber-optics operatory lights, scalers, curing lights, and electric micro motors
- Size: 7.5" - 5.5" - 2.5". Easily fits in all Forest floor boxes
- Quick disconnect connectors for all of Forest ancillary products
- Available in both 230 volt and 110 volt
- Auto-resettable Klixon switch: inbuilt over current and over temperature protection, eliminating the need for manually changing fuses or circuit breakers

Please note: built in air-electric switch, power source turns on/off w/pilot air only



Handpiece Illumination System, 5 hole with NEW integrated LED bulb

Factory option

1700-021-CS

Handpiece Illumination System, 6 Pin

Factory option

1700-019-CS

- Dual channel outputs individually adjustable from 2.5VDC to 4.0VDC
- Accommodates one additional handpiece tubing
- Illumination delay (adjustable between 0-15 seconds) can be set individually for two handpieces
- Automatic activation of handpiece illumination w/the dental unit foot control
- Simple installation into any new or existing dental unit
- Indicator LED's confirms correct operation of the control module
- Powers all brands of optic handpieces and air-motors
- ISO standard handpiece terminal w/stainless steel nut
- Handpiece tubing, 7' long (4' to 5' installed)
- Power transformer, output 24VAC, 1.05 Amps
- Installation hardware and instructions



Handpiece Illumination Tubing, 5 hole with NEW integrated LED bulb

Factory option

0017-615-CS

Used to replace power optic handpiece tubing

- ISO standard power optic handpiece terminal w/stainless steel nut
- Handpiece tubing, 7' (4' to 7' factory installed)



Handpiece Illumination Tubing, 6 pin

Factory option

0017-603-CS

Used to replace power optic handpiece tubing

- ISO standard power optic handpiece terminal w/stainless steel adapter and nut
- Handpiece tubing, 7' (4' to 7' factory installed)



Standard Deluxe Air/Water Syringe with Quick-release Tip

Factory option 1109-010

with Sheathed Tubing Assembly, 5'

Sterling	1107-755
Shadow	1107-755-SH
Black	1107-755-BLK

Straight deluxe Air/Water Syringe with Quick-release Tip

Factory option 1109-025

with Sheathed Tubing Assembly, 4'

Sterling	1109-585-STG
Shadow	1109-585-SHW
Black	1109-585-BLK

Standard Deluxe One Button Syringe with Quick-release Tip

Factory option 1109-035

with Sheathed Tubing Assembly, 5'

Sterling	1109-485-STG
Shadow	1109-485-SHW
Black	1109-485-BLK



Additional HVE with Tubing 5' and Autoclavable Oral Dry Cup

Factory option 1400-113-CS

- Standard vacuum valve

Oral dry cup only 1400-010



Hands-Free HVE Chair-Mounted Flexible HVE Holder w/ HVE and 6' Tubing

Factory options

White Holder Arm w/Sterling HVE Tubing	1400-279-STG
White Holder Arm w/Grey HVE Tubing	1400-279-GRY
Black Holder Arm w/Shadow HVE Tubing	1400-279-SHW
Black Holder Arm w/Black HVE Tubing	1400-279-BLK
Oral dry cup (pictured above)	1400-010



Dual Articulating Quad Holder Assembly

Factory option	0024-227-CS
4 standard, 4 universal, White	0024-227
4 standard, 4 universal, Forest Grey	0024-227-G
4 standard, 4 universal, Designer Friendly® (specify color)	0024-227-C
8 universal, White	0024-224
8 universal, Forest Grey	0024-224-G
8 universal, Designer Friendly (specify color)	0024-224-C



Support Shelf

Factory option	0023-665
Forest Grey	0023-665-G

- Provides adequate support for use of foot control while allowing cart mobility.

Forest 6400 Packages

Forest 6400 Ordering Guides

Forest 3900 Packages

Forest 3900 Ordering Guide

Forest Fusion Packages

Forest Fusion Ordering Guide

Surface Treatments

Chairs and Stools

Lights and Monitor Mounts

Delivery Systems 7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps



Contoured Utility Center Housing and Cover

Factory option	3900-101
Forest Grey	3900-101-G
Shadow	3900-101-SH
Designer Friendly® (specify color)	3900-101-C



Deluxe Utility Center Housing and Cover

Forest Grey	1115-007-G
Shadow	1115-007-SH
Black	1115-007-BLK
White	1115-007-WH
• Rear umbilical access	

Complete Deluxe Utility Center

Forest Grey	1115-005-G
Shadow	1115-005-SH
Black	1115-005-BLK
White	1115-005-WH
• Rear umbilical access	
• Air and water pressure regulator w/filter	
• Master air and water shut-off valves	



Standard Utility Center Housing and Cover

Factory option	1115-020
Forest Grey	1115-020-G
Designer Friendly (specify color)	1115-020-C



Recessed utility center housing and cover

1115-305

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



Tray Holder Assembly, Use with 1095/1295 Rear Pivot Assistant's Instrumentation

Factory option	1110-680
Forest Grey	1110-680-G
Designer Friendly® (specify color)	1110-680-C
Includes:	
• Pivoting powder coated tray holder, 9 3/4" x 13 1/2"	
• Clamp assembly adapts to 1 1/2" horizontal or vertical post	



Slide with Swivel 7 5/8" Extension, use w/1284, 2184, 2400, MT's

Factory option (2400/MT's)	1700-140
Forest Grey (2400/MT's)	1700-140-G
Factory option (1284, 2184) (Does not include swivel bracket)	1700-142



Folding Arm Assembly, use with 2400, MT's

Factory option	1110-115-CS
Forest Grey	1110-115-G
Includes:	
• 8" folding arm with 1/2" pin and umbilical retainer	
• 4" x 4" wall mounting plate which accepts 1/2" pin	
• 12" folding arm with 1/2" pin and umbilical retainer	
• Control mount with 1/2" pin	



Universal Heavy Duty Mounting Plate

10" x 24" beneath counter	0010-084-BLK
11" x 14 1/2" beneath counter	0010-086-BLK
10" x 24" in lieu of 11" x 14 1/2"	0010-084-CS



26" Rigid with 24" Flexible Arm In Lieu of Top Post Mount

Factory option	1110-771-RM
----------------	-------------

22" Rigid with 18" Flexible Arm In Lieu of Top Post Mount

Factory option	1110-536-RM
----------------	-------------

26" Rigid With 18" Flexible Arm In Lieu Of Top Post Mount

Factory option	1110-537-RM
• Designed for use with 2" diameter post.	



Floating Work Surface Arm Assembly, 20"

Forest Grey	0023-815-G
Designer Friendly (specify color)	0023-815-C

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified

Factory Options and Parts

General Information
Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



Flow Control

Quick Disconnects

1/4" Female and Male QD with Flow Control

Factory Installed 0016-454-CS

1/4" Female QD with flow control, black/grey/shadow/white

Dealer Installed 0016-454-BLK /-G /-SH /-WH

1/4" Female and Male QD With Flow Control, Arm and Control Head Mount

Factory Installed 0016-455-CS

1/4" Female QD with Flow Control, Arm and Control Head Mount, Black /Grey

Dealer Installed 0016-455-BLK /-G

3/8" Female and Male QD With Flow Control, Factory Installed

Factory Installed 0016-456-CS

3/8" Female QD with Flow Control, Black /Grey/Shadow/White

Dealer Installed 0016-456-BLK /-G /-SH /-WH



Water Shut-Off
and Regulator

Routing Valve Kits

— Check with Local Authorities About Installation Requirements

Factory installed with water shut-off and regulator 1115-315-CS

Factory installed without water shut-off and regulator 1115-316-CS

Dealer Installed Includes:

Toggle Valve/Bracket/Tubing/Water Shut-Off/Regulator

black / grey / shadow / white 1115-310-BLK / -G / -SH / -WH

Dealer Installed Includes: Toggle Valve, Bracket, And Tubing

black / grey / shadow / white 1115-312-BLK / -G / -SH / -WH



1100-431

1100-423

High Density Bottle

1 1/2" liter with quick disconnect adapter 1100-431

1 1/2" liter with quick disconnect upgrade kit 1100-423





For use with clean water system

- Opaque bottle
- High Density Polyethylene (HDPE)



General Information - Dimensions

Dimensions

Pivot Delivery Systems



Fixed with Sidebox Delivery Systems



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


Fixed with Sidebox Delivery Systems




Fixed without Sidebox Delivery Systems



1636/1637 Basic Cuspidor



3800 Chair



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


3900 Chair




6176 Doctor's Comfort Stool



6159C Doctor's Deluxe Stool




6156C Doctor's Economy Stool




Utility Centers


Contoured




Deluxe




Standard




Recessed




6187 Assistant's Comfort Stool



6160C/6161C Assistant's Deluxe Stool



6157C Assistant's Economy Stool



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions



Chair-Mounted Light



Ceiling-Mounted Light



Track Mount Light/Monitor



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


Standard & Tall Wall & Cabinet Mounted Light Tall Wall & Cabinet Mounted Light




Wall & Cabinet Mount




Optional Cabinet Mount



1300RM



1300TPM



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide



Utility
Warranty

Sales Reps



General Information

Dimensions




7020PRO-BC Beneath Counter



7021PRO-BC Beneath Counter



- Large contoured work surface standard for 7020/7021
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"





7020/7021PRO-BC Mounting Plates

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps



General Information

Dimensions


7020PRO-WM Wall Mount




7021PRO-WM Wall Mount




- Large contoured work surface standard for 7020/7021
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"



**7020/7021 WM
Mounting Plate**



**7020/7021 WM
Leveling Mounting Board**



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide



Utility
Warranty

Sales Reps



General Information

Dimensions

7020PRO-UC Under Cabinet Floor Mount





7021PRO-UC Under Cabinet Floor Mount




- Large contoured work surface standard for 7020/7021
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"

7020/7021PRO-UC Under Cabinet Floor Mount Leveling/Mounting Plate



1300-025-CS Monitor Assembly Option





Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series

Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps



General Information

Dimensions


7020PRO-MC Duo Cart




7021PRO-MC Duo Cart



- Large contoured work surface standard for 7020/7021
- Work surface 9/16" thick
- Shown with organization features
- Height adjustable work surface, 5"



1300-025-CS Monitor Assembly Option



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


7010BC Beneath Counter




- Work surface 9/16" thick (semi-round), 7010
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly
- Height adjustable work surface, 5"




7010BC Mounting Plates 0010-086




0010-084-CS Option




7010WM Mounting Plate




7010WM Leveling/Mounting Board



7010WM Wall Mount



1300-025-CS Monitor Assembly Option



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


- Work surface 9/16" thick (semi-round), 7010
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly
- Height adjustable work surface, 5"




7010UC Under Cabinet Floor Mount




1300-025-CS Monitor Assembly Option



7010UC Under Cabinet Floor Mount Leveling/Mounting Plate



7010MC Duo Cart



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


3486BC and WM Rear Beneath Counter Mount




3486BC Heavy Duty Mounting Plate




3486WM Wall Mount Hub




3486WM Leveling/Mounting Board




3486BC and WM Mounting Arm



- Semi-round work surface standard for 2184BC/2186BC/3486
- Work surface 9/16" thick (semi-round and rectangular)
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly



3486MC Mobile Cart



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


3486UC Under Cabinet Floor Mount Assistant's




- Semi-round work surface standard for 2184BC/2186BC/3486
- Work surface 9/16" thick (semi-round and rectangular)
- Shown with optional storage organizer and swivel tray
- Shown with additional optional asst's quad holder assembly




2184 Rear Beneath Work Surface Duo Unit




2184 MOUNTING PLATE




3486UC Under Cabinet Floor Mount Leveling/Mounting Plate




1300-025-CS Monitor Assembly Option



2184/2500 Wall Cut-Out



1700-142 SLIDE MOUNT OPTION, 1284/2184




Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps


General Information

Dimensions


2186 Rear Beneath Work Surface Duo Unit




2186 Control Mounting




2184BC Beneath Counter Duo Unit




2184BC/2186BC Beneath Counter Duo Unit




2186BC Rear Beneath Counter Duo Unit




2500 PANEL MOUNT



2520 PANEL MOUNT



2520 Wall Cut-Out



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps

General Information


Dimensions

3235 Cabinet Side Delivery




Minimum cabinet opening 18" deep x 30" wide


3230 Wall Side Delivery




2400 Rear Mount Beneath Counter




2400 Mounting Plate




3230/3235 Leveling Mounting Plate




1110-115-CS Folding Arm System Option




MT600/MT610 Compact Controls



MT600/610BS Leveling/Mounting Plate



1700-140 Slide Mount With Swivel Option



Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide


Utility
Warranty

Sales Reps


General Information

Dimensions


1284 Rear, Beneath Work Surface




1286 Rear, Beneath Shelf




1700-142 Slide Mount Option, 1284/2184




1284/1286 Mounting Plate



5350/5355 Standard IC Cart



5350/5355 Work Surface



Proper maintenance will prolong the life of the dental equipment. It is important to adhere to routine practices for patient safety and avoid premature aging of the dental equipment. Following the step by step details outlined in the Forest Operator Guides included with the product will provide valuable information to achieve optimum success maintaining patient safety and the lifespan of the treatment room equipment.

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Warranty

DENTALEZ and its employees are proud of the products we provide to the dental community. We stand behind these products with a warranty against defects in material and workmanship as provided below.

In the event that you experience difficulty with the application or operation of any of our products, please contact our technical service department at our expense at (800)423-3555.

If we cannot resolve the issue by telephone, we will arrange for a representative to contact you or suggest that the product be repaired by a dealer service technician.

If product return is required, we will provide you with a Return Authorization number and shipping instructions to return the product to the proper facility. If the product is under warranty we will ask you to provide proof of purchase such as a copy of your invoice. Please be sure to include the Return Authorization number on the package you are returning. It is the responsibility of the purchaser/end user to immediately cease using the product upon discovery that the product is not in proper working order and to dispose of any defective product that is not returned to DENTALEZ.

Freight costs for product returns are the responsibility of the customer. Products under warranty will be repaired or replaced, at our sole discretion, and returned at our expense. We are not responsible for shipping damages. We will, however, help you file a claim with the freight carrier.

DENTALEZ warrants all equipment and parts to be free of defects in material and workmanship, under normal usage, under the following terms:

Forest

	Limited Warranty Period*
Patient chairs / Provider stools	5 years from date of installation
Chair / Stool upholstery	2 years from date of installation
Delivery systems	5 years from date of installation
7000 series delivery system chassis	20 years from date of installation
Handpiece illumination systems	6 months from date of installation
Power suppliers, Syringe buttons, HVE and SE valves	1 year from date of installation
Water bottles, gas springs, and handpiece tubing	2 years from date of installation
Air and Water tubing connectors (T's), Distribution block, Master control block, Flow valves	20 years from date of installation
Operatory lights	5 years from date of installation

Please note the following additional terms of our warranty and return policy:

- This warranty covers manufacturing defects only and does not cover defects resulting from abuse, accidents, misuse, improper handling, cleaning, care or maintenance. Damage, including, without limitation, discoloration of upholstery, resulting from contact, such as dye transfer, from garments worn by dental staff/team, the use of chemicals, cleaners, disinfectants or sterilization is not covered under this warranty nor is normal wear and tear. Failure to follow or observe our operating/use, maintenance and/or installation instructions voids this warranty as does the use of non-authorized parts and repairs made by an unauthorized repair facility.
- Liability is limited to repair or replacement of the defective part(s) or products at our sole discretion. All other liabilities, in particular liability for damages, including, without limitation, special, exemplary, consequential or incidental damages are excluded.**
- THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. NO EMPLOYEE, REPRESENTATIVE OR DEALER IS AUTHORIZED TO CHANGE THIS WARRANTY IN ANY WAY OR TO GRANT ANY OTHER WARRANTY.

CANCELLATION POLICY:

Standard orders: Written cancellation request required within 3 days of order entry date, otherwise a 25% restock fee will be charged.
Expedited ship date requests are non-cancellable.

CUSTOM ORDERS:

- Custom configurations are those that deviate from Forest's published product offering; these orders are non-cancellable and non-returnable.
- May be subject to an additional "Deviation" charge and/or a non-refundable engineering charge.
- Custom powder coat, custom plastic surface treatments and/or custom treatments such as multiple color combinations or custom colors may result in an upcharge of 10-15% per model and are non-cancellable and non-returnable.

WARRANTY REPLACEMENTS:

Parts replaced on a product that is in warranty will be warranted for the duration of that product's original warranty

PRODUCT RETURN:

Opened products or product returns more than a year old cannot be returned for credit. There will be a 25% (\$25.00 minimum) restocking charge on all items authorized for return.

Ultraleather is a registered trademark of Ultrafabrics, Inc., Formica is a registered trademark of The Diller Corporation, Corian is a registered trademark of DuPont, Wilsonart is a registered trademark of Wilsonart LLC., Meganite is a registered trademark of Lai, James, Naugahyde is a registered trademark of United States Rubber Products, Inc., and Active tilt is a registered trademark of Crown Seating.

Designer Friendly is a registered trademark of DTE Oregon, Inc.

Forest 6400 Packages	
Forest 6400 Ordering Guides	
Forest 3900 Packages	
Forest 3900 Ordering Guide	
Forest Fusion Packages	
Forest Fusion Ordering Guide	
Surface Treatments	
Chairs and Stools	
Lights and Monitor Mounts	
Delivery Systems 7000 Series	
Delivery Systems	
Simplified	
Factory Options and Parts	
General Information	
Equipment Warranty	
Compressors	
Dry Vacuums	
Tanks	
Wet Vacuums	
Utility Room Accessories	
Utility Quick Reference Guide	
Utility Warranty	
Sales Reps	



Connecting innovation to care

DENTALEZ is a proven leader, dedicated to providing innovative and trusted products, service, and support. We empower healthcare professionals to maximize their operations for increased patient care.

Powered by an easy-to-use and intuitive dashboard, you will know about potential utility preventative maintenance – anytime, anywhere, so you can be prepared for anything.


With Aeras Intel, you will have one less thing to worry about.

AERAS, Inc. is a separate corporation from DENTALEZ, Inc.

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps

Ramvac

UTILITY



- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps



PROTECT YOUR REVENUE WITH LESS DOWNTIME!




Sound shields are now options for all compressor models to reduce noise in your utility room.



Keep your practice running with the reliability of oil lubricated Ramvac Dry Vacuums.



Monitor ambient room conditions and maintenance needs to make sure your office is running efficiently long term.



Every Ramvac Compressor & Dry Vacuum with Aeras Intel comes standard with smart monitoring.

Register after purchase to take advantage & make your practice smarter!

AERAS, Inc. is a separate corporation from DENTALEZ, Inc.

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

UTILITY

Compressors	128	Forest 6400 Packages
Ramvac Compressors	128	Forest 6400 Ordering Guides
Dry Vacuums	131	Forest 3900 Packages
Dry Vacuums	131	Forest 3900 Ordering Guide
Badger LF	133	Forest Fusion Packages
Tanks	134	Forest Fusion Ordering Guide
15 Gallon and 30 Gallon Tanks	134	Surface Treatments
InfiniTank	134	Chairs and Stools
Wet Vacuums	135	Lights and Monitor Mounts
Barracuda	135	Delivery Systems 7000 Series
Maintenance	137	Delivery Systems
Compressor and Vacuum Maintenance	137	Simplified
Utility Room Accessories	138	Factory Options and Parts
Remote Panels	138	General Information
Master Water Control Valve	138	Equipment Warranty
Steel Frame Stands	138	Compressors
Amalgam Separator	139	Dry Vacuums
VACHECK / FLOWCHECK	139	Tanks
SlugBuster	139	Wet Vacuums
Quick Reference Guide	140	Utility Room Accessories
Utility Warranty	142	Utility Quick Reference Guide
DSX Sales Team	143	Utility Warranty



KNOW BEFORE AT NO CHARGE!




Nomenclature Key

Series Name	Max Number of Users	Max Input Power
C = Ramvac® Compressor	3 = up to 3 users	12 = 115 volts or 230v
SS = Sound Shield included	5 = up to 5 users	2 = 230 volts
	8 = up to 8 users	
	12 = up to 12 users	
	18 = up to 18 users	


Sizing Guide

Model Number	Description	Number of Users	Number of Heads	Sound Shield
C3-12	Ramvac 3 Compressor	3	2	8658-016 (x2)
C3-12SS	Ramvac 3 Compressor w/Sound Shield	3	2	Included
C5-12	Ramvac 5 Compressor	5	2	Integrated
C8-2	Ramvac 8 Compressor	8	3	Integrated
C12-2	Ramvac 12 Compressor	12	2	8625-022 (x2)
C12-2SS	Ramvac 12 Compressor w/Sound Shield	12	2	Included
C18-2	Ramvac 18 Compressor	18	3	8625-022 (x3)
C18-2SS	Ramvac 18 Compressor w/Sound Shield	18	3	Included

Ramvac Compressor without sound shield



Ramvac Compressor with 1/2 hp sound shield



NOTES:

Note for all compressors: Pressure is 85-115 PSI. High-Pressure Mode is 105-125 PSI.
Pump Up Cycle Time, 85 to 115 PSI (Sec) is 40-60.

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Compressors

Electrical Requirements

Model Number	Description	Voltage	Amps	Breaker Size	Total HP
C3-12	Ramvac 3 Compressor	115v / 230v	12 / 6.1	15 amp	1.0
C3-12SS	Ramvac 3 Compressor w/Sound Shield	115v / 230v	12 / 6.1	15 amp	1.0
C5-12	Ramvac 5 Compressor	115v / 230v	17.6 / 9	20 / 15 amp	1.5
C8-2	Ramvac 8 Compressor	230v	16	20 amp	2.25
C12-2	Ramvac 12 Compressor	230v	21	30 amp	4
C12-2SS	Ramvac 12 Compressor w/Sound Shield	230v	21	30 amp	4
C18-2	Ramvac 18 Compressor	230v	30.5	40 amp	6
C18-2SS	Ramvac 18 Compressor w/Sound Shield	230v	30.5	40 amp	6

Weight and Dimensions

Model Number	Description	Weight (lbs)	Dimensions W x D x H (in)	Dimensions w Sound Shield W x D x H (in)	Tank Size (Gallon)
C3-12	Ramvac 3 Compressor	136	30" x 21" x 26"	NA	12
C3-12SS	Ramvac 3 Compressor w/Sound Shield	136	NA	30" x 20" x 26"	12
C5-12	Ramvac 5 Compressor	128	30" x 21" x 26"	NA	12
C8-2	Ramvac 8 Compressor	155	30" x 21" x 27"	NA	12
C12-2	Ramvac 12 Compressor	230	35" x 25" x 31"	NA	20
C12-2SS	Ramvac 12 Compressor w/Sound Shield	230	NA	35" x 25" x 31"	20
C18-2	Ramvac 18 Compressor	335	41" x 29" x 32"	NA	30
C18-2SS	Ramvac 18 Compressor w/Sound Shield	335	NA	41" x 29" x 32"	30

Ramvac Compressors



Ramvac 3 Compressor

Part Number C3-12



Ramvac 3 Compressor with Sound Shield

Part Number C3-12SS



Ramvac 5 Compressor

Part Number C5-12



Ramvac 8 Compressor

Part Number C8-2

NOTES:

Note for all compressors: Pressure is 85-115 PSI. High-Pressure Mode is 105-125 PSI.
Pump Up Cycle Time, 85 to 115 PSI (Sec) is 40-60.

TABLE OF CONTENTS

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Compressors

Ramvac Compressors



Ramvac 12 Compressor

Part Number C12-2



Ramvac 12 Compressor with Sound Shield

Part Number C12-2SS



Ramvac 18 Compressor

Part Number C18-2



Ramvac 18 Compressor with Sound Shield

Part Number C18-2SS

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors**
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

See page 137 for maintenance parts.



KNOW BEFORE AT NO CHARGE!

Nomenclature Key

Series Name	Max Number of Users	Max Input Power	Tank Included
RV = Ramvac® Dry Vac	4 = up to 4 users	12 = 115v or 230v	T = tank included
	5 = up to 5 users	2 = 230v	
	7 = up to 7 users	23 = 230v 3 phase	
	10 = up to 10 users	43 = 460v 3 phase	
	12 = up to 12 users		

Sizing Guide: 1 HVE = 1 user, 1 SE = 1/5 user, Nitrous scavenger = 1/4 user**Sizing Guide**

Model Number	Description	Users	Breaker Size
RV4-12T	Ramvac 4 Dry Vacuum System (Bulldog QT 2)	4	15 amp
RV5-2T	Ramvac 5 Dry Vacuum System (Bison 3)	5	20 amp
RV5-23T	Ramvac 5 Dry Vacuum System (Bison 3)	5	15 amp
RV7-2T	Ramvac 7 Dry Vacuum System (Bison 5)	7	20 amp
RV7-23T	Ramvac 7 Dry Vacuum System (Bison 5)	7	15 amp
RV10-2T	Ramvac 10 Dry Vacuum System (Bison 7)	10	20 amp
RV10-23T	Ramvac 10 Dry Vacuum System (Bison 7)	10	15 amp
RV12-23T	Ramvac 12 Dry Vacuum System (Bison 9)	12	20 amp
RV12-43T	Ramvac 12 Dry Vacuum System (Bison 9)	12	15 amp

All Ramvac Vacuum Systems
(pump and tank) except
RV4-12T include a platform
(575056) for stacking.

Electrical Requirements

Model Number	Description	Output Power	Input Power	Max Amps
RV4-12T	Ramvac 4 Dry Vacuum System (Bulldog QT 2)	1 hp (0.745kW)	115v or 230v, 1 ph 60 Hz	9 @ 115v, 4.5 @ 230v
RV5-2T	Ramvac 5 Dry Vacuum System (Bison 3)	3 hp (2.2kW)	230v, 1 ph 60 Hz	16.0 @ 230v
RV5-23T	Ramvac 5 Dry Vacuum System (Bison 3)	3 hp (2.2kW)	230v, 3 ph 60 Hz	8.4 @ 230v
RV7-2T	Ramvac 7 Dry Vacuum System (Bison 5)	3 hp (2.2kW)	230v, 1 ph 60 Hz	16.0 @ 230v
RV7-23T	Ramvac 7 Dry Vacuum System (Bison 5)	3 hp (2.2kW)	230v, 3 ph 60 Hz	8.4 @ 230v
RV10-2T	Ramvac 10 Dry Vacuum System (Bison 7)	3 hp (2.2kW)	230v, 1 ph 60 Hz	16.0 @ 230v
RV10-23T	Ramvac 10 Dry Vacuum System (Bison 7)	3 hp (2.2kW)	230v, 3 ph 60 Hz	8.4 @ 230v
RV12-23T	Ramvac 12 Dry Vacuum System (Bison 9)	5 hp (3.7kW)	230v, 3 ph 60 Hz	13.2 @ 230v
RV12-43T	Ramvac 12 Dry Vacuum System (Bison 9)	5 hp (3.7kW)	460v, 3 ph 60 Hz	6.6 @ 460v

Weight and Dimensions

Model Number	Description	Weight (lbs) Pump only	Dimensions (Pump only) W x D x H (in)	Stacked Dimensions (Pump and Tank) W x D x H (in)
RV4-12T	Ramvac 4 Dry Vacuum System (Bulldog QT 2)	200 lbs	25" x 17" x 28"	29" x 19" x 48"
RV5-2T	Ramvac 5 Dry Vacuum System (Bison 3)	345 lbs	26" x 19" x 31"	33" x 21" x 59"
RV5-23T	Ramvac 5 Dry Vacuum System (Bison 3)	345 lbs	26" x 19" x 31"	33" x 21" x 59"
RV7-2T	Ramvac 7 Dry Vacuum System (Bison 5)	345 lbs	26" x 19" x 31"	33" x 21" x 59"
RV7-23T	Ramvac 7 Dry Vacuum System (Bison 5)	345 lbs	26" x 19" x 31"	33" x 21" x 59"
RV10-2T	Ramvac 10 Dry Vacuum System (Bison 7)	395 lbs	26" x 22" x 31"	33" x 22" x 59"
RV10-23T	Ramvac 10 Dry Vacuum System (Bison 7)	395 lbs	26" x 22" x 31"	33" x 22" x 59"
RV12-23T	Ramvac 12 Dry Vacuum System (Bison 9)	395 lbs	26" x 22" x 31"	33" x 22" x 71"
RV12-43T	Ramvac 12 Dry Vacuum System (Bison 9)	395 lbs	26" x 22" x 31"	33" x 22" x 71"

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface Treatments

Chairs and Stools

Lights and Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

Dry Vacuums

Ramvac Dry Vacuums

**Ramvac 4 Dry Vacuum System**

Part Number

RV4-12T

**Ramvac 5 Dry Vacuum System**

Part Number

RV5-2T

**Ramvac 7 Dry Vacuum System**

Part Number

RV7-2T

**Ramvac 10 Dry Vacuum System**

Part Number

RV10-2T

**Ramvac 12 Dry Vacuum System**

Part Number

RV12-23T

Description	Part Number	Description	Part Number
Vacuum Only			
Ramvac 4 Dry Vacuum (Bulldog QT 2) 4 User/115V & 230V Single Phase	RV4-12	Ramvac System (with 30-gal tank and platform)	
Ramvac 5 Dry Vacuum (Bison 3) 5 User 230V Single Phase	RV5-2	Ramvac 5 Dry Vacuum System (Bison 3) 5 User 230V Single Phase	RV5-2T30
Ramvac 7 Dry Vacuum (Bison 5) 7 User 230V Single Phase	RV7-2	Ramvac 7 Dry Vacuum System (Bison 5) 7 User 230V Single Phase	RV7-2T30
Ramvac 7 Dry Vacuum (Bison 5) 7 User 230V Three Phase	RV7-23	Ramvac 10 Dry Vacuum System (Bison 7) 10 User 230V Single Phase	RV10-2T30
Ramvac 10 Dry Vacuum (Bison 7) 10 User 230V Single Phase	RV10-2	Ramvac 12 Dry Vacuum System (Bison 9) 12 User 230V Three Phase	RV12-23T
Ramvac 10 Dry Vacuum (Bison 7) 10 User 230V Three Phase	RV10-23	Ramvac 12 Dry Vacuum System (Bison 9) 12 User 460V Three Phase	RV12-43T
Ramvac 12 Dry Vacuum (Bison 9) 12 User 230V Three Phase	RV12-23	See page 137 for maintenance parts.	
Ramvac 12 Dry Vacuum (Bison 9) 12 User 460V Three Phase	RV12-43		

NOTES:

RV5-2T, RV7-2T, RV10-2T and RV12-23T include a platform.

- Forest 6400 Packages
- Forest 6400 Ordering Guides
- Forest 3900 Packages
- Forest 3900 Ordering Guide
- Forest Fusion Packages
- Forest Fusion Ordering Guide
- Surface Treatments
- Chairs and Stools
- Lights and Monitor Mounts
- Delivery Systems 7000 Series
- Delivery Systems
- Simplified
- Factory Options and Parts
- General Information
- Equipment Warranty
- Compressors
- Dry Vacuums
- Tanks
- Wet Vacuums
- Utility Room Accessories
- Utility Quick Reference Guide
- Utility Warranty
- Sales Reps

Badger™ LF



The efficient, LubeFree dry vacuum system

- Compact LubeFree rotary vane pump with integrated 5 gallon tank
- Dental office can transition easily from wet to dry vacuum
- Open frame design means quick maintenance
- Saves 200,000 gallons of water annually
- 5 gallon polyethylene tank
- 2 year system warranty and 5 year pump warranty

Badger LF Dry Vacuum

Part Number BADGLF2

Badger Specifications

Power Requirements		Pump Specifications	
Users	4	HP	1.5 HP 1,200 rpm
Voltage	208-230 VAC, Single Phase	Weight	Unit: 175lbs Shipping: 195lb
Full Load	6.2A	Maximum Vacuum	15 in Hg
Frequency	60 Hz	Breaker Size	15amp
Exhaust Vent	2"		
Dimensions:	W25"xD18"xH24"		

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Tanks

Ramvac Tanks



15 Gallon Tank**

Part Number A575015



30 Gallon Tank**

Part Number A575030

- Store liquids and debris removed from the treatment room
- Tanks must be sized to correctly match the volume of liquids and debris removed during a vacuum operating cycle.

Specifications	15 Gallon	30 Gallon
Treatment Rooms*	1-6	7-12
W	29"	29"
D	19"	19"
H	24"	36"
Description	Part Number	
Additional Options		
Ramvac Vacuum Sizes 5-12 Platform	575056	
Suspension Set	575215	
12" Height Extender	575212	
Auxiliary Vacuum Adapter	590504	

Ramvac Infini Tanks



Ramvac Infini Tank

Part Number A575080

InfiniTanks

- Simple maintenance free design
- Drains without vacuum interruption
- Significantly reduces the size and number of separating tanks
- Completely automated to provide worry free operation

Specifications	InfiniTank	InfiniTank-High Flow
Treatment Rooms*	1-20	1-20
W	29"	55"
D	21"	27"
H	38"	62"

50 Gallon Cylindrical Tank



50 Gallon Cylindrical Tank

Part Number 590053

50 Gallon Cylindrical Tank

- Easy to Install
- Reliable
- Fiberglass construction won't corrode
- Clean-drain tapered bottom
- Clog resistant drain valve
- Omni-directional drain & pump connections

Specifications

Treatment Rooms*	13-20
Includes	2" & 3" inlet & vacuum fittings, 3 each galvanized pipe legs, and drain valve
OD	26"
H	62"

NOTES:

* This assumes 1 chair per room. For bays, let "rooms"= chairs.

** Ramvac tanks include level sensor and cable interface to Ramvac Vacuum Control with Aeras.

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps

Wet Vacuums

Barracuda™

• 100% Dual Redundancy

- Brass impeller and housing with sealed motor bearings for long life and chemical resistance
- 1 & 2 HP Water Injected Pump models available to meet a practice's specific requirements
- Optional water recirculator reduces water consumption by about 75%
- Two independent pumps to handle peak loads and back-up protection (Dual only)
- Maintains vacuum in peak work periods (Dual only)
- Provides back-up protection to minimize downtime (Dual only)
- Now with a 5 year limited warranty from date of installation.

Barracuda Single

Product Name	Barracuda CV-101FS	Barracuda CV-101FSW	Barracuda CV-102FS	Barracuda CV-102FSW
Recirculator	N/A	Included	N/A	Included
Horsepower	1x1 HP	1x1 HP	1x2 HP	1x2 HP
Volts	115 / 208-230	115 / 208-230	208-230	208-230
Amps (each motor)	13.2 / 7.2-6.6	13.2 / 7.2-6.6	11.9-10.7	11.9-10.7
Phase	Single	Single	Single	Single
RPM	3450	3450	3450	3450
Ambient Temperature	32°C / 104°F	32°C / 104°F	32°C / 104°F	32°C / 104°F
Cycle	60 Hz	60 Hz	60 Hz	60 Hz
Width	12"	12"	14"	14"
Depth	9"	9"	11"	11"
Height	16"	16"	20"	20"
Shipping Weight	68 lbs.	68 lbs.	78 lbs.	78 lbs.
Typical application	2 Users	2 Users	3 Users	3 Users
Breaker Size	20 amp	20 amp	20 amp	20 amp
Max sealed vacuum	Approximately 20-25 inches Hg Adjustable			

Barracuda Dual

Product Name	Barracuda MC-201FS	Barracuda MC-201FSW	Barracuda MC-202FS	Barracuda MC-202FSW
Recirculator	N/A	Included	N/A	Included
Horsepower	2 x 1 HP	2 x 1 HP	2 x 2 HP	2 x 2 HP
Volts	208-230	208-230	208-230	208-230
Amps (each motor)	7.5 each motor	7.5 each motor	15 each motor	15 each motor
Phase	Single	Single	Single	Single
RPM	3450	3450	3450	3450
Ambient Temperature	32°C / 104°F	32°C / 104°F	32°C / 104°F	32°C / 104°F
Cycle	60 Hz	60 Hz	60 Hz	60 Hz
Width	24"	26"	24"	26"
Depth	17"	20"	17"	20"
Height	24"	24"	24"	24"
Shipping Weight	175 lbs.	175 lbs.	215 lbs.	215 lbs.
Typical application	4 Users	4 Users	6 Users	6 Users
Breaker Size*	15 amp	15 amp	20 amp	20 amp
Max sealed vacuum	Approximately 20-25 inches Hg Adjustable			



Barracuda CV-101FS



Barracuda MC-202FSW

NOTES:

Dimensions do not include recirculator

*Each motor.

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Wet Vacuums

Barracuda Single



Barracuda CV-101FS

Part Number 64586121

Barracuda CV-101FSW

Part Number 64586121W

Barracuda CV-102FS

Part Number 64586122

Barracuda CV-102FSW

Part Number 64586122W

Barracuda Dual



Barracuda MC-201FS

Part Number 64586110D

Barracuda MC-201FSW

Part Number 6458611DW

Barracuda MC-202FS

Part Number 64586211D

Barracuda MC-202FSW

Part Number 6458621DW



NEW STYLE WITH RELAY BOX

For use with wet pumps with serial numbers starting with "V"



OLD STYLE NO RELAY BOX

For use with wet pumps with serial numbers starting with "S"

Description

Part Number

Water Recirculators

WR-100 - Water Recirculator for CV-101	64592011
WR-100 - Water Recirculator for CV-102	64592012
WR-100 - Water Recirculator for MC-201 & MC-202 Pumps	64592005
Cleaner for Water Recirculator (40oz. container with 1oz. scoop)	64568129

Water Separator

AM-200 - Air/Water or Nitrous Oxide Separator	64656200
---	----------

Maintenance / Replacement Pumps

MC-201FS 1 HP	64586811R
MC-202FS 2 HP	64586812R
MC-201FS with Water Recirculator	64586813R
MC-202FS with Water Recirculator	64586814R
MC-201 Old Style-208/230 Volt (no relay box)	64586100R
MC-202 Old Style-208/230 Volt (no relay box)	64586200R

Accessories

Top shelf with a 300 lb. Capacity for Stacking Compressors	64568202
--	----------

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Maintenance and Accessories

COMPRESSORS

Compressor		Sound Shield/Covers	Filter Kit	Dryer	Fresh Air Intake	Buck Boost
Part #	Users	Part #	Part #	Part #	Part #	Part #
C3-12	3	8658-016(X2)	8801-080	004018SP	8658-009	64637006
C3-12SS	3	Included	8801-079	004018SP	8658-014	64637006
C5-12 ACR Serial Prefix	5	8658-016(X3)	8801-080	004018SP	8658-010	64637006
C5-12 ACS Serial Prefix	5	Integrated	8801-080	004019SP	8658-009	64637006
C8-2 ACR Serial Prefix	8	8625-022 (x2)	8801-087	004020SP	8658-012	64637006
C8-2 ACS Serial Prefix	8	Integrated	8801-080	004019SP	8658-010	64637006
C12-2	12	8625-022 (x2)	8801-087	004020SP	8658-012	64637013
C12-2SS	12	Included	8801-101	004020SP	8658-019	64637013
C18-2	18	8625-022 (x3)	8801-087	005413SP	8658-013	64637013
C18-2SS	18	Included	8801-100	004020SP	8658-020	64637013
C3-2*	3	8543-022	005147SP	004018SP	005002	64637006
C4-1*	4	8543-021	005147SP	004018SP	005002	64637013
C4-2*	4	8543-022	005147SP	004018SP	005002	64637006
C4-12*	4	8658-016(X3)	8801-080	004018SP	8658-010	64637006
C4-12SS*	4	Included	8801-079	004018SP	8658-015	64637006
C5-12SS*	5	Included	8801-079	004018SP	8658-015	64637006
C6-2*	6	8543-023	005350SP	004020SP	8658-012	64637006
C7-2*	7	8543-023	005148SP	004020SP	005002	64637013
C8-2SS*	8	Included	8801-101	004020SP	8658-019	64637006
C11-2*	11	8543-024	005149SP	005413SP	005014	64637013
OSP22*	3	8543-022	005147SP	004018SP	005002	64637006
OSP13*	4	8543-021	005147SP	004018SP	005002	64637013
OSP23*	4	8543-022	005147SP	004018SP	005002	64637006
OSP24*	6	8543-023	005350SP	004020SP	005002	64637006
OSP25*	7	8543-023	005148SP	004020SP	005002	64637013
OSP28*	11	8543-024	005149SP	005413SP	005014	64637013
1022*		N/A	004015	004018SP	N/A	N/A
1013D*		N/A	004016	004018SP	N/A	N/A
1023D*		N/A	004016	004018SP	N/A	N/A
1025D*		N/A	004016	004018SP	N/A	N/A
1027D20*		N/A	004016	004019SP	N/A	N/A
1029T*		N/A	004017	004020SP	N/A	N/A
1029T30*		N/A	004017	004020SP	N/A	N/A

*Discontinued

VACUUMS

Vacuum		Service Filter Set		Oil		Buck Boost	
Part #	Users	Part #	Part #	Part #	Part #	Part #	Part #
RV4-12	4	003740		7809-001		64637006	
RV5-2	5	003750		7809-001		64637006	
RV7-2	7	003750		7809-001		64637006	
RV10-2	10	003750		7809-001		64637006	
RV12-23	12	003750		7809-001		64637006	
BADGLF2 prior to BG1307xxx	4	N/A		N/A		64637006	
BADGLF2 after BG1307xxx	4	N/A		N/A		64637006	
All Wet Vacuums	2-6	64568145		N/A		64637006	

Models 550, 1400 and 2800 use DELVAC 15W40 Pack of 4 gallons (part # 310703)

4 user vacuum requires 3 qts of oil, some require 4 qts of oil

NOTES: Please refer to your User Guide or contact Ramvac Technical Support to determine your Ramvac's oil capacity.
Remote air intakes come with sound covers, not with sound shield.

TABLE OF CONTENTS

Forest 6400 Packages	
Forest 6400 Ordering Guides	
Forest 3900 Packages	
Forest 3900 Ordering Guide	
Forest Fusion Packages	
Forest Fusion Ordering Guide	
Surface Treatments	
Chairs and Stools	
Lights and Monitor Mounts	
Delivery Systems 7000 Series	
Delivery Systems	
Simplified	
Factory Options and Parts	
General Information	
Equipment Warranty	
Compressors	
Dry Vacuums	
Tanks	
Wet Vacuums	
Utility Room Accessories	
Utility Quick Reference Guide	
Utility Warranty	
Sales Reps	

Accessories

Remote Panels



Remote Panel #2

Part Number 8-600-132

Remote Panel #4

Part Number 8-600-134

Remote Panel #1

Part Number 8-600-131

Remote Panel #3

Part Number 8-600-133

Remote Panel #5

Part Number 8-600-135

Master Water Control Valves



Master Water Control Valve

Adjustable between 1/2" and 3/4"

Part Number MWC075

Master Water Control Valve with Regulator

Adjustable between 1/2" and 3/4" with regulator

Part Number MWC075R

Master Water Control Valve 1"

Part Number MWC001

NOTES:

Master Water Control Valves now come standard with 24V Solenoid valve and transformer for easy installation to Ramvac master control panels. 115V versions are available; please call 1-866-DTE-INFO for ordering information and availability.

Description	Part Number
Replacement Cartridges	
AP110 Water Filter for 1/2" and 3/4" Master water controls (2 pack, for older WC series)	64520006
AP801 Water Filter for 1" Master water controls (2 pack, for older WC series)	64520012
Filter Cartridge Kit for new MWC series (2 pack)	8801-040

Master Water Control Valve Replacement Kits

Solenoid Valve Kit, 3/4"

8801-041

Solenoid Valve Kit, 1"

8801-042

Transformer Kit

8801-043

O-Ring Replacement Kit

8801-044

Adjustable Utility Stand



Adjustable Utility Stand, shown at 62, 50 and 40 inches.

Adjustable Utility Stand

36" W x 24" D x 40"-62" H

Part Number 4060

NOTES:

The Ramvac 18 Compressor (C18-2) with Sound Cover does not fit underneath stand. The RV12 with 30 gallon tank does not fit underneath stand 4060.

Forest 6400 Packages

Forest 6400 Ordering Guides

Forest 3900 Packages

Forest 3900 Ordering Guide

Forest Fusion Packages

Forest Fusion Ordering Guide

Surface Treatments

Chairs and Stools

Lights and Monitor Mounts

Delivery Systems 7000 Series

Delivery Systems

Simplified

Factory Options and Parts

General Information

Equipment Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room Accessories

Utility Quick Reference Guide

Utility Warranty

Sales Reps

Accessories

Amalgam Separator



NXT Hg5 Amalgam Separator

(1-10 users)

Part Number 8707-002

NXT Hg5 MINI Amalgam Separator

(1-4 users)

Part Number 8707-003

NXT Hg5 High Volume Amalgam Separator (1-10 users)

Part Number 8707-004

Description	Part Number
-------------	-------------

Additional Options

Bracket Mount to Ramvac Tank 8562-001

Vacuum Line Cleaner



SlugBuster™

SlugBuster Powder Starter Pack

Includes a 2 gallon mixing jug, two 40 oz. containers of SlugBuster powder with scoop

Part Number 900143

SlugBuster Powder 6 Pack

Includes six 40 oz. containers of SlugBuster powder with scoop

Part Number 900146

SlugBuster Liquid Single Gallon

Includes a single gallon of SlugBuster concentrate and hand pump dispenser

Part Number 900241

NXT Amalgam Separator

- Compact design from Solmetex with internal manifold, ideal for hard to fit areas
- Integrates seamlessly with all Ramvac wet and dry vacuums
- Environmentally friendly packaging
- Simple mail-back recycling program
- Meets all "EPA Dental Rule" requirements
 - Clear system, easy to inspect
 - Two-minute container change out
 - No tools necessary
- ISO 11143 certified by NSF International

NOTES:

Replacement canisters are supplied by Solmetex, to order, please contact them directly. Used collection containers must be delivered to an EPA Certified Recycler. Ramvac is not a recycler. Do NOT send used collection containers to Ramvac.

VACHECK® / FLOWCHECK®

Check the vacuum performance for your treatment room.



VACHECK

Part Number 800500



FLOWCHECK

Part Number 800501

VACHECK for Simple Pass-Fail" indication.(Patent 5,779,472)

FLOWCHECK for measurement of flow in SCFM. (Patent 6,925,886 B2)

NOTES:

Please refer to your User Guide or contact Ramvac Technical Support to determine your Ramvac's oil capacity.

Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps

Utility Room Quick Reference Guide

Part No	Description	Part No	Description	Forest 6400 Packages
Ramvac Compressors		Dry Vacuums		Forest 6400 Ordering Guides
C3-12	Ramvac 3 Compressor (Osprey 22)	BADGLF2	Badger LF Dry Vacuum	Forest 3900 Packages
C3-12SS	Ramvac 3 Compressor w/Sound Shield	Tanks		Forest 3900 Ordering Guide
C5-12	Ramvac 5 Compressor (Osprey 23)	A575015	15-Gallon Tank	Forest Fusion Packages
C8-2	Ramvac 8 Compressor	A575030	30-Gallon Tank	Forest Fusion Ordering Guide
C12-2	Ramvac 12 Compressor (Osprey 28)	A575080	Ramvac Infini Tank	Surface Treatments
C12-2SS	Ramvac 12 Compressor w/Sound Shield	590053	50 Gallon Cylindrical Tank	Chairs and Stools
C18-2	Ramvac 18 Compressor	575056	Vacuum Platform	Lights and Monitor Mounts
C18-2SS	Ramvac 18 Compressor w/Sound Shield	575215	Suspension Set	Delivery Systems 7000 Series
Ramvac Vacuums		575212	12" Height Extender	Delivery Systems
Vacuum System (with 15 gallon tank and platform)		590504	Auxiliary Vacuum Adapter	Simplified
RV4-12T	Ramvac 4 Dry Vacuum System (Bulldog QT 2) 4 User/115V & 230V Single Phase	Wet Vacuums		Factory Options and Parts
RV5-2T	Ramvac 5 Dry Vacuum System (Bison 3) 5 User 230V Single Phase	64586121	Barracuda Single, CV-101FS	General Information
RV7-2T	Ramvac 7 Dry Vacuum System (Bison 5) 7 User 230V Single Phase	64586121W	Barracuda Single, CV-101FSW	Equipment Warranty
RV7-23T	Ramvac 7 Dry Vacuum System (Bison 5) 7 User 230V Three Phase	64586122	Barracuda Single, CV-102FS	Compressors
RV10-2T	Ramvac 10 Dry Vacuum System (Bison 7) 10 User 230V Single Phase	64586122W	Barracuda Single, CV-102FSW	Dry Vacuums
RV10-23T	Ramvac 10 Dry Vacuum System (Bison 7) 10 User 230V Three Phase	64586110D	Barracuda Dual, MC-201FS	Tanks
Vacuum System (with 30 gallon tank and platform)		64586111DW	Barracuda Dual, MC-201FSW	Wet Vacuums
RV5-2T30	Ramvac 5 Dry Vacuum System (Bison 3) 5 User 230V Single Phase	64586211D	Barracuda Dual, MC-202FS	Utility Room Accessories
RV7-2T30	Ramvac 7 Dry Vacuum System (Bison 5) 7 User 230V Single Phase	64586211DW	Barracuda Dual, MC-202FSW	Utility Quick Reference Guide
RV10-2T30	Ramvac 10 Dry Vacuum System (Bison 7) 10 User 230V Single Phase	64592011	WR-100, Water Recirculator for CV-101	Utility Warranty
RV12-23T	Ramvac 12 Dry Vacuum System (Bison 9) 12 User 230V Three Phase	64592012	WR-100, Water Recirculator for CV-102	Sales Reps
RV12-43T	Ramvac 12 Dry Vacuum System (Bison 9) 12 User 460V Three Phase	64592005	WR-100, Water Recirculator for MC-201 & MC-202 Pumps	
Vacuum Only		64568129	Cleaner for Water Recirculator (40 oz. container)	
RV4-12	Ramvac 4 Dry Vacuum (Bulldog QT 2) 4 User/115V & 230V Single Phase	64656200	AM-200, Air/Water or Nitrous Oxide Separator	
RV5-2	Ramvac 5 Dry Vacuum (Bison 3) 5 User 230V Single Phase	Maintenance / Replacement Pumps		
RV7-2	Ramvac 7 Dry Vacuum (Bison 5) 7 User 230V Single Phase	64586811R	MC-201FS 1 HP	
RV7-23	Ramvac 7 Dry Vacuum (Bison 5) 7 User 230V Three Phase	64586812R	MC-202FS 2 HP	
RV10-2	Ramvac 10 Dry Vacuum (Bison 7) 10 User 230V Single Phase	64586813R	MC-201FS with Water Recirculator	
RV10-23	Ramvac 10 Dry Vacuum (Bison 7) 10 User 230V Three Phase	64586814R	MC-202FS with Water Recirculator	
RV12-23	Ramvac 12 Dry Vacuum (Bison 9) 12 User 230V Three Phase	64586100R	MC-201 Old Style - 208/230 Volt (no relay box)	
RV12-43	Ramvac 12 Dry Vacuum (Bison 9) 12 User 460V Three Phase	64586200R	MC-202 Old Style - 208/230 Volt (no relay box)	
		64568202	Top shelf with a 300 lb. capacity for stacking compressors	

Utility Room Quick Reference Guide

Part No Description

Remote Panels

- 8-600-131 Remote Panel #1
- 8-600-132 Remote Panel #2
- 8-600-133 Remote Panel #3
- 8-600-134 Remote Panel #4
- 8-600-135 Remote Panel #5

Master Water Control Valves

- MWC075 Adjustable between 1/2" and 3/4"
- MWC075R Adjustable between 1/2" and 3/4" with regulator

Master Water Control Valve 1"

Replacement Cartridges

- 64520006 AP110 Water Filter for 1/2" and 3/4" Master water controls (2 pack, for older WC series)
- 64520012 AP801 Water Filter for 1" Master water controls (2 pack, for older WC series)
- 8801-040 Filter Cartridge Kit for new MWC series (2 pack)

Replacement Kits

- 8801-041 Solenoid Valve Kit, 3/4"
- 8801-042 Solenoid Valve Kit, 1"
- 8801-043 Transformer Kit
- 8801-044 O-Ring Replacement Kit

Buck/Boost Transformers

- Facility voltage range w/o transformer
115VAC | 230VAC up to 40 AMP Circuit
- 64637013 94-98VAC | 206-212VAC
- Facility voltage range w/o transformer
115VAC | 230VAC up to 20 AMP Circuit
- 64637006 94-98VAC | 206-212VAC
- 590504 Auxillary Vacuum Adapter

Adjustable Utility Stand

- 4060 Adjustable Utility Stand
36"W x 24"D x 40"-62"H

VACHECK / FLOWCHECK

- 800500 VACHECK
- 800501 FLOWCHECK

NXT Amalgam Separator - Hg5 Models

- 8707-002 NXT Hg5 Amalgam Separator (1-10 users)
- 8707-003 NXT Hg5 MINI Amalgam Separator (1-4 users)
- 8707-004 NXT Hg5 High Volume Amalgam Separator (11-20 users)
- 8562-001 Bracket Mount to Ramvac Otter Tank

Part No Description

SlugBuster

- 900143 Powder Starter Pack - 2 gallon mixing jug, two 40oz. containers of SlugBuster powder with scoop
- 900146 Powder 6 Pack - six 40oz. containers of SlugBuster powder with scoop
- 900241 Liquid Single Gallon - single gallon of SlugBuster concentrate and hand pump dispenser
- 900243 Liquid Starter Pack - 2 gallon mixing jug, 1 gallon of SlugBuster concentrate, & hand pump dispenser
- 900246 Liquid 4 Pack Case - 4 single gallons of SlugBuster concentrate, & 4 hand pump dispensers

Compressor Maintenance

- 8801-080 Filter Maintenance Kit for Dual Voltage Compressor, C3-12 w/o Sound Shield
- 8801-079 Filter Maintenance Kit for Dual Voltage Compressor, C3-12SS, C4-12SS w/Sound Shield
- 8801-100 Filter Maintenance Kit for Compressor C8-2SS, C12-2SS
- 8801-101 Filter Maintenance Kit for Compressor C18-2SS
- 005147SP Filter Maintenance Kit for 3-4 user
- 005350SP Filter Maintenance Kit for 6 user
- 005148SP Filter Maintenance Kit for 7 user
- 8801-087 Filter Maintenance Kit for C6, C8, C12, C18
- 8658-009 Fresh Air Intake Kit, Ramvac Dual Votage C3-12
- 8658-010 Fresh Air Intake Kit, Ramvac Dual Votage C4-12
- 8658-012 Fresh Air Intake Kit for C6, C8, C12
- 8658-013 Fresh Air Intake Kit for C18
- 8658-019 Fresh Air Intake Kit for C8-2SS and C12-2SS
- 8658-020 Fresh Air Intake Kit for C18-2SS
- 005002 Fresh Air Intake Kit for 3-7 user
- 005014 Fresh Air Intake Kit for 11 user
- 004018SP Service Part, Dryer HR2 for 3 & 4 user
- 004019SP Service Part, Dryer HR3 for 7 user
- 004020SP Service Part, Dryer HR4 for 6 & 7 user
- 005413SP Service Part, Dryer HR4L for 11 user
- 004015 Maintenance Kit for 1022
- 004016 Maintenance Kit for 1013D, 1023D, 1025D & 1027D20
- 004017 Maintenance Kit for 1029T and 1029T30

Vacuum Maintenance

- 003583 Air Filter Set for 4 user
- 310705 Oil Filter Element
- 003740 Service Filter Set for 4 user
- 003700 Air Filter Set for 5-12 user
- 003750 User Service Filter Set for 5-12 user
- 7809-001 Ramvac Synthetic Oil, 6 x 1qt. bottles
- 710720 "Standard" Air Filter Set
- 310703 Pump Oil, DELVAC 15W40 pack of 4 gallon for models 550, 1400 and 2800
- 004691SP 25 Micron Filter Element (for units prior to serial number: BG1307xxx)
- 005301SP Badger Service Filter Set (for units later than serial number: BG1307xxx) (requires both filters)

Forest 6400

Packages

Forest 6400

Ordering Guides

Forest 3900

Packages

Forest 3900

Ordering Guide

Forest Fusion

Packages

Forest Fusion

Ordering Guide

Surface

Treatments

Chairs

and Stools

Lights and

Monitor Mounts

Delivery Systems

7000 Series

Delivery Systems

Simplified

Factory Options

and Parts

General

Information

Equipment

Warranty

Compressors

Dry Vacuums

Tanks

Wet Vacuums

Utility Room

Accessories

Utility Quick

Reference Guide

Utility

Warranty

Sales Reps

Ramvac Warranty

Ramvac Product Support Services

DENTALEZ and its employees are proud of the products we provide to the dental community. We stand behind these products with a warranty against defects in material and workmanship as provided below and have our own in-house repair facility to service our products. In the event that you experience difficulty with the application or operation of any of our products, please contact our customer service department at our expense at (866)DTE-INFO.

If we cannot resolve the issue by telephone, we will arrange for a representative to contact you or suggest that the product be returned to our factory for inspection.

If product return or repair is required, we will provide you with a Return Authorization number and shipping instructions to return the product to the proper facility. If the product is under warranty, we will ask you to provide proof of purchase such as a copy of your invoice. Please be sure to include the Return Authorization number on the package you are returning. Products returned without a return authorization number cannot be repaired.

Freight costs for product returns are the responsibility of the customer. Products under warranty will be repaired or replaced, at our sole discretion, and returned at our expense. Products outside the warranty limits will be repaired and returned with costs invoiced to the customer. We are not responsible for shipping damages. We will, however, help you file a claim with the freight carrier. Written repair estimates are available.

DENTALEZ warrants all equipment and parts to be free of defects in material and workmanship, under normal usage, under the following terms:

Ramvac

	Limited Warranty Period*
Ramvac Compressors	6 Years from date of installation
Ramvac Dry Vacuums	6 years from date of installation / 10 years for pump only
Badger LF Vacuum	2 years from date of installation / 5 years for pump only
Wet Vacuum Systems	5 Years from date of installation
Master Water Control	2 Years from date of installation
NXT Hg5 Amalgam Separator	Warranty from Solmetex, 2 years

Please note the following additional terms of our warranty and return policy:

- This warranty covers manufacturing defects only and does not cover defects resulting from abuse, accidents, misuse, improper handling, cleaning, care or maintenance. Damage resulting from the use of chemicals, cleaners, disinfectants or sterilization is not covered under this warranty nor is normal wear and tear. Failure to follow or observe our operating/use, maintenance and/or installation instructions voids this warranty as does the use of non-authorized parts and repairs made by an unauthorized repair facility.
- Liability is limited to repair or replacement of the defective part(s) or products at our sole discretion. All other liabilities, in particular liability for damages, including, without limitation, special, exemplary, consequential or incidental damages are excluded.
- **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. NO EMPLOYEE, REPRESENTATIVE OR DEALER IS AUTHORIZED TO CHANGE THIS WARRANTY IN ANY WAY OR TO GRANT ANY OTHER WARRANTY.**

Warranty Repairs:

Parts repaired or replaced on a product that is in warranty will be warranted for the duration of that product's original warranty.

Product Returns:

Opened products or product returns more than a year old cannot be returned for credit. There will be a 15% (\$25.00 minimum) restocking charge on all items authorized for return.

NOTES:

*Provided conditions defined in instruction manual are met. DENTALEZ is a registered trademark of DENTALEZ, Inc. Ramvac, Bison, Bulldog, FLOWCHECK, InfiniTank, Ramclean and VACHECK are registered trademarks and Badger, Barracuda, Osprey, and SlugBuster are trademarks of Ramvac Dental Products, Inc. PowerScrub is a trademark of Solmetex, LLC.

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps



DSXe - Equipment (Equipment Sales Specialists) - Forest, Ramvac, Porter, Radiance and Imagine+

Steve Sinclair

Divisional Sales Manager,
East Equipment Sales
ssinclair@dentalez.com
480.639.7064



Rob Depalma

Northeast
rdepalma@dentalez.com
860.993.3518



Mike Geren

Southeast
mgeren@dentalez.com
615.533.8385



Marc Bruckner

New York
mbruckner@dentalez.com
201.321.3906



Mark DeVore

Mid-Atlantic
mdevore@dentalez.com
410.320.6374



Troy Parr

Ohio
tparr@dentalez.com
513.550.3006



Dean Constantinou

Eastern Canada
dconstantinou@dentalez.com
647.203.1126



Phillip Eycue

North & South Carolina
peycke@dentalez.com
919.606.2811



Todd Smith

Michigan
tsmith@dentalez.com
616.581.4879



Doug Dasher

Florida
ddasher@dentalez.com
954.818.2618



Nick Carter

Divisional Sales Manager,
West Equipment Sales
ncarter@dentalez.com
817.583.3145



Mark Azurin

Southern California
mazurin@dentalez.com
562.485.8172



Ryan Behbahani

Pacific Northwest
rbehbahani@dentalez.com
206.595.0952



Joe Dardis

Minnesota
jarddis@dentalez.com
calendly.com/jarddis
651.894.2507



Brian Kennedy

North Texas
bkennedy@dentalez.com
214.663.2465



Jeff Clark

Northern California
jclark@dentalez.com
916.622.4493



John Gustafson

Chicagoland
jgustafson@dentalez.com
847.331.0195



Steve McClure

South Texas
s McClure@dentalez.com
210.445.4191



Sean Clifford

Mountain
s clifford@dentalez.com
303.335.6234




Dylan Hickey

Western Canada
dhickey@dentalez.com
778.846.5149



Heather Murray

Heartland
h murray@dentalez.com
636.342.0306





Forest 6400
Packages

Forest 6400
Ordering Guides

Forest 3900
Packages

Forest 3900
Ordering Guide

Forest Fusion
Packages

Forest Fusion
Ordering Guide

Surface
Treatments

Chairs
and Stools

Lights and
Monitor Mounts

Delivery Systems
7000 Series

Delivery Systems

Simplified

Factory Options
and Parts

General
Information

Equipment
Warranty

Compressors

Dry Vacuums

Tanks


Wet Vacuums

Utility Room
Accessories

Utility Quick
Reference Guide

Utility
Warranty

Sales Reps



DSXi - Institutional (Institutional Sales Specialists) Star, Ramvac, Forest, Columbia, Porter, Radiance and Imagine+

Dustin Stember

National Sales Manager, Institutional
dstember@dentalez.com
913.972.4963



Ashley Buechner

Mid-Atlantic/Southeast
abuechner@dentalez.com
919.397.3503



Chad Price

Southwest/Texas
cprice@dentalez.com
805.245.4315



David Armstrong

Northeast/Ohio
jarmstrong@dentalez.com
440-479-3860



Matt Mecuk

Heartland/Michigan
mecukm1@yahoo.com
323.348.8536



Jim Zann

California/Pacific Northwest
jzann@dentalez.com
916.835.2804



DENTALEZ CORPORATE
301 Lindenwood Drive, Ste.100
Malvern, PA 19355

T: 800.293.5206
F: 610.725.9894

DENTALEZ.COM

Visit DSXsales.com to find your local
representative. For Forest and Ramvac orders
and/or questions please call our Customer
Experience Representatives at 866.DTE.INFO.

2025 DENTALEZ, Inc. All rights reserved. DENTALEZ is
a registered trademark and DENTALEZ INC. Ramvac is a
registered trademark of Ramvac Dental Products Inc.

Forest 6400 Packages
Forest 6400 Ordering Guides
Forest 3900 Packages
Forest 3900 Ordering Guide
Forest Fusion Packages
Forest Fusion Ordering Guide
Surface Treatments
Chairs and Stools
Lights and Monitor Mounts
Delivery Systems 7000 Series
Delivery Systems
Simplified
Factory Options and Parts
General Information
Equipment Warranty
Compressors
Dry Vacuums
Tanks
Wet Vacuums
Utility Room Accessories
Utility Quick Reference Guide
Utility Warranty
Sales Reps